Biteblock, single use

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English

INTENDED USE

To protect the tracheal tube, laryngeal tube and laryngeal mask from occlusion during ventilation of a patient.

Clinical benefit: protection against occlusion of the ventilation tube

Patient target group: children to adults

Place of use: clinic and preclinic

INDICATIONS

- When using a tracheal tube, laryngeal tube or laryngeal mask.

No other indications are known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel.

- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU Member State (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.

- The product is intended for single use and must not be reused and/or reprocessed. The function of the product is impaired by reprocessing. Any reuse entails a potential infection hazard.

- The product is sterile (ethylene oxide).

- The product must not be used if the packaging is damaged or the expiration date has expired.

VISUAL INSPECTION

- ▶ Check the packaging for damages.
- ▶ Check the product for damages and loose particles.

The faulty product must be disposed of (see Section "Disposal").

USE

To use the bite block, select the appropriate size (#) according to the table.

| # | Tracheal Tube | Laryngeal Tube LTS-D | Laryngeal Mask |
|---|---------------|----------------------|----------------|
| 4 | 6.5 - 8.0 mm | - | #2, #2.5 |
| 5 | 8.5 - 11.0 mm | #2, #2.5 | - |
| 6 | - | #3, #4, #5 | #3, #4, #5 |

- ▶ Place the bite block to the ventilation tube in situ.

Option

- ▶ Place the elastic necktape around the neck of the patient and hook it onto the upper pin.
- ▶ Fix the ventilation tube to the bite block by looping the elastic necktape several times around the ventilation tube and bar and attaching the end of the elastic necktape to the lower pin.

SHELF LIFE

The shelf life of the product is 5 years.

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

The used or defective products must be disposed of in accordance with the applicable national and international statutory regulations.

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