








ET Tube with Stabilizer

for intubating Laryngeal Tube iLTS-D


DE	Gebrauchsanweisung	Seite 5 - 8
EN	Instructions for use.....	Page 9 - 12
CS	Návod k použití.....	Strana 13 - 16
ES	Instrucciones de uso.....	Página 17 - 20
FR	Manuel d'utilisation	Page 21 - 24
IT	Istruzioni d'uso.....	Pagina 25 - 28
PT	Instruções de utilização	Pág. 29 - 32
TR	Kullanım kılavuzu	Sayfa 33 - 36
	General symbol description.....	Page 38 - 40

Deutsch

PRODUKTBESCHREIBUNG

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inflates and deflates the cuff

Exit mark

indicates the exit of the ET Tube tip from the iLTS-D ventilation channel

Pilot balloon

with check valve and luer connection

Valve opener

Armoured Tube


Cuff

- high volume low pressure cuff
- seals the trachea
- reduces the risk of aspiration

Murphy Eye

Stabilizer

stabilises ET Tube during iLTS-D removal

INTENDED USE

The ET Tube is used for endotracheal intubation via a fiberoptic scope through the ventilation channel of the intubating Laryngeal Tube iLTS-D.

INDICATIONS

- Difficult Airway Management.
- To secure and maintain the airway in routine and emergency situations.

CONTRAINDICATIONS

- Fiberoptic view of the glottis opening is not possible despite repositioning of the iLTS-D.

Further contraindications are not known.

SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- The product must only be used by medical trained personnel with sufficient knowledge in the use of the product
- The product is not suitable for MRI!
- Maximal use of the product is 7 days.
- The product is suitable for conducting oxygen, air and anaesthetic gases. Anaesthetic gases used must not interact with the product materials (see material data).
- Laser or electrosurgical equipment may cause the product to be flammable.
- The product is for single use and must not be reprocessed. Reprocessing leads to impairment of the function of the product. Any reuse entails a potential risk of cross-contamination.
- The product is sterile (ethylene oxide). The product must not be used if the packaging is damaged or the shelf life has elapsed.

USE

- ▶ Check the size of the inserted iLTS-D and choose the suitable ET Tube size (see table "Product Specifications").
- ▶ Remove the valve opener.

Perform the following visual and functional check:

- ▶ Check the ET Tube for damage, patency and loose particles.
- ▶ Inflate the cuffs and check for leaks.
- ▶ Ensure that the ISO connector is removable.

CAUTION

- ▶ Re-attach the ISO connector slightly as afterwards the connector needs to be removed again.

Damaged products must not be used.

Have a spare tube ready and prepared for immediate use.

- ▶ Deflate the cuff completely and lubricate with a suitable water soluble lubricant.
- ▶ Place the ET Tube over a fiberoptic scope with a suitable diameter.
- ▶ Look for sufficient pre-oxygenation.
- ▶ Disconnect the iLTS-D from the ventilation system directly before insertion of the fiberoptic scope.

Insertion ET Tube:

- ▶ Insert the fiberoptic scope through the ventilation channel of the iLTS-D into the trachea.
- ▶ Reposition the iLTS-D in case that there is no view to the glottis opening. In addition to the change of depth insertion, it may be helpful to rotate the iLTS-D to the left or right, in order to optimise the fiberoptic view on the glottis. Pay attention to sufficient oxygenation.

CAUTION

Never use force to insert the ET tube! During insertion, the cuff must not be damaged.

- ▶ Railroad the ET Tube along the fiberoptic scope into the trachea (picture 1). The Murphy Eye should be positioned to the right of the patient. Once the exit mark of the ET Tube (see product description) passes the connector of the iLTS-D, it indicates that the tip of the ET Tube exits the iLTS-D. For tracheal placement the ET Tube must be advanced further.
- ▶ Inflate the cuff with the recommended pressure (see table "Product Specifications", picture 2).

CAUTION

Oxygen, air, anaesthetic gases or repositioning of the patient can change the cuff pressure.

- ▶ Check the cuff pressure regularly and adjust accordingly.
- ▶ Remove the fiberoptic scope.
- ▶ Connect the ET Tube to a ventilation system (picture 3).
- ▶ Use an appropriate procedure to check ventilation.

Removal iLTS-D:

- ▶ Deflate the iLTS-D cuffs completely.
- ▶ Disconnect the ventilation system from the ET Tube.
- ▶ Remove the ISO connector of the ET Tube temporarily (picture 4).
- ▶ Hold the iLTS-D with one hand and position the stabilizer on the ET Tube (picture 5).
- ▶ Carefully remove the iLTS-D while maintaining the ET Tube in position with the stabilizer (picture 6).

CAUTION

Do not actively push the ET Tube forward.

- ▶ Once the ET Tube is visible, remove the stabilizer, hold the ET Tube with one hand and remove the iLTS-D completely (picture 7).

CAUTION

Don't damage the inflation line of the ET Tube during iLTS-D removal.

- ▶ Firmly fix the ISO connector to the ET Tube (picture 8).
- ▶ Connect the ET Tube to the ventilation system and start ventilation (picture 9).
- ▶ Use an appropriate procedure to check ventilation.
- ▶ Fix the tube, if necessary.

PRODUCT SPECIFICATIONS

ET Tube size	I.D. 5.5 mm O.D. 8.3 mm	I.D. 7.5 mm O.D. 10.3 mm
Suitable for iLTS-D size	# 2.5 / 3	# 4 / 5
Recommended cuff pressure	22 - 32 cmH ₂ O	

SCOPE OF DELIVERY

Set	Set content			
	iLTS-D	ET Tube with Stabilizer	Syringe 60 ml	Syringe 100 ml
ET Tube with Stabilizer Set of 10 # 5.5	-	●	-	-
ET Tube with Stabilizer Set of 10 # 7.5	-	●	-	-
iLTS-D Intubation Set # 2.5 / 3	●	●	●	-
iLTS-D Intubation Set # 4 / 5	●	●	-	●

MATERIAL DATA

Tube	PVC (polyvinyl chloride), stainless steel
Cuff Inflation Line Pilot Balloon	PVC (polyvinyl chloride)
Connector Stabilizer	PP (polypropylene)
Check valve with luer connection	PVC (polyvinyl chloride), PBT (polybutylene terephthalate), nitrile, stainless steel

SHELF LIFE

Expiry date: see product label.


STORAGE AND TRANSPORT CONDITIONS

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

Used or damaged products must be disposed of in accordance with the applicable national and international legal regulations.

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