## iLTS-D

---

intubating Laryngeal Tube -

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
## PRODUCT DESCRIPTION

## INTENDED USE

The iLTS-D is a supraglottic airway with drain tube. The iLTS-D offers the additional possibility of endotracheal intubation via the ventilation channel.

## INDICATIONS

- As alternative adjunct to temporarily secure the airway during difficult airway management and cardiopulmonary resuscitation in clinical as well as preclinical environments (current guidelines of the European Resuscitation Council, American Heart Association and American Society of Anesthesiologists).
- For mechanical ventilation or spontaneous breathing of patients during short anaesthesia procedures.

## CONTRAINDICATIONS

- Complete obstruction of the upper airway.
- In elective use: patients who have not fasted as well as patients with retained gastric contents.

No other contraindications are known.

## SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- The product must only be used by trained medical personnel with adequate skills in the use of the product.
- During MRI treatment, the pilot balloon of the tube can cause artefacts around the area where it is positioned. A non-clinical test verified that the tube can be used in an MRI environment.
- The product is suitable for conducting oxygen, air and anaesthetic gases. Anaesthetic gases used must not interact with the product materials (see material data).
- Laser or electrosurgical equipment may cause the product to be flammable.
- Tube and syringe (optional accessories) are for single use and must not be reprocessed. Reprocessing leads to impairment of the function of the products. Any reuse entails a potential risk of cross-contamination.
- The tube is sterile (ethylene oxide). The product must not be used if the packaging is damaged or the shelf life has elapsed.

## USE

- ▶ Choose the correct tube size (see table "Product Specifications").
- ▶ Remove the valve opener.

### Perform the following visual and functional check:

- ▶ Check the tube for damage, patency and loose particles.
- ▶ Inflate the cuffs and check for leaks.  
Damaged products must not be used!  
Have a spare tube ready and prepared for immediate use.
- ▶ Deflate the cuffs completely (picture 1).
- ▶ Lubricate tip and posterior wall of the tube with a suitable water-soluble lubricant.
- ▶ Pre-oxygenation is advisable in all cases.
- ▶ Head position  
The ideal head position for insertion is the "sniffing position". However, the tube's shape also allows it to be inserted with the head in a neutral position.

## CAUTION

Achieve an appropriate depth of anaesthesia before inserting the tube!

## Insertion:

- ▶ Open the mouth with the free hand and apply chin lift (jaw thrust) to create space in the pharynx (lift base of the tongue – picture 2).
- ▶ Hold the tube like a pen in the area with the black lines (teeth marks).
- ▶ Insert the tube into the mouth and advance the tube tip along the base of the tongue in direction of the oesophagus (while advancing tilt the tube towards the nose – picture 3). The tongue must not fold back during the downward movement of the tube.
- ▶ Insert the tube without resistance to the upper teeth mark of the respective size (picture 4).

**CAUTION**

Never use force to insert the tube! During insertion, the cuffs must not be damaged. Resistance during insertion can indicate a malposition of the tube.

- ▶ Use a syringe to inflate the cuffs with the recommended volume (see table "Product Specifications" - picture 5).  
Due to the high-volume low-pressure cuffs, a seal can already be achieved with low cuff pressures (< 60 cmH<sub>2</sub>O).
- ▶ Connect the tube to a ventilation system (picture 6).

**CAUTION**

Oxygen, air, anaesthetic gases or repositioning of the patient can change the cuff pressure.

- ▶ Check regularly for the airway leak pressure and adjust accordingly.  
Excessive cuff pressure can lead to swelling of the tongue.
- ▶ Use an appropriate procedure to check ventilation and exclude gastric insufflation.
- ▶ If no sufficient ventilation can be verified, withdraw the tube with inflated cuffs between the teeth marks until ventilation is easy and free flowing (large tidal volume with minimal airway pressure).
- ▶ Check the cuff pressure again and adjust to the Airway Leak Pressure.
- ▶ Fix the tube, if necessary.

The drain tube (see table) allows the insertion of a gastric or suction catheter.

The ventilation channel allows endotracheal intubation via a fiberoptic scope. For intubation, it is recommended to use the VBM ET Tube for iLTS-D. If ET Tubes of other manufacturers are used, ensure that the ET Tube connector is removable and the outer diameter is not exceeded (see table "Product Specifications"). VBM does not take responsibility for the use of unsuitable ET Tubes.

Further information regarding intubation with the VBM ET Tube can be found in the instructions of use of the ET Tube.


**Removal:**

- ▶ Deflate the cuffs completely (picture 7).

Once the tube is in the correct position, it is well tolerated until the return of the protective reflexes.

Remove the tube in an area where suction equipment and the ability for rapid intubation are present.

**PRODUCT SPECIFICATIONS**

iLTS-D Tube size	# 2.5 / 3	# 4 / 5
Colour coding / teeth marks	# 2.5 orange # 3 yellow	# 4 red # 5 purple
Size	# 2.5 125 - 150 cm # 3 < 155 cm	# 4 155 - 180 cm # 5 > 180 cm
Recommended cuff volume	# 2.5 40 - 45 ml # 3 50 - 60 ml	# 4 70 - 80 ml # 5 80 - 90 ml
Maximum cuff pressure	60 cmH <sub>2</sub> O	
Drain tube	16 Fr	18 Fr
Max. ET Tube size	I.D. 6.5 mm (max. O.D. 9.0 mm)	I.D. 8.0 mm (max. O.D. 11.0 mm)
Recommended VBM ET Tube	I.D. 5.5 mm	I.D. 7.5 mm
Outer diameter iLTS-D 	A) 23.5 mm B) 15.5 mm	A) 29.5 mm B) 19 mm

## SCOPE OF DELIVERY

Set	Set content			
	iLTS-D	ET-Tube with Stabilizer	Syringe 60 ml	Syringe 100 ml
iLTS-D Single Set # 2.5 / 3	●	-	●	-
iLTS-D Single Set # 4 / 5	●	-	-	●
iLTS-D Intubation Set # 2.5 / 3	●	●	●	-
iLTS-D Intubation Set # 4 / 5	●	●	-	●
iLTS-D Set of 10 # 2.5 / 3	●	-	-	-
iLTS-D Set of 10 # 4 / 5	●	-	-	-

## MATERIAL DATA

Tube Cuffs Inflation line Pilot Balloon	PVC (polyvinyl chloride)
Connector	PC (polycarbonate)
Check valve and luer connection	PVC (polyvinyl chloride), PBT (polybutylene terephthalate), nitrile, stainless steel
Syringe	PP (polypropylene), IR (isoprene)

### SHELF LIFE

Expiry date: see product label.


### STORAGE AND TRANSPORT CONDITIONS

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

### DISPOSAL

Used or damaged products must be disposed of in accordance with the applicable national and international legal regulations.

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