
Catheter Mount Tubing, reusable


	TECHNICAL DATA	3
DE	GÄNSEGURGEL, WIEDERVERWENDBAR Gebrauchsanweisung	4 - 6
EN	CATHETER MOUNT TUBING, REUSABLE Instructions for use	7 - 9
BG	ГОФРИРАНА ТРЪБА, ЗА МНОГОКРАТНА УПОТРЕБА Инструкции за употреба	10 - 12
CS	SPOJKA KATÉTRU, OPAKOVANĚ POUŽITELNÁ Návod k použití	13 - 15
DA	KATETERSLANGE, GENANVENDELIG Bruganvisning	16 - 18
EL	ΣΩΛΗΝΩΣΗ ΤΟΠΟΘΕΤΗΣΗΣ ΚΑΘΕΤΗΡΑ, ΕΠΑΝΑΧΡΗΣΙΜΟΠΟΙΟΥΜΕΝΗ Οδηγίες χρήσης	19 - 21
ES	TUBO CORRUGADO, REUSABLE Instrucciones de uso	22 - 24
ET	Y-ÜHENDUSDETAIL, TAASKASUTATAV Kasutamisyjuhís	25 - 27
FI	HENGITYSLETKUN YHDISTÄJÄ, UDELLEENKÄYTETTÄVÄ Käyttööhje	28 - 30
FR	CATHÉTER MOUNT, RÉUTILISABLE Manuel d'utilisation	31 - 33
HR	REBRASTI NASTAVAK VENTILACIJSKOG CRIJEVA, ZA VIŠEKRAТNU PRIMJENU Upute za uporabu	34 - 36
HU	GÉGECSŐ, ÚJRAFELHASZNÁLHATÓ Használati utasítás	37 - 39
IT	CATETERE MOUNT, RIUTILIZZABILE Istruzioni per l'uso	40 - 42
LT	TARPINIS PRIJUNGIMO VAMZDELIS, DAUGKARTINIS Naudojimo instrukcija	43 - 45
LV	TRAHEOSTOMAS CAURULE, IZMANTOJAMA ATKĀRTOTI Lietošanas instrukcija	46 - 48
NL	ZWANENHALS, HERBRUIKBAAR Gebruiksaanwijzing	49 - 51
NO	KATETERKOBLING, TIL GJENBRUK Bruksanvisning	52 - 54
PL	ADAPTER TYPU GĘSIA SZYJA, WIELOKROTNEGO UŻYTKU Instrukcja użycia	55 - 57
PT	TUBULAÇÃO DE MONTAGEM DE CATETER COM CONECTORES, REUTILIZÁVEL Instruções de utilização	58 - 60
RO	MONTURĂ FLEXIBILĂ PENTRU CATETER, REUTILIZABILĂ Instrucţiuni de utilizare	61 - 63
RU	ПЕРЕХОДНИК, МНОГОРАЗОВЫЙ Инструкция по применению	64 - 66
SK	PREDLŽOVACIA TRUBICA, OPAKOVANE POUŽITELNÁ Návod na použitie	67 - 69
SL	CEV ZA VSTAVITEV KATETRA, ZA PONOVARNO UPORABO Navodila za uporabo	70 - 72
SV	FÖRBINDESELSESLANG, ÅTERANVÄNDBAR Bruksanvisning	73 - 75
TR	KATETER MONTAJ BORUSU, TEKRAR KULLANILABİLİR Kullanım kılavuzu	76 - 78
	SYMBOL DESCRIPTION	79 - 80

TECHNICAL DATA


REF 65-50-130, 65-50-200, 65-55-130, 65-55-200

REF 65-60-130, 65-60-200, 65-80-130, 65-80-200

REF 65-68-130, 65-68-200, 65-88-130, 65-88-200

REF	A	B	C	D	E	F
65-50-130	13 cm	17 cm	10 mm	15 mm I.D. / 22 mm O.D.	15 mm O.D.	-
65-50-200	20 cm	24 cm	10 mm	15 mm I.D. / 22 mm O.D.	15 mm O.D.	-
65-55-130	13 cm	17 cm	10 mm	15 mm I.D. / 22 mm O.D.	22 mm I.D.	-
65-55-200	20 cm	24 cm	10 mm	15 mm I.D. / 22 mm O.D.	22 mm I.D.	-
65-60-130	13 cm	18 cm	10 mm	15 mm I.D. / 22 mm O.D.	15 mm O.D.	-
65-60-200	20 cm	25 cm	10 mm	15 mm I.D. / 22 mm O.D.	15 mm O.D.	-
65-80-130	13 cm	18 cm	10 mm	15 mm I.D. / 22 mm O.D.	22 mm I.D.	-
65-80-200	20 cm	25 cm	10 mm	15 mm I.D. / 22 mm O.D.	22 mm I.D.	-
65-68-130	13 cm	18 cm	9 mm	15 mm I.D. / 22 mm O.D.	15 mm O.D.	5 mm
65-68-200	20 cm	25 cm	9 mm	15 mm I.D. / 22 mm O.D.	15 mm O.D.	5 mm
65-88-130	13 cm	18 cm	9 mm	15 mm I.D. / 22 mm O.D.	22 mm I.D.	5 mm
65-88-200	20 cm	25 cm	9 mm	15 mm I.D. / 22 mm O.D.	22 mm I.D.	5 mm

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English

INTENDED USE

The Catheter Mount Tubing is used to connect a breathing system to a patient's airway.

Clinical benefit: The Catheter Mount Tubing facilitates ventilation due to the extension of the breathing system.

Patient target group: Children, adults

Place of use: Hospital and prehospital

INDICATIONS

- To enhance flexibility of the breathing system during ventilation of a patient.

No other indications are known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU Member State (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- The product is not intended for infants.
- The connectors of the Catheter Mount Tubing must not be disassembled.
- Before every use, the product must be inspected visually for damages (cracks, breakage, etc.) (see "Checks"). A faulty product must not be used.
- The maximum application time is 30 days. After this time, the product must be replaced and reprocessed.
- The product is suitable for conducting oxygen, air, anaesthetic gases and their gas mixtures.
- Due to a possibility to cause static charge, the product must not be used with or in striking distance of flammable anaesthetic gases.
- The reusable products are supplied non-sterile and must be cleaned, disinfected and sterilised prior to every use (see "Reprocessing").

KEY FOR TECHNICAL DATA

A	Length of connecting tubing
B	Length of Catheter Mount Tubing
C	Minimum inner diameter
D	Connector to the patient
E	Connector to the breathing system
F	Bronchoscope / suction cap

USE

WARNING

The user must ensure efficient ventilation during application.

- Connect the Catheter Mount Tubing (E) to the breathing system.
- Connect the Catheter Mount Tubing (D) to the patient's airway (e. g. face mask or tube).

- After use, carefully disconnect the Catheter Mount Tubing from the breathing system and patient's airway.

The manufacturer offers an uncoupling wedge (REF 60-60-999) to facilitate disconnection of stuck connectors.

Option: Bronchoscope and suction

- A bronchoscope or suction catheter can be inserted via the bronchoscope / suction cap (5 mm) of the Catheter Mount Tubing (F).
- After use, close the bronchoscope / suction cap.

REPROCESSING

(CLEANING, DISINFECTION, STERILISATION)

GENERAL INFORMATION

- The manufacturer is not responsible for damages resulting from improper reprocessing.
- The user is obliged to validate their procedures or the devices and accessories and to comply with the validated parameters during every reprocessing.
- It is recommended to use an automated procedure due to a significantly higher efficiency of cleaning and disinfection.
- The effectiveness has been verified by an independent and accredited test laboratory.
- To achieve effective reprocessing, heavy soiling must not be allowed to dry on the product and must be removed immediately after use.
- Silicone products must not come into contact with oil and grease.

WARNING

The connectors of the Catheter Mount Tubing must not be disassembled.

CLEANING / DISINFECTION

Automated cleaning / disinfection

Notes regarding the use of washer disinfectors

- Use de-ionised water.
- Use a washer disinfectant in compliance with ISO 15883.
- For the validation a Miele G7882 equipped with cabinet E435/3 was used.

Notes regarding the use of cleaning agents and disinfectants

- When alkaline cleaning agents are used, the agent must be subsequently neutralised.
- Do not use drying agents.

Procedure validated by the manufacturer

- Connect products to the injector nozzle (remove plug and suction cap).

CAUTION

Ensure effective rinsing of all hollow spaces.

- Start program with the following parameters:
 - Pre-rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 1 min.
 - Cleaning at 55 °C (131 °F) with an exposure time of 5 min using de-ionised water and "Sekumatic® ProClean" cleaning agent (dose: 0.5 % (5 ml/l)).
 - Neutralisation with "Sekumatic® FNZ" at 20 °C (68 °F) with an exposure time of 2 min (dose: 0.1 % (1 ml/l)).
 - Rinse with de-ionised water (max. 100 CFU/ml) at 20 °C (68 °F) and an exposure time of 2 min.

- e. Thermal disinfection at 93 °C (199 °F) for 5 min with de-ionised water.
 - f. Dry at 100 °C (212 °F).
3. Check for visible contamination. If necessary, repeat the reprocessing.
 4. Check the products according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

Manual cleaning / disinfection

Notes regarding the use of cleaning agents and disinfectants

- Prepare the cleaning and disinfection solutions prior to each manual cycle.
- Use de-ionised water.

Procedure validated by the manufacturer

1. Prepare a 2 % (30 ml/l) cleaning and disinfection solution of "Sekusept® Aktiv" with de-ionised water at 20 °C (68 °F). Stir several times in the next 15 min. After 15 min the cleaning and disinfection solution is ready for use.
2. Clean the products with soft sponges in the cleaning and disinfection solution, the products must be completely immersed for at least 20 sec. Areas difficult to access must be cleaned with soft brushes. Hard brushes and other materials which may damage the surface of the products must not be used. Clean the lumen by stroking soft nylon brushes at least 10 times in and out the lumen.
3. Place the products into the cleaning and disinfection solution and leave for 15 min. Ensure all hollow spaces are covered with solution.
4. Remove cleaning and disinfection solution by rinsing thoroughly with de-ionised water (max. 100 CFU/ml). Remaining residues can shorten the product shelf life or can lead to material damage.
5. Dry the products. Avoid any accumulation of water.
6. Check for visible contamination. If necessary, repeat the reprocessing.
7. Check the products according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

CHECKS

After disinfection, the products must be inspected for damages (cracks, breakage, etc.), foreign bodies and consistency. Faulty products must be disposed of (see "Disposal").

PACKAGING

After the cleaning and disinfection procedure, the products must be packed into a sterile barrier system which is suitable for steam sterilisation. The sterile barrier system must be in compliance with ISO 11607-1. For the validation a single pouch sterile barrier system was used.

STERILISATION

Notes regarding sterilisation

- Products or sterilisation packages must be protected from mechanical damage.
- The sterilizer should meet the requirements of DIN EN 285 or DIN EN 13060.

Sterilisation validated by the manufacturer

- The products must be steam sterilised by a fractionated vacuum process.
- Exposure time is 5 min with a sterilisation temperature of 134 °C (273 °F).

- Drying time is 10 min.
- For the validation a Systec HX-320 autoclave was used.

REUSABILITY

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use. Considering the maximum shelf life of 5 years, reusable products can be reprocessed up to 100 times if appropriately used and under adherence of the cleaning, disinfection and sterilisation instructions. Every further reuse is the responsibility of the user (see section "Checks").

WARNING

By using the products on patients with suspected prion disease there may be a high risk of transmission. In such case the physician decides if the product is either disposed of (see section "Disposal") or reprocessed according to national directives.

SHELF LIFE

Product shelf life is 5 years and can be reprocessed within the shelf life up to 100 times.

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must have run through the complete reprocessing process to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

DISPOSAL

The product must be disposed of in accordance with the applicable national and international legal regulations.

MATERIAL DATA

Connecting tubing, Bronchoscope / suction cap	Silicone
Connectors	Polysulfone (PSU)
O-ring	Ethylene propylene diene rubber (EPDM)

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