





Classic Infusor

Intended to be left blank.

DE	Classic Infusor Gebrauchsanweisung 5 - 6	LT	Classic Infusor Naudojimo instrukcija 24 - 25
EN	Classic Infusor Instructions for use 6 - 7	LV	Classic Infusor Lietošanas instrukcija 26 - 27
BG	Classic Infusor Инструкции за употреба 8 - 9	NL	Classic Infusor Gebruiksaanwijzing 27 - 28
CS	Classic Infusor Návod k použití 9 - 10	NO	Classic Infusor Bruksanvisning 29 - 30
DA	Classic Infusor Bruganvisning 11 - 12	PL	Classic Infusor Instrukcja użycia 30 - 31
EL	Classic Infusor Οδηγίες χρήσης 12 - 13	PT	Classic Infusor Instruções de utilização 32 - 33
ES	Classic Infusor Instrucciones de uso 14 - 15	RO	Classic Infusor Instrucțiuni de utilizare 33 - 34
ET	Classic Infusor Kasutamisinjuhised 15 - 16	RU	Classic Infusor Инструкция по использованию 35 - 36
FI	Classic Infusor Käyttöohje 17 - 18	SK	Classic Infusor Návod na použitie 36 - 37
FR	Classic Infusor Manuel d'utilisation 18 - 19	SL	Classic Infusor Navodila za uporabo 38 - 39
HR	Classic Infusor Upute za uporabu 20 - 21	SV	Classic Infusor Bruksanvisning 39 - 40
HU	Classic Infusor Használati utasítás 21 - 22	TR	Classic Infusor Kullanım kılavuzu 41 - 42
IT	Classic Infusor Istruzioni per l'uso 23 - 24		Symbol Description 43 - 44

Deutsch

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English

INTENDED USE

To pressurise fluid bags by external pressure.

Clinical benefit:

- Rapid compensation of fluid/volume deficiency
- Irrigation of the operation field
- Support for continuous blood pressure monitoring

Patient target group: all patients with a medical indication.

Place of use: hospital and prehospital

INDICATIONS

Pressure infusion, irrigation purposes, transfusion procedures and invasive pressure monitoring.

No other indications are known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

-  Read and follow the instructions for use carefully before using the product and keep them for future reference.
-  This product must only be used by medically trained personnel.

- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.

- Before every use, the product must be inspected visually for damage (cracks, breakage, etc.) (see Section "Visual Inspection"). A faulty product must not be used.

- Application time: maximum of 5 days, depending on the indication.

- It is not permitted to make any changes to the product.

- Use original VBM individual parts only. Using individual parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.

- The product must be disinfected prior to use on any new patient (see Section "Wipe disinfection").

- The product is not sterile.

- The product is not suitable for MRI.

- The product must not be used if the expiration date has expired.

PRODUCT DESCRIPTION

a - Hanger to hang up on the IV-stand

b - Hook for fixation of the fluid bag

c - Netting

d - Three-way stopcock

e - Hand inflator

f - Zipper

VISUAL INSPECTION

► Check the product for damage.

A faulty product must be disposed of (see "Disposal").

APPLICATION

NOTE

- For optimal pressure transmission the fluid bag must be covered completely by the netting (c) of the Classic Infuser.
- The material expands during inflation, which may initially cause a slight drop in pressure.
- To avoid damaging the Classic Infuser, do not exceed the maximum pressure of 300 mmHg.
- The pressure drop due to emptying of the fluid bag must be monitored constantly and readjusted if necessary.

Classic Infusor 500 ml and 1000 ml (Figs. 1 and 2)

- ▶ Insert the fluid bag from below and fix it to one of the two hooks (b) depending on the shape of the bag.
- ▶ Hang the Classic Infusor on the IV-stand.
- ▶ Adjust the Classic Infusor to the desired pressure using the hand inflator (e) and the three-way stopcock (d).

Classic Infusor 3000 / 5000 ml (Figs. 3 and 4)

- ▶ Open the zipper (f).
 - ▶ Insert the fluid bag from above.
- The connections of the fluid bag must face downwards.
- ▶ Close the zipper.
 - ▶ Adjust the Classic Infusor to the desired pressure using the hand inflator (e) and the three-way stopcock (d).

Connect the Classic Infusor with Luer lock (REF 58-01-300) to the automatic pressure infusor for inflation.

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- Inflation ①
- Deflation ②
- Airtight ③

AFTER USE

- ▶ Deflate the Classic Infusor via the three-way stopcock (d) or via the automatic pressure infusor (3000 / 5000 ml only).
- ▶ Remove the fluid bag.

Before inserting the new fluid bag, the Classic Infusor must be fully deflated.

WIPE DISINFECTION

WARNING

The risk of contamination increases if the product is heavily soiled or if a bioburden is suspected. In this case the physician decides if the product is either disposed of (see Section "Disposal") or disinfected according to national directives.

CAUTION

- The product must not be mechanically or manually reprocessed or sterilised.
- Do not immerse the product/device in liquids.
- In case of unsuccessful wipe disinfection or heavy contamination, the product must be disposed of as specified in "Disposal".

Wipe disinfection must be performed with commercially available surface disinfectants based on alcohol, QACs (quaternary ammonium compounds) or hydrogen peroxide. When choosing a product for disinfection, a disinfectant with appropriate ranges of action must be used: bactericidal, yeasticidal and virucidal. After wipe disinfection, inspect the product for visible contamination. If necessary, repeat wipe disinfection. After the wipe disinfection, check the product (see "Visual inspection").

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must be cleaned and disinfected to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture and can be reprocessed 100 times within the period of the shelf life.

Use-by date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

PRODUCT SPECIFICATIONS

REF	Fluid bag	Model	Weight (incl. packaging)	Packaging unit
58-02-050	500 ml	with hand inflator, three-way stopcock and colour coded pressure indicator	380 g (13.4 oz)	1
58-08-050	500 ml	with hand inflator, three-way stopcock and manometer	430 g (15.2 oz)	1
58-02-100	1000 ml	with hand inflator, three-way stopcock and colour coded pressure indicator	400 g (14.1 oz)	1
58-08-100	1000 ml	with hand inflator, three-way stopcock and manometer	450 g (15.9 oz)	1
58-01-300	3000 / 5000 ml	with female Luer lock	500 g (17.6 oz)	1
58-08-300	3000 / 5000 ml	with hand inflator, three-way stopcock and manometer	630 g (22.2 oz)	1

MATERIAL DATA

REF	58-02-050 58-02-100	58-08-050 58-08-100	58-01-300	58-08-300
Fabric	PA (TPU laminated)			
Netting, fleece disc, connector Luer lock	PA			
Hook	HDPE			-
Loop	PET			
Eyelet	Nickel-plated brass			
Tubing	Silicone			
Connector	TPU			
Piping	PA			
Piping netting	-		PA	
Zipper	-		PES	
Inflation bulb	PVC		-	PVC
Connecting tubing	TPU		-	TPU
Three-way stopcock	PC / HDPE		-	PC / HDPE
Colour coded pressure indicator	ABS	-	-	-
Manometer glass	-	PC	-	PC
Manometer protective ring	-	EPDM	-	EPDM

REPLACEMENT PARTS (FIG. 6)

REF	Description	Weight (incl. packaging)	Packaging unit
20-17-301	Manometer 300 mmHg for fixation to the pressure infusion cuff	75 g (2.7 oz)	1
51-02-666	Hand inflator with three-way stopcock and colour coded pressure indicator	82 g (2.9 oz)	1
51-03-666	Hand inflator with three-way stopcock	61 g (2.2 oz)	1

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