
Intended to be left blank.

Intended to be left blank.

ScalpelCric

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INTENDED USE

Surgical cricothyrotomy set for securing the airway in case of upper airway obstruction or as ultima ratio if all other attempts to ventilate the patient have failed. ScalpelCric can be used with resuscitation bags or other standard ventilation systems.

Clinical benefit: Ventilation of a patient

Patient target group: Adults

Location of use: Hospital and prehospital including military use

INDICATION

- Emergency airway access via the cricothyroid membrane. Life-threatening dyspnoea that cannot be controlled in any other way.

CONTRAINDICATION

- Not known.

SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product and keep them with the product.

- The product must only be used by a physician or by authorized and trained medical personnel with adequate skills in the use of the product.
- Only use the product for the patient target group.
- Only use the products for patients for whom a 6 mm tube is suitable for cricothyrotomy.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before use, the product must be inspected acc. to a functional check (see "Functional check"). A faulty product must not be used.
- The product serves to temporarily secure the airway. If possible it must be replaced by a permanent airway. Maximum recommended period of application: 4 hours.
- Tube and extension tubing are suitable for conducting oxygen and air.
- The enclosed extension tubing is not suitable for use with mechanical ventilation systems.
- Avoid damage to the cuff of the tube by sharp-edged objects.
- Do not kink the inflation line of the tube.
- Use the scalpel, bougie and tube with appropriate force.
- It is not permitted to make any changes to the product.
- The product is not suitable for MRI.
- The product is intended for single use and must not be reused and/or reprocessed. The function of the product is impaired by reprocessing. Any reuse entails a potential infection hazard.
- The product is sterile (ethylene oxide).
- The product must not be used if the packaging is damaged or the expiration date has expired.

FUNCTIONAL CHECK

- ▶ Check the products for damage and loose particles.
- ▶ If required, shape the bougie to the desired form.

▶ Inflate the cuff with the syringe and check for leakage.

A faulty product must be disposed of (see "Disposal").

▶ Remove all air from the cuff with the syringe.

USE

▶ Hyperextend the patient's head.

▶ **1** Locate the cricothyroid membrane by the palpable depression between the thyroid cartilage and the cricoid cartilage and stabilise this site with index finger and thumb.

If the cricothyroid membrane is **palpable**, proceed as described in **2a**. If the cricothyroid membrane is **impalpable**, proceed as described in **2b**.

▶ **2a** Make a horizontal stab incision through the skin and cricothyroid membrane using the scalpel. Proceed as described in **3**.

▶ **2b** Make a midline vertical skin incision using the scalpel.

▶ Separate tissues with fingers of both hands.

▶ Identify and stabilize the cricothyroid membrane again.

▶ Make a horizontal stab incision through the cricothyroid membrane.

▶ **3** Keep the scalpel perpendicular to the skin and turn the sharp edge 90° caudally.

▶ **4** Maintain gentle traction, pulling the scalpel towards you in order to open the incision and to keep it open for the insertion of the bougie.

▶ Slide the bougie down the scalpel blade into the trachea.

▶ **5** Rotate the bougie caudally and advance gently up to 10 - 15 cm into the trachea (according to the distance marks).

▶ Remove the scalpel.

▶ **6** Hold the bougie and railroad the tube gently over the bougie into the trachea, until the flange is flush with the patient's neck. Lubricate the tube if required.

▶ Remove the bougie. Make sure that the position of the tube is not changed.

▶ If required the tube can be repositioned after loosening the flange.

The distance marks on the tube indicate the insertion depth.

▶ Then tighten the flange to prevent axial dislocation of the tube.

▶ Inflate the cuff to the minimum pressure required to seal the trachea using the syringe and check the cuff at regular intervals.

CAUTION

Refrain from excessive rotating movements and traction. Avoid accidental disconnection or kinking of the inserted tube.

▶ Connect the tube to a ventilation system.

If necessary, use the enclosed extension tubing when using manual ventilation systems.

During mechanical ventilation ensure for proper strain-relief.

▶ Ventilate the patient. Verify correct placement and successful ventilation at regular intervals by a suitable method.

▶ Place the necktape around the patient's neck and connect it firmly to the flange.

▶ Check the tube for secure hold.

▶ Tracheal suctioning is possible through the connector of the extension tubing with a suction catheter \leq 16 Fr.

▶ Prior to extubation deflate the cuff completely with a syringe.

The tube is radiopaque.

SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture.

Expiration date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

CAUTION

- The product may be contaminated with potentially infectious substances of human origin.
- Injury hazard, the scalpel is sharp.

The used or faulty product must be disposed of in accordance with the applicable national and international regulations.

WEIGHT

REF	30-08-447-1
Weight (incl. packaging)	0.21 kg (7.4 oz)
Packaging unit	1

MATERIAL DATA

Scalpel	Blade: Carbon steel; handle: PA, SAN	
Bougie	PE	
Tube	Tube: PVC; connector: PSU; flange: PP	
Extension tubing	PP	
Swivel connector	Connector: PC; suction cap: S; o-ring: EPDM	
Syringe	PP, IR	
Necktape	Foam	Inside: PU; outside: PA
	Hook and loop fastener	PES, PP

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