
Intended to be left blank.

Intended to be left blank.

## ScalpelCric®

---

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## English

### INTENDED USE

Surgical cricothyrotomy set for securing the airway in case of upper airway obstruction or as ultima ratio if all other attempts to ventilate the patient have failed. Emergency airway access via the cricothyroid membrane, using scalpel, bougie and tube.

Clinical benefit: Ventilation of a patient

Patient target group: Adult

Place of use: Hospital and prehospital including military use

### INDICATION

- Life-threatening dyspnoea that cannot be controlled in any other way.

No other indications are known.

### CONTRAINDICATION

- Not known.

### SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- The product must only be used by a physician or by authorized and trained medical personnel with adequate skills in the use of the product.
- Only use the product for the patient target group.
- Only use the products for patients for whom a 6.0 mm tube is suitable for cricothyrotomy.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before use, the product must be inspected visually and a functional check must be performed (see chapters "Visual inspection and functional check"). A faulty product must be disposed of.
- The product serves to temporarily secure the airway. If possible it must be replaced by a permanent airway. Maximum recommended period of application: 4 hours.
- Tube and extension tubing are suitable for conducting oxygen and air.
- The enclosed extension tubing is not suitable for use with mechanical ventilation systems.
- Avoid damage to the cuff of the tube by sharp-edged objects.
- Do not kink the inflation line of the tube.
- It is not permitted to make any changes to the product..
- The product is MRI unsafe. The tube can be used under certain conditions (see chapter "MRI note for tube in situ").
- The product is intended for single use and must not be reused and/or reprocessed. The function of the product is impaired by processing. Any reuse entails a potential infection hazard.
- The product is sterile (ethylene oxide).
- The product must not be used if the packaging is damaged or the expiration date has expired.

## VISUAL INSPECTION AND FUNCTIONAL CHECK

- ▶ Check the product for damage and loose particles.
- ▶ If required, shape the bougie to the desired form.
- ▶ Inflate the cuff with the syringe and check for leakage.
- ▶ Remove all air from the cuff with the syringe.

A faulty product must be disposed of (see "Disposal").

## USE

- ▶ Hyperextend the patient's head.
- ▶ **1** Locate the cricothyroid membrane by the palpable depression between the thyroid cartilage and the cricoid cartilage and stabilise this site with index finger and thumb.

If the cricothyroid membrane is **palpable**, proceed as described in **2a**. If the cricothyroid membrane is **impalpable**, proceed as described in **2b**.

### CAUTION

- Applying too much force during the application or too deep insertion of the scalpel can cause injury to the posterior tracheal wall.

- ▶ **2a** Make a horizontal stab incision through the skin and cricothyroid membrane using the scalpel. Proceed as described in **3**.
- ▶ **2b** Make a midline vertical skin incision using the scalpel.
- ▶ Separate tissues with fingers of both hands.
- ▶ Identify and stabilize the cricothyroid membrane again.
- ▶ Make a horizontal stab incision through the cricothyroid membrane.
- ▶ **3** Keep the scalpel perpendicular to the skin and turn the sharp edge 90° caudally.
- ▶ **4** Maintain gentle traction, pulling the scalpel towards you in order to open the incision and to keep it open for the insertion of the bougie.
- ▶ Slide the bougie down the scalpel blade into the trachea.
- ▶ **5** Rotate the bougie caudally and advance gently up to 10 - 15 cm into the trachea (according to the centimetre marks).
- ▶ Remove the scalpel.
- ▶ **6** If necessary, lubricate the tube and ensure that the cuff is completely deflated.
- ▶ Hold the bougie and railroad the tube gently over the bougie into the trachea, until the flange is flush with the patient's neck.
- ▶ Remove the bougie. Make sure that the position of the tube is not changed.


### CAUTION

- The tube must be held in position as long as it is not sufficiently fixed.

- ▶ Inflate the cuff to the minimum pressure required to seal the trachea using the syringe. Check the cuff at regular intervals and readjust if necessary. The cuff pressure may change during repositioning or transportation of the patient.
- ▶ Connect the tube to a ventilation system.  
If necessary, use the enclosed extension tubing when using manual ventilation systems.  
During mechanical ventilation ensure for proper strain relief.
- ▶ Ventilate the patient. Verify correct placement and successful ventilation at regular intervals by a suitable method.

- ▶ Place the necktape around the patient's neck and connect it firmly to the flange.
- ▶ Check the tube for secure hold.
- ▶ Tracheal suctioning is possible through the connector of the extension tubing with a suction catheter ≤ 16 Fr.
- ▶ Prior to removal of the tube deflate the cuff completely with a syringe.

## MRI note for tube in situ:

- ▶ Fix the pilot balloon  to the patient using adhesive tape.
- During MRI examination, the pilot balloon can cause artefacts around the area where it is positioned.
- An MR field strength of up to 3 T is permitted.

## SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture.

Expiration date: See product label

## STORAGE AND TRANSPORT CONDITIONS

### CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

## DISPOSAL

### CAUTION

- The product may be contaminated with potentially infectious substances of human origin.
- Injury hazard, the scalpel is sharp.

The used or faulty product must be disposed of in accordance with the applicable national and international regulations.

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- ▶ Kiinnitä ohjainballonki © kiinnitysteipillä potilaaseen.