
Silicone Resuscitator Bag

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English

INTENDED USE / INDICATION

Any emergency with limited respiratory function or signs of hypoxia.
The silicone resuscitator bag is intended to be used for the following patients:

Patient	Body weight
Infants	< 5 kg
Children	5 - 20 kg
Adults	> 20 kg

Place of use: No restrictions

CONTRAINDICATION

In patients who do not have an empty stomach and with increased risk of aspiration, mask ventilation must not be employed except in emergency situations.

⚠ SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- Intended for use by qualified personnel.
- A functional check must be performed each time the product is used or after every reprocessing (see section "Functional check"). Incorrect assembly of the valves (1 and 2 (see section "Assembly")) may have an impact on performance.
- Do not use the product in a toxic atmosphere.
- Do not use oil, grease or any hydrocarbon-based substance on the product so as to avoid possible flammability during use.
- Do not disassemble the pressure relief valve (6).
- The product must be cleaned, disinfected and sterilised prior to reuse (see section "Reprocessing").
- Reservoir bag and oxygen tubing are for single use and must not be reprocessed.
- Damaged products must not be used.
- Use original VBM parts only. Using parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.

FUNCTIONAL CHECK

Prior to every use, the valve functions must be checked as follows to ensure proper functionality of the resuscitator bag:

A reservoir bag is needed to complete the test procedures described below:

1.0 Intake valve (2):

- a) ▶ Firstly compress the resuscitator bag with one hand and then close the patient valve (1) of the resuscitator bag with the other hand. Release your grip on the resuscitator bag.

Rapid re-expansion of the resuscitator bag confirms efficient air intake via the intake valve (2).

- b) ▶ Close the patient valve (1) of the resuscitator bag and try to compress the resuscitator bag.

If the resuscitator bag cannot be compressed with reasonable force, or if resuscitator bag compression forces the air out between your hand and the neck of the resuscitator bag, the intake valve (2) efficiently prevents backward leakage of air.

2.1 Patient valve (pressure relief valve in closed position) (1):

- a) ▶ Hold a reservoir bag over the patient valve (1) and use your thumb to press on the reservoir bag connector.

- ▶ Ensure a tight seal between the patient valve (1) and reservoir bag.
- ▶ Compress the resuscitator bag with your other hand several times. Ensure that the lip valve (4) opens during compression.

Filling of the reservoir bag confirms that the patient valve (1) efficiently directs air to the patient.

- b) ▶ Compress the filled reservoir bag attached to the patient valve (1) and inspect the flutter valve (5).

Rising of the flutter valve (5) confirms that air is properly released into the atmosphere, rather than returning into the resuscitator bag.

2.2 Patient valve (pressure relief valve in opened position)

- ▶ Close the patient valve (1) with your thumb and compress the resuscitator bag several times. The resuscitator bag can be compressed and the visible and audible opening of the pressure relief valve (6) confirms its proper operation.

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3.0 Reservoir valve (in the intake valve) (2):

- a) ▶ Proceed as described above in 2.1a) in order to fill the reservoir bag with ambient air.
- ▶ Attach the filled reservoir bag to the intake valve (2) and press on the reservoir bag.
- The visible rise of the diaphragm (7) when compressing the reservoir bag confirms that the reservoir bag efficiently vents excessive gas to the atmosphere.
- b) ▶ Proceed as described above in 2.1a) in order to fill a reservoir bag with ambient air.
- ▶ Attach the filled reservoir bag to the intake valve (2).
 - ▶ Perform several compression-release cycles on the resuscitator bag until the reservoir bag is flat and empty. Rapid re-expansion of the resuscitator bag confirms efficient air intake via the diaphragm (8).

USE

- ▶ Optional: Connect reservoir bag and oxygen tubing to the resuscitator bag.
- ▶ Connect the mask to the patient valve.
- ▶ Connect the oxygen tubing to a regulated oxygen source.
- ▶ Adjust the oxygen flow. The reservoir bag expands completely during inflow and nearly collapses as the compressed resuscitator bag refills during exhalation.
- ▶ Prior to connecting to a patient, check the function of the resuscitator bag and make sure that all connections are in proper position. Observe the intake valve, reservoir bag and patient valve in all phases of ventilation; leakages must not occur.
- ▶ Press the mask tightly to the patient's face so that it forms an airtight seal. If the mask does not fit properly, the ventilation may be affected.
- ▶ Fix the mask by hand.
- ▶ Compress the resuscitator bag to deliver a breath. Observe the patient's chest wall rise to confirm inspiration.
- ▶ Release the resuscitator bag to allow patient exhalation. Observe the chest wall fall to confirm exhalation.
- ▶ If there is any contamination in the resuscitator bag during operation, clean the contamination immediately.

REPROCESSING (CLEANING, DISINFECTION, STERILISATION)

General information

- The reusable products are supplied non-sterile and must be unpacked, cleaned, disinfected and sterilised prior to first use.
- Products must be reprocessed prior to each use.
- The manufacturer is not responsible for damage resulting from improper reprocessing.
- The user is obliged to validate their procedures or the devices and accessories and to comply with the validated parameters during every reprocessing.
- It is recommended to use an automated procedure because automated cleaning and disinfection are significantly more effective.
- The effectiveness has been verified by an independent and accredited test laboratory.
- To achieve effective reprocessing, heavy soiling must not be allowed to dry on the product and must be removed immediately after use.
- Silicone products must not come into contact with oil and grease.

Cleaning / disinfection

Automated cleaning / disinfection

Notes regarding the use of washer disinfectors

- Use de-ionised water.

Notes regarding the use of cleaning agents and disinfectants

- When alkaline cleaning agents are used, the agent must be subsequently neutralised.
- Do not use drying agents.

Procedure validated by the manufacturer:

1. Disassemble the resuscitator bag and its accessories into the parts according to the diagram ("Assembly").
2. Position the products in the injector carriage:

Product/accessories	Positioning in the injector carriage
Resuscitator bag	Injector jet
Mask	
Patient valve	Mesh tray (disassemble products into parts; secure with a cover net)
Intake valve	

CAUTION

Ensure effective rinsing of all hollow spaces.

3. Start programme with the following parameters:
 - a. Pre-rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 1 min.
 - b. Cleaning at 55 °C (131 °F) with an exposure time of 5 min using de-ionised water and "Sekumatic® ProClean" cleaning agent (dose: 0.5% (5 ml/d)).
 - c. Neutralisation with "Sekumatic® FNZ" at 20 °C (68 °F) with an exposure time of 2 min (dose: 0.1% (1 ml/d)).
 - d. Rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 2 min.
 - e. Thermal disinfection at 93 °C (199.4 °F) for 5 min with de-ionised water.
 - f. Dry at 100 °C (212 °F).
4. Check for visible contamination. If necessary, repeat the reprocessing.
5. All parts must be dry and free of any damage (cracks, breakages, etc.) before reassembling. Reassemble the resuscitator bag and its accessories according to the diagram ("Assembly").

CAUTION

There must be only "ONE" lip valve in the patient valve.

6. Check the products according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

Manual cleaning / disinfection

Notes regarding the use of cleaning agents and disinfectants

- Prepare the cleaning and disinfection solutions prior to each manual cycle.
- Use de-ionised water.

Procedure validated by the manufacturer:

1. Prepare a 2 % (30 ml/d) cleaning and disinfection solution of "Sekusept® Aktiv" with de-ionised water at 20 °C (68 °F). After 15 min the cleaning and disinfection solution can be used.
2. Disassemble the resuscitator bag and its accessories into the parts according to the diagram ("Assembly").
3. Clean the products with soft sponges in the cleaning and disinfection solution. Areas difficult to access must be cleaned with soft brushes. Hard brushes and other materials which may damage the surface of the products must not be used.
4. Place the products into the cleaning and disinfection solution and leave for 15 min. Ensure all hollow spaces are covered with solution.
5. Remove cleaning and disinfection solution by rinsing thoroughly with de-ionised water. Remaining residues can shorten the product shelf life or can lead to material damage.
6. Dry the products. Avoid any accumulation of water.
7. Check for visible contamination. If necessary, repeat the reprocessing.
8. All parts must be dry and free of any damage (cracks, breakages, etc.) before reassembling. Reassemble the resuscitator bag and its accessories according to the diagram ("Assembly").

CAUTION

There must be only "ONE" lip valve in the patient valve.

9. Check the products according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

Checks

After disinfection and before every use, the products must be inspected for damage (cracks, breakages, etc.). Damaged products must not be used.

Packaging

- Connection pieces must not be sterilised while connected to other connection pieces/products.

After the cleaning and disinfection procedure, the products must be packed into a package which is suitable for steam sterilisation.

Sterilisation

Notes regarding sterilisation:

- Products or sterilisation packages must be protected from mechanical damage.

Sterilisation validated by the manufacturer

- The products must be steam sterilised by a fractional vacuum process.
- Exposure time is 5 min with a sterilisation temperature of 134 °C (273.2 °F).

Reusability

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use. Considering the maximum shelf life of 5 years, reusable products can be reprocessed up to 100 times if used properly and if the cleaning, disinfection and sterilisation instructions are complied with. Every further reuse is the responsibility of the user (see section "Checks").

WARNING

By using the products on patients with suspected prion disease there may be a high risk of transmission. In such case the physician decides if the product is either disposed of (see section "Disposal") or reprocessed according to national directives.

A list with compatible cleaning agents is available at www.vbm-medical.de in the download section.

SHELF LIFE

The shelf life of the product is 5 years and it can be reprocessed 100 times within the shelf life.

STORAGE AND TRANSPORT CONDITIONS

Operating temperature: -18 °C ... +50 °C (-0.4 °F ... 122 °F)

Storage temperature: -40 °C ... +60 °C (-40 °F ... 140 °F)

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must have run through the complete reprocessing process to exclude any risk for the staff of the manufacturer. The manufacturer reserves the right to refuse soiled or contaminated products for reasons of safety.

MAINTENANCE

The medical device may only be serviced by persons who have the technical knowledge, qualification and the required means to properly carry out the service. After servicing the medical device, the structural and functional characteristics essential for safety and functionality must be checked, insofar as they can be influenced by the service procedure.

DISPOSAL

Used or damaged products must be disposed of in accordance with the applicable national and international regulations.

TECHNICAL DATA

	Infants	Children	Adults
Inspiratory resistance	≤ 5 cmH ₂ O @ 50 l/min		
Expiratory resistance	≤ 5 cmH ₂ O @ 50 l/min		
Pressure relief valve	Gas release when pressure ≥ 40 cmH ₂ O		
Dead space	< 6 ml		
Expected delivery volume	up to 125 ml	up to 250 ml	up to 675 ml
Product dimension assembled (L x W x D)	460 x 140 x 120 mm	550 x 140 x 145 mm	630 x 140 x 150 mm
Mass	< 500 grams		

Conversion of pressure units: 1 hPa = 1.01973 cmH₂O = 0.75006 mmHg

SCOPE OF DELIVERY

	Infants	Children	Adults
Resuscitator bag	280 ml	450 ml	1500 ml
Mask	# 1	# 3	# 5
Reservoir bag	600 ml	2500 ml	2500 ml
Oxygen tubing	200 cm		

MATERIAL DATA

Item	Material
Patient valve	PSF (polybisphenylsulfone), silicone
Resuscitator bag	Silicone
Intake valve	PSF (polybisphenylsulfone), silicone
Reservoir bag	PVC (polyvinyl chloride)
Mask	Silicone
Oxygen tubing	PVC (polyvinyl chloride)

Oxygen flow rate	Percent oxygen delivered		
	Device under test Infants	Device under test Children	Device under test Adults
2 l/min	92 %	69 %	50 %
4 l/min	96 %	95 %	74 %
6 l/min	98 %	98 %	96 %
8 l/min	99 %	99 %	98 %
10 l/min	99 %	99 %	99 %
15 l/min	99 %	99 %	99 %
All with assembled oxygen reservoir bag.			

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