

Silicone Resuscitator Bag

reusable

ASSEMBLY	3	LT Naudojimo instrukcija.....	38 - 40
DE Gebrauchsanweisung	4 - 6	LV Lietošanas instrukcija	41 - 43
EN Instructions for use.....	7 - 9	NL Gebruiksaanwijzing	43 - 45
BG Инструкции за употреба	10 - 12	NO Bruksanvisning.....	46 - 48
CS Návod k použití.....	13 - 15	PL Instrukcja użycia.....	48 - 50
DA Brugsanvisning.....	15 - 17	PT Instruções de utilização	51 - 53
EL Οδηγίες χρήσης.....	18 - 20	RO Instrucțiuni de utilizare	53 - 55
ES Instrucciones de uso	21 - 23	RU Инструкция по применению.....	56 - 58
ET Kasutamishuhtis	23 - 25	SK Návod na použitie	59 - 61
FI Käyttöohje	26 - 28	SL Navodila za uporabo	61 - 63
FR Mode d'emploi.....	28 - 30	SV Bruksanvisning.....	64 - 66
HR Upute za uporabu	31 - 33	TR Kullanma kılavuzu	66 - 68
HU Használati utasítás	33 - 35	SYMBOL DESCRIPTION	69 - 71
IT Istruzioni per l'uso.....	36 - 38		

ASSEMBLY

Deutsch

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Englisch

INTENDED USE

For manual ventilation of patients.

Clinical benefit: Ventilation and oxygen supply of patients:

Patient	Body weight
Infant	< 5 kg
Child	5 - 20 kg
Adult	> 20 kg

Place of use: Hospital and prehospital

INDICATIONS

Any emergency with limited respiratory function or signs of hypoxia.

No other indications are known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel.

- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.

- Each time the product is used, visually inspect the product for damage (cracks, breakage, etc.) and carry out a functional check (see chapter "Functional check"). A faulty product must be disposed of.
- Incorrect assembly of the valves (1 and 2 (see section "Assembly")) may have an impact on performance.

- Do not disassemble the pressure relief valve (6).

- The reusable product must be processed completely prior to every use (see chapter "processing").

- Reservoir bag and oxygen tubing are for single use and must not be processed. The function of these products is impaired by processing. Reuse carries the potential risk of infection.

- Use original VBM parts only. Using parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.

- The product is MR unsafe.

FUNCTIONAL CHECK

Prior to every use, the valve functions must be checked as follows to ensure proper functionality of the resuscitator bag:

A reservoir bag is needed to complete the test procedures described below:

1.0 Intake valve (2):

- Firstly compress the resuscitator bag with one hand and then close the Non-Rebreathing Valve (1) of the resuscitator bag with the other hand. Release your grip on the resuscitator bag.

Rapid re-expansion of the resuscitator bag confirms efficient air intake via the intake valve (2).

- Close the Non-Rebreathing Valve (1) of the resuscitator bag and try to compress the resuscitator bag.

If the resuscitator bag cannot be compressed with reasonable force, or if resuscitator bag compression forces the air out between your hand and the neck of the resuscitator bag, the intake valve (2) efficiently prevents backward leakage of air.

2.1 Non-Rebreathing Valve (pressure relief valve in closed position) (1):

- Hold a reservoir bag over the Non-Rebreathing Valve (1) and use your thumb to press on the reservoir bag connector.
 - Ensure a tight seal between the Non-Rebreathing Valve (1) and reservoir bag.
- Compress the resuscitator bag with your other hand several times. Ensure that the lip valve (4) opens during compression.

Filling of the reservoir bag confirms that the Non-Rebreathing Valve (1) efficiently directs air to the patient.

- Compress the filled reservoir bag attached to the Non-Rebreathing Valve (1) and inspect the flutter valve (5).

Rising of the flutter valve (5) confirms that air is properly released into the atmosphere, rather than returning into the resuscitator bag.

2.2 Non-Rebreathing Valve (pressure relief valve in opened position)

- Close the Non-Rebreathing Valve (1) with your thumb and compress the resuscitator bag several times. The resuscitator bag can be compressed and the visible and audible opening of the pressure relief valve (6) confirms its proper operation.

3.0 Reservoir valve (in the intake valve) (2):

- Proceed as described above in 2.1a) in order to fill the reservoir bag with ambient air.

- Attach the filled reservoir bag to the intake valve (2) and press on the reservoir bag.

The visible rise of the diaphragm (7) when compressing the reservoir bag confirms that the reservoir bag efficiently vents excessive gas to the atmosphere.

- Proceed as described above in 2.1a) in order to fill a reservoir bag with ambient air.

- Attach the filled reservoir bag to the intake valve (2).

- Perform several compression-release cycles on the resuscitator bag until the reservoir bag is flat and empty. Rapid re-expansion of the resuscitator bag confirms efficient air intake via the diaphragm (8).

USE

- Optional: Connect reservoir bag and oxygen tubing to the resuscitator bag.

- Connect the mask to the Non-Rebreathing Valve.

- Connect the oxygen tubing to a regulated oxygen source.

- Adjust the oxygen flow. The reservoir bag expands completely during inflow and nearly collapses as the compressed resuscitator bag refills during exhalation.

- Prior to connecting to a patient, check the function of the resuscitator bag and make sure that all connections are in proper position. Observe the intake valve, reservoir bag and Non-Rebreathing Valve in all phases of ventilation; leakages must not occur.

- Apply the Mask firmly over the patient's mouth and nose so that it forms an airtight seal of the upper airway. If the mask does not fit properly, the ventilation may be affected.

- Hold the Mask in position during ventilation.

- Compress the resuscitator bag to deliver a breath. Observe the patient's chest wall rise to confirm inspiration.

- Release the resuscitator bag to allow patient exhalation. Observe the chest wall fall to confirm exhalation.

- If there is any contamination in the resuscitator bag during operation, clean the contamination immediately.

PROCESSING (CLEANING / DISINFECTION, STERILISATION)

General Information

- The reusable product is supplied non-sterile and must be processed completely prior to every use.

- The manufacturer is not responsible for damage resulting from improper processing.

- The user is obliged to validate their procedures or the devices and accessories and to comply with the validated parameters during every processing.

- The effectiveness has been verified by an independent test laboratory.

- To achieve effective processing, heavy soiling must not be allowed to dry on the product and must be removed immediately after use.

- Silicone products must not come into contact with oil and grease.

CLEANING / DISINFECTION

Automated cleaning / disinfection

Notes

- Use de-ionised water.
- Use a washer disinfectant which complies with ISO 15883.
- A Miele PG8581 with Cabinet Miele A 105 was used for validation.
- Do not use drying agents.

Procedure validated by the manufacturer

- Disassemble the breathing bag and accessories into parts as shown in the diagram ("Assembly").
- Position the product in the injector carriage:

Product/accessories	Positioning in the injector carriage
Resuscitator bag Mask	Connected to an injector jet
Non-Rebreathing Valve Intake valve	In a mesh tray (disassemble products into parts; secure with a cover net)

CAUTION

- Ensure effective rinsing of all hollow spaces by using a suitable injector jet.

3. Start program with the following parameters:

- Pre-rinse with water at 20 °C (68 °F) and an exposure time of 1 min.
 - Cleaning at 55 °C (131 °F) with an exposure time of 5 min using de-ionised water and a cleaning agent based on alkalinity donors, surfactants and enzymes. The cleaning agent Neodisher® MediClean Forte (Dr. Weigert GmbH; dose: 0.5 % (5 ml/l)) was used for validation.
 - Rinse with de-ionised water (max. 100 CFU/ml) at 20 °C (68 °F) and an exposure time of 1 min.
 - Thermal disinfection at 90 °C (194 °F) - 95 °C (203 °F) for 5 min with de-ionised water.
 - Dry at 100 °C (212 °F) for 10 min.
4. Check for visible contamination. If necessary, repeat processing.
5. All parts must be dry and free of any damage (cracks, breakages, etc.) before reassembling. Reassemble the resuscitator bag and its accessories according to the diagram ("Assembly").

CAUTION

There must be only one lip valve in the Non-Rebreathing Valve.

6. Check the product according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

CHECKS

- After disinfection, the product must be inspected for damage (cracks, breakage, etc.).

A faulty product must be disposed of (see chapter "Disposal").

PACKAGING

- Connection pieces must not be sterilised while connected to other connection pieces/products.
- After the cleaning and disinfection procedure, the product must be packed into sterile barrier systems suitable for steam sterilisation. The sterile barrier systems must comply with ISO 11607-1. Sterilising See-Through Heat Sealable Rolls (K953776) were used for validation.

STERILISATION

Notes

- Product and sterilisation packaging must be protected from mechanical damage.
- Use a steriliser which complies with the requirements of DIN EN 285 or DIN EN 13060.

Sterilisation validated by the manufacturer

- Steam sterilise the product using a fractionated vacuum process.
- Exposure time is 5 min with a sterilisation temperature of 134 °C (273 °F).
- The drying time is 10 min.
- A Tuttnauer ELARA 11-D (K090783) autoclave was used for validation.

SHELF LIFE

WARNING

By using the products on patients with suspected prion disease there may be a high risk of transmission. In such case the physician decides if the product is either disposed of (see chapter "Disposal") or processed according to national directives.

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use.

The shelf life of the product documented by the manufacturer is 5 years from the date of manufacture. The product can be processed up to 100 times within the shelf life. Every further use of the product is the responsibility of the user and the product must be inspected beforehand in accordance with chapter "Checks".

Use-by date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must have run through the complete processing process to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

MAINTENANCE

The medical device may only be serviced by persons who have the technical knowledge, qualification and the required means to properly carry out the service. After servicing the medical device, the structural and functional characteristics essential for safety and functionality must be checked, insofar as they can be influenced by the service procedure.

In order for the manufacturer to carry out repairs quickly, the medical device must be returned with a description of the fault that is as precise as possible.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

CAUTION

The product may be contaminated with potentially infectious substances of human origin.

PRODUCT SPECIFICATIONS

REF	Patient	Resuscitator bag	Mask	Reservoir bag	Oxygen tubing	Packaging unit
88-10-100	Adult	1500 ml	-	-	-	1
88-10-190			#5	2500 ml	200 cm	
88-10-200	Child	450 ml	-	-	-	
88-10-290			#3	2500 ml	200 cm	
88-10-300	Infant	280 ml	-	-	-	
88-10-390			#1	600 ml	200 cm	

TECHNICAL DATA

Patient	Infant	Child	Adult
Inspiratory resistance	≤ 5 cmH ₂ O @ 60 l/min		
Expiratory resistance	≤ 5 cmH ₂ O @ 60 l/min		
Pressure relief valve	40 cmH ₂ O		
Non-Rebreathing Valve	22 mm O.D. / 15 mm I.D.		
Dead space	≤ 5 ml + 10% of the supplied volume	≤ 5 ml + 10% of the supplied volume	≤ 5 ml + 10% of the supplied volume
Delivered volume	≥ 20 ml	≥ 150 ml	≥ 600 ml

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