
Paediatric Mask

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Einsteinstrasse 1 | 72172 Sulz a.N. | Germany | Tel.: +49 7454 9596-0
Fax: +49 7454 9596-33 | e-mail: info@vbm-medical.de | www.vbm-medical.de

DE	Gebrauchsanweisung	4
EN	Instructions for use	8
BG	Инструкции за употреба	12
CS	Návod k použití	16
DA	Brugsanvisning	19
EL	Οδηγίες χρήσης	23
ES	Instrucciones de uso	27
ET	Kasutusjuhend	31
FI	Käyttöohje	34
FR	Manuel d'utilisation	38
HR	Upute za uporabu	42
HU	Használati utasítás	46
IT	Istruzioni per l'uso	49
LT	Naudojimo instrukcija	53
LV	Lietošanas instrukcija	57
NL	Gebruiksaanwijzing	60
NO	Bruksanvisning	64
PL	Instrukcja użycia	68
PT	Instruções de utilização	71
RO	Instrucţiuni de utilizare	75
RU	Инструкция по применению	79
SK	Návod na použitie	83
SL	Navodila za uporabo	87
SV	Bruksanvisning	90
TR	Kullanım kılavuzu	94
	Symbol Description	98

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EN Instructions for use

INTENDED USE

To seal mouth and nose during ventilation/oxygenation of a patient.

Clinical benefit: Adequate ventilation/oxygenation of a patient by sealing the upper airway.

Patient target group: Newborn, infant and child

Place of use: Hospital and prehospital

INDICATIONS

- Temporary manual or mechanical ventilation/oxygenation.

No other indications are known.

CONTRAINDICATIONS

In non-fasted patients with increased risk of aspiration, mask ventilation must not be employed except in emergency situations.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Each time the product is used, visually inspect the product for damage (cracks, breakage, etc. (see section "Checks"). A faulty product must be disposed of.
- Insufficient ventilation can lead to permanent damage in patients.
- Maximum application time: ≤ 2 hours
- If ventilation/oxygenation is not possible with this product, another adequate method of ventilation must be initiated.
- It is not permitted to make any changes to the product.
- The product is suitable for conducting oxygen, air, anaesthetic gases and their gas mixtures.
- The reusable product is supplied non-sterile and must be cleaned, disinfected and sterilised prior to every use (see section "Reprocessing").
- The product is MR safe.

PREPARATION PRIOR TO USE

- ▶ Select the suitable mask size (see table „Product specifications“).

APPLICATION

WARNING

- In the event of injuries or critical deformities of the face, it is at the discretion of the user to what extent mask ventilation is possible.
- The user must ensure efficient ventilation/oxygenation during application.

- ▶ Connect the mask to a resuscitator bag or ventilation system.
- ▶ Extend the head to keep the airways open

CAUTION

In case of cervical spine injuries, consider whether or not the head should be extended.

- ▶ Apply the mask firmly over the patient's mouth and nose so that it forms an airtight seal of the upper airway. If the mask does not fit properly, the ventilation may be affected.
- ▶ Keep the mask in position during ventilation.

FOLLOWING APPLICATION

- ▶ Remove and disconnect the mask.

PROCESSING (CLEANING, DISINFECTION, STERILISATION)

General information

- The reusable product is supplied non-sterile and must be processed completely prior to every use.
- The manufacturer is not responsible for damage resulting from improper processing.
- The user is obliged to validate their procedures or the devices and accessories and to comply with the validated parameters during every processing.
- The effectiveness has been verified by an independent test laboratory.
- To achieve effective processing, heavy soiling must not be allowed to dry on the product and must be removed immediately after use.
- Silicone products must not come into contact with oil and grease.

CLEANING / DISINFECTION

Automated cleaning / disinfection

Notes

- Use de-ionised water.
- Use a washer disinfectant which complies with ISO 15883.
- A Miele PG8581 with Cabinet Miele A 105 was used for validation.
- Do not use drying agents.

Procedure validated by the manufacturer

1. Position the product in the injector carriage:

Product / Parts	Position in injector carriage
Hook ring	In a mesh tray (secure with a cover net)
Mask	Connected to an injector jet

CAUTION

- Ensure effective rinsing of all hollow spaces by using a suitable injector jet.

2. Start program with the following parameters:
 - a. Pre-rinse with water at 20 °C (68 °F) and an exposure time of 1 min.
 - b. Cleaning at 55 °C (131 °F) with an exposure time of 5 min using de-ionised water and a cleaning agent based on alkalinity donors, surfactants and enzymes. The cleaning agent Neodisher® MediClean Forte (Dr. Weigert GmbH; dose: 0.5 % (5 ml/l)) was used for validation.
 - c. Rinse with de-ionised water (max. 100 CFU/ml) at 20 °C (68 °F) and an exposure time of 1 min.
 - d. Thermal disinfection at 90 °C (194 °F) - 95 °C (203 °F) for 5 min with de-ionised water.
 - e. Dry at 100 °C (212 °F) for 10 min.
3. Check for visible contamination. If necessary, repeat processing.
4. Check the product according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

CHECKS

- ▶ After disinfection, the product must be inspected for damage (cracks, breakage, etc.). A faulty product must be disposed of (see chapter "Disposal").

PACKAGING

- ▶ After the cleaning and disinfection procedure, the product must be packed into a sterile barrier system suitable for steam sterilisation. The sterile barrier systems must comply with ISO 11607-1. Striking See-Through Heat Sealable Rolls (K953776) were used for validation.

STERILISATION

Notes

- Product and sterilisation packaging must be protected from mechanical damage.
- Use a steriliser which complies with the requirements of DIN EN 285 or DIN EN 13060.

Sterilisation validated by the manufacturer

- Steam sterilise the product using a fractionated vacuum process.
- Exposure time is 5 min with a sterilisation temperature of 134 °C (273 °F).
- The drying time is 10 min.
- A Tuttnauer ELARA 11-D (K090783) autoclave was used for validation.

SHELF LIFE

WARNING

By using the products on patients with suspected prion disease there may be a high risk of transmission. In such case the physician decides if the product is either disposed of (see chapter "Disposal") or processed according to national directives.

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use.

The shelf life of the product documented by the manufacturer is 5 years from the date of manufacture. The product can be processed up to 100 times within the shelf life.

Every further use of the product is the responsibility of the user and the product must be inspected beforehand in accordance with chapter "Checks".

Use-by date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must have run through the complete processing process to exclude any risk for the manufacturer's staff and third parties. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

CAUTION

- The product may be contaminated with potentially infectious substances of human origin.

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