## Mask

---

for infant and small child

<b>DE</b>	<b>MASKE</b> Gebrauchsanweisung .....Seite 3 - 4	<b>LT</b>	<b>KAUKĖ</b> Naudojimo instrukcija..... 29 - 30 psl.
<b>EN</b>	<b>MASK</b> Instructions for use.....Page 5 - 6	<b>LV</b>	<b>MASKA</b> Lietošanas instrukcija ..... 31. - 32. lpp.
<b>BG</b>	<b>MACKA</b> Инструкции за употреба .....страница 7 - 8	<b>NL</b>	<b>MASKER</b> Gebruiksaanwijzing .....pagina 33 - 34
<b>CS</b>	<b>MASKA</b> Návod k použití.....Strana 9 - 10	<b>NO</b>	<b>MASKE</b> Bruksanvisning.....Side 35 - 36
<b>DA</b>	<b>MASKE</b> Brugsanvisning.....Side 11 - 12	<b>PL</b>	<b>MASKA</b> Instrukcja użycia..... Strona 37 - 38
<b>EL</b>	<b>ΜΑΣΚΑ</b> Οδηγίες χρήσης.....Σελίδα 13 - 14	<b>PT</b>	<b>MÁSCARA</b> Instruções de utilização .....Páginas 39 - 40
<b>ES</b>	<b>MASCARILLA</b> Instrucciones de uso .....Página 15 - 16	<b>RO</b>	<b>MASCĂ</b> Instrucțiuni de utilizare .....pagina 41 - 42
<b>ET</b>	<b>MASK</b> Kasutamisinjuhised.....Lehekülg 17 - 18	<b>RU</b>	<b>MACKA</b> Инструкция по применению .....Страница 43 - 44
<b>FI</b>	<b>MASKI</b> Käyttöohje .....Sivu 19 - 20	<b>SK</b>	<b>MASKA</b> Návod na použitie .....Strana 45 - 46
<b>FR</b>	<b>MASQUE</b> Mode d'emploi.....Page 21 - 22	<b>SL</b>	<b>MASKA</b> Navodila za uporabo ..... Stran 47 - 48
<b>HR</b>	<b>MASKA</b> Upute za uporabu ..... stranica 23 - 24	<b>SV</b>	<b>MASK</b> Bruksanvisning.....Sid. 49 - 50
<b>HU</b>	<b>MASZK</b> Használati utasítás.....25 - 26. oldal	<b>TR</b>	<b>MASKE</b> Kullanma kılavuzu ..... Sayfa 51 - 52
<b>IT</b>	<b>MASCHERA</b> Istruzioni per l'uso.....Pagina 27 - 28		<b>SYMBOL DESCRIPTION</b> .....Page 53 - 56

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# English

## INTENDED USE / INDICATION

During manual or mechanical ventilation, the mask is used to seal the mouth and nose area.

Place of use: No restrictions

## CONTRAINDICATION

In patients who do not have an empty stomach and with increased risk of aspiration, mask ventilation must not be employed except in emergency situations.

## SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- Intended for use by qualified personnel.
- A check must be performed each time the product is used or after every reprocessing (see section "Checks"). Damaged products must not be used.
- Insufficient ventilation can lead to permanent damage in patients.
- If ventilation is not possible with this product, another adequate type of ventilation must be initiated.
- Do not use the product in a toxic atmosphere.
- Do not use oil, grease or any hydrocarbon-based substance on the product so as to avoid possible flammability during use.
- The product must be cleaned, disinfected and sterilised prior to reuse (see section "Reprocessing").
- Use original VBM parts only. Using parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.

## USE

- ▶ Select the suitable mask size (see table "Technical data").
- ▶ Connect the mask to a resuscitator bag or ventilation system.

### CAUTION

In case of cervical spine injuries, consider whether or not the head should be extended.

- ▶ Extend the head so that the airway remains patent!

- ▶ Press the mask tightly to the patient's face so that it forms an airtight seal. If the mask does not fit properly, the ventilation may be affected.
- ▶ Fix the mask by hand. Option: Fix the mask using a fixation tape and a hook ring.

## REPROCESSING (CLEANING, DISINFECTION, STERILISATION)

### General information

- The reusable products are supplied non-sterile and must be unpacked, cleaned, disinfected and sterilised prior to first use.
- Products must be reprocessed prior to each use.
- The manufacturer is not responsible for damage resulting from improper reprocessing.
- The user is obliged to validate their procedures or the devices and accessories and to comply with the validated parameters during every reprocessing.
- It is recommended to use an automated procedure because automated cleaning and disinfection are significantly more effective.

- The effectiveness has been verified by an independent and accredited test laboratory.
- To achieve effective reprocessing, heavy soiling must not be allowed to dry on the product and must be removed immediately after use.
- Silicone products must not come into contact with oil and grease.

## Cleaning / disinfection

### Automated cleaning / disinfection

#### Notes regarding the use of washer disinfectors

- Use de-ionised water.

#### Notes regarding the use of cleaning agents and disinfectants

- When alkaline cleaning agents are used, the agent must be subsequently neutralised.
- Do not use drying agents.

#### Procedure validated by the manufacturer:

1. Position the products in the injector carriage:

Product / Parts	Position in injector carriage
Hook ring	Mesh tray (secure with a cover net)
Mask	Injector jet

### CAUTION

Ensure effective rinsing of all hollow spaces.

2. Start programme with the following parameters:
  - a. Pre-rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 1 min.
  - b. Cleaning at 55 °C (131 °F) with an exposure time of 5 min using de-ionised water and "Sekumatic® ProClean" cleaning agent (dose: 0.5% (5 ml/l)).
  - c. Neutralisation with "Sekumatic® FNZ" at 20 °C (68 °F) with an exposure time of 2 min (dose: 0.1% (1 ml/l)).
  - d. Rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 2 min.
  - e. Thermal disinfection at 93 °C (199.4 °F) for 5 min with de-ionised water.
  - f. Dry at 100 °C (212 °F).
3. Check for visible contamination. If necessary, repeat the reprocessing.
4. Check the products according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

### Manual cleaning / disinfection

#### Notes regarding the use of cleaning agents and disinfectants

- Prepare the cleaning and disinfection solutions prior to each manual cycle.
- Use de-ionised water.

#### Procedure validated by the manufacturer:

1. Prepare a 2 % (30 ml/l) cleaning and disinfection solution of "Sekusept® Aktiv" with de-ionised water at 20 °C (68 °F). After 15 min the cleaning and disinfection solution can be used.

- Clean the products with soft sponges in the cleaning and disinfection solution. Areas difficult to access must be cleaned with soft brushes. Hard brushes and other materials which may damage the surface of the products must not be used.
- Place the products into the cleaning and disinfection solution and leave for 15 min. Ensure all hollow spaces are covered with solution.
- Remove cleaning and disinfection solution by rinsing thoroughly with de-ionised water. Remaining residues can shorten the product shelf life or can lead to material damage.
- Dry the products. Avoid any accumulation of water.
- Check for visible contamination. If necessary, repeat the reprocessing.
- Check the products according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

### Checks

After disinfection and before every use, the products must be inspected for damage (cracks, breakages, etc.).

Damaged products must not be reused.

### Packaging

After the cleaning and disinfection procedure, the products must be packed into a package which is suitable for steam sterilisation.

### Sterilisation

#### Notes regarding sterilisation:

- Products or sterilisation packages must be protected from mechanical damage.

#### Sterilisation validated by the manufacturer

- The products must be steam sterilised by a fractional vacuum process.
- Exposure time is 5 min with a sterilisation temperature of 134 °C (273.2 °F).

### Reusability

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use. Considering the maximum shelf life of 5 years, reusable products can be reprocessed up to 100 times if used properly and if the cleaning, disinfection and sterilisation instructions are complied with.

Every further reuse is the responsibility of the user (see section "Checks").

### WARNING

By using the products on patients with suspected prion disease there may be a high risk of transmission. In such case the physician decides if the product is either disposed of (see section "Disposal") or reprocessed according to national directives.

A list with compatible cleaning agents is available at [www.vbm-medical.de](http://www.vbm-medical.de) in the download section.

### SHELF LIFE

The shelf life of the product is 5 years and it can be reprocessed 100 times within the shelf life.

### STORAGE AND TRANSPORT CONDITIONS

Operating temperature: -18 °C ... +50 °C (-0.4 °F ... 122 °F)

Storage temperature: -40 °C ... +60 °C (-40 °F ... 140 °F)

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

### SERVICE

Prior to the return of medical devices for complaints/repairs, the products must have run through the complete reprocessing process to exclude any risk for the staff of the manufacturer. The manufacturer reserves the right to refuse soiled or contaminated products for reasons of safety.

### DISPOSAL

The used or damaged products must be disposed of in accordance with the applicable national and international statutory regulations.

### TECHNICAL DATA

REF	88-11-000	88-11-001	88-11-002	88-60-266
Description	Small infant mask	Large infant mask	Small child mask	Hook ring
Size	#0	#1	#2	-
Material	Silicone			PSF (polybisphenylsulfone)
Connections	15 mm O.D.	15 mm O.D.	22 mm I.D.	Ø 22 mm

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