## PEEP-Valve & PEEP-Adapter


---

reusable

<b>ASSEMBLY</b> .....	Page 3
<b>DE</b> PEEP-VENTIL & PEEP-ADAPTER, WIEDERVERWENDBAR Gebrauchsanweisung .....	Seite 4 - 5
<b>EN</b> PEEP-VALVE & PEEP-ADAPTER, REUSABLE Instructions for use.....	Page 6 - 7
<b>BG</b> PEEP ВЕНТИЛ И PEEP АДАПТЕР, ЗА МНОГОКРАТНА УПОТРЕБА Инструкции за употреба .....	страница 8 - 9
<b>CS</b> VENTIL PEEP A ADAPTÉR PEEP, OPAKOVANÉ POUŽITELNÉ Návod k použití.....	Strana 10 - 11
<b>DA</b> PEEP-VENTIL & PEEP-ADAPTER, KAN GENANVENDES Brugsanvisning.....	Side 12 - 13
<b>EL</b> ΒΑΛΒΙΔΑ PEEP ΚΑΙ ΠΡΟΣΑΡΜΟΓΕΑΣ PEEP, ΕΠΑΝΑΧΡΗΣΙΜΟΠΟΙΟΥΜΕΝΑ Οδηγίες χρήσης.....	Σελίδα 14 - 15
<b>ES</b> VÁLVULA PEEP Y ADAPTADOR PEEP, REUTILIZABLES Instrucciones de uso.....	Página 16 - 17
<b>ET</b> PEEP KLAPP JA PEEP ADAPTER, TAASKASUTATAVAD Kasutamisejuhised.....	Lehekülg 18 - 19
<b>FI</b> PEEP-VENTTIILI JA PEEP-SOVITIN, UUDELEENKÄYTETTÄVÄ Käyttöohje .....	Sivu 20 - 21
<b>FR</b> VALVE DE PEP ET ADAPTATEUR DE PEP, RÉUTILISABLE Mode d'emploi.....	Page 22 - 23
<b>HR</b> PEEP VENTIL I PEEP ADAPTER ZA VIŠEKRATNU UPORABU Upute za uporabu .....	stranica 24 - 25
<b>HU</b> PEEP SEZLEP ÉS PEEP ADAPTER, ÚJRAFELHASZNÁLHATÓ Használati utasítás.....	26 - 27. oldal
<b>IT</b> VALVOLA PEEP E ADATTATORE PEEP, RIUTILIZZABILI Istruzioni per l'uso.....	Pagina 28 - 29
<b>LT</b> PEEP VOŽTUVAS IR PEEP ADAPTERIS, DAUGKARTINIO NAUDOJIMO Naudojimo instrukcija.....	30 - 31 psl.
<b>LV</b> PEEP VĀRSTS UN PEEP ADAPTERIS, IZMANTOJAMI ATKĀRTOTI Lietošanas instrukcija .....	32. - 33. lpp.
<b>NL</b> PEEP-VENTIEL & PEEP-ADAPTER, HERBRUIKBAAR Gebruiksaanwijzing.....	pagina 34 - 35
<b>NO</b> PEEP-VENTIL OG PEEP-ADAPTER, GJENBRUKBAR Bruksanvisning.....	Side 36- 37
<b>PL</b> ZAWÓR PEEP I ADAPTER PEEP, DO PONOWNEGO UŻYCIA Instrukcja użycia.....	Strona 38 - 39
<b>PT</b> VÁLVULA PEEP E ADAPTADOR PEEP, REUTILIZÁVEL Instruções de utilização .....	Página 40 - 41
<b>RO</b> VALVĂ PEEP ȘI ADAPTOR PEEP, REUTILIZABIL Instrucțiuni de utilizare .....	pagina 42 - 43
<b>RU</b> PEEP-КЛАПАН И PEEP-АДАПТЕР, МНОГОРАЗОВЫЕ Инструкция по применению .....	Страница 44 - 45
<b>SK</b> VENTIL PEEP A ADAPTÉR PEEP, OPĀTOVNE POUŽITELNÉ Návod na použitie.....	Strana 46 - 47
<b>SL</b> VENTIL PEEP IN ADAPTER PEEP ZA PONOVRNO UPORABO Navodila za uporabo .....	Stran 48 - 49
<b>SV</b> PEEP-VENTIL OCH PEEP-ADAPTER, ÅTERANVÄNDBARA Bruksanvisning.....	Sid. 50 - 51
<b>TR</b> PEEP VALFİ & PEEP ADAPTÖRÜ, YENİDEN KULLANILABİLİR Kullanma kılavuzu .....	Sayfa 52 - 53
<b>SYMBOL DESCRIPTION</b> .....	Page 54 - 56

# ASSEMBLY

## PEEP-VALVE

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# English

## INTENDED USE

The PEEP-Valve is intended to generate a positive end-expiratory pressure.

## INDICATION

Can prevent a collapse of the alveoli and thus prevent atelectasis.

## CONTRAINDICATION

None known.

## SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- Intended for use by qualified personnel.
- A check must be performed each time the product is used or after every reprocessing (see section "Checks"). Damaged products must not be used.
- Only use the PEEP-Valve if the patient's condition can be precisely monitored with blood gas analyses and haemodynamic measurements.
- The use of the PEEP-Valve on resuscitator bags and ventilation systems is only permitted if its use has been explicitly specified by the manufacturer.
- Do not use the product in a toxic atmosphere.
- Do not use oil, grease or any hydrocarbon-based substance on the product so as to avoid possible flammability during use.
- The product must be cleaned, disinfected and sterilised prior to reuse (see section "Reprocessing").
- Use original VBM parts only. Using parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.

## USE

- ▶ Connect the PEEP-Valve (using the PEEP-Adapter if necessary) and set the desired PEEP-level by turning the adjustment knob on the PEEP-Valve.

## REPROCESSING (CLEANING, DISINFECTION, STERILISATION)

### General information

- The reusable products are supplied non-sterile and must be unpacked, cleaned, disinfected and sterilised prior to first use.
- Products must be reprocessed prior to each use.
- The manufacturer is not responsible for damage resulting from improper reprocessing.
- The user is obliged to validate their procedures or the devices and accessories and to comply with the validated parameters during every reprocessing.
- It is recommended to use an automated procedure because automated cleaning and disinfection are significantly more effective.
- The effectiveness has been verified by an independent and accredited test laboratory.
- To achieve effective reprocessing, heavy soiling must not be allowed to dry on the product and must be removed immediately after use.
- Silicone products must not come into contact with oil and grease.

- Prior to automated / manual cleaning and sterilisation, set the PEEP-Valve to the lowest PEEP-level (turn the adjustment knob upwards).

## Cleaning / disinfection

### Automated cleaning / disinfection

#### Notes regarding the use of washer disinfectors

- Use de-ionised water.

#### Notes regarding the use of cleaning agents and disinfectants

- When alkaline cleaning agents are used, the agent must be subsequently neutralised.
- Do not use drying agents.

### Procedure validated by the manufacturer

1. Disassemble products into three parts in accordance with the diagram (see section "Assembly").
2. Position products in suitable mesh trays in the injector carriage. Secure products with a cover net.

### CAUTION

Ensure effective rinsing of all hollow spaces.

3. Start programme with the following parameters:
  - a. Pre-rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 1 min.
  - b. Cleaning at 55 °C (131 °F) with an exposure time of 5 min using de-ionised water and "Sekumatic® ProClean" cleaning agent (dose: 0.5 % (5 ml/l)).
  - c. Neutralisation with "Sekumatic® FNZ" at 20 °C (68 °F) with an exposure time of 2 min (dose: 0.1 % (1 ml/l)).
  - d. Rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 2 min.
  - e. Thermal disinfection at 93 °C (199.4 °F) for 5 min with de-ionised water.
  - f. Dry at 100 °C (212 °F).
4. Check for visible contamination. If necessary, repeat the reprocessing.
5. All parts must be dry and free of any damage (cracks, breakages, etc.) before reassembling. Reassemble products in accordance with the diagram (see section "Assembly").
6. Check products according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

### Manual cleaning / disinfection

#### Notes regarding the use of cleaning agents and disinfectants

- Prepare the cleaning and disinfection solutions prior to each manual cycle.
- Use de-ionised water.

### Procedure validated by the manufacturer

1. Prepare a 2 % (30 ml/l) cleaning and disinfection solution of "Sekusept® Aktiv" with de-ionised water at 20 °C (68 °F). After 15 min the cleaning and disinfection solution can be used.
2. Disassemble products into three parts in accordance with the diagram (see section "Assembly").
3. Clean the products with soft sponges in the cleaning and disinfection solution. Areas difficult to access must be cleaned with soft brushes. Hard brushes and other materials which may damage the surface of the products must not be used.

- Place the products into the cleaning and disinfection solution and leave for 15 min. Ensure all hollow spaces are covered with solution.
- Remove cleaning and disinfection solution by rinsing thoroughly with de-ionised water. Remaining residues can shorten the product shelf life or can lead to material damage.
- Dry the products. Avoid any accumulation of water.
- Check for visible contamination. If necessary, repeat the reprocessing.
- All parts must be dry and free of any damage (cracks, breakages, etc.) before reassembling. Reassemble products in accordance with the diagram (see section "Assembly").
- Check products according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

### Checks

After disinfection and before every use, the products must be inspected for damage (cracks, breakages, etc.) and correct functioning.

Test the PEEP-Valve (2 - 10 cmH<sub>2</sub>O) at 5 and 10 cmH<sub>2</sub>O.

Test the PEEP-Valve (5 - 20 cmH<sub>2</sub>O) at 10 and 20 cmH<sub>2</sub>O.

The tolerance is  $\pm 2$  cmH<sub>2</sub>O, and the values must be within this tolerance for 3 seconds.

Damaged products must not be reused.

### Packaging

- Do not sterilise the assembled PEEP-Valve and PEEP-Adapter with other products.

After the cleaning and disinfection procedure, the products must be packed into a package which is suitable for steam sterilisation.

### Sterilisation

#### Notes regarding sterilisation

- Products or sterilisation packages must be protected from mechanical damage.

#### Sterilisation validated by the manufacturer

- The products must be steam sterilised by a fractional vacuum process.
- Exposure time is 5 min with a sterilisation temperature of 134 °C (273.2 °F).

### Reusability

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use. Considering the maximum shelf life of 5 years, reusable products can be reprocessed up to 100 times if used properly and if the cleaning, disinfection and sterilisation instructions are complied with.

Every further reuse is the responsibility of the user (see section "Checks").

### WARNING

By using the products on patients with suspected prion disease there may be a high risk of transmission. In such case the physician decides if the product is either disposed of (see section "Disposal") or reprocessed according to national directives.

A list with compatible cleaning agents is available at [www.vbm-medical.de](http://www.vbm-medical.de) in the download section.

### SHELF LIFE

The shelf life of the product is 5 years and it can be reprocessed 100 times within the shelf life.

### STORAGE AND TRANSPORT CONDITIONS

Operating temperature: -18 °C ... +50 °C (-0.4 °F ... 122 °F)

Storage temperature: -40 °C ... +60 °C (-40 °F ... 140 °F)

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

### SERVICE

Prior to the return of medical devices for complaints/repairs, the products must have run through the complete reprocessing process to exclude any risk for the staff of the manufacturer. For safety reasons, the manufacturer preserves the right to refuse soiled or contaminated products.

### DISPOSAL

The used or damaged products must be disposed of in accordance with the applicable national and international statutory regulations.

### TECHNICAL DATA

REF	88-13-006	88-13-007	88-13-003
Description	PEEP-Valve		PEEP-Adapter
Colour	Blue		-
Pressure range	2 - 10 cmH <sub>2</sub> O	5 - 20 cmH <sub>2</sub> O	-
Material	PSF (polybisphenylsulfone)		PSF (polybisphenylsulfone) Silicone (sealing ring)
Connections	30 mm I.D. / 22 mm O.D.		22 mm I.D. / 30 mm O.D.

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