

PEEP-Valve & PEEP-Adapter


reusable

ASSEMBLY

PEEP-VALVE

PEEP-ADAPTER

Ø 22 mm I.D. / 30 mm O.D.

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Einsteinstrasse 1 | 72172 Sulz a.N. | Germany | Tel.: +49 7454 9596-0
Fax: +49 7454 9596-33 | e-mail: info@vbm-medical.de | www.vbm-medical.de

DE	Gebrauchsanweisung	5
EN	Instructions for use	9
BG	Инструкции за употреба	13
CS	Návod k použití	18
DA	Brugsanvisning	21
EL	Οδηγίες χρήσης	25
ES	Instrucciones de uso	30
ET	Kasutusjuhend	34
FI	Käyttöohje	38
FR	Manuel d'utilisation	42
HR	Upute za uporabu	46
HU	Használati utasítás	50
IT	Istruzioni per l'uso	54
LT	Naudojimo instrukcija	58
LV	Lietošanas instrukcija	61
NL	Gebruiksaanwijzing	65
NO	Bruksanvisning	69
PL	Instrukcja użycia	73
PT	Instruções de utilização	77
RO	Instrucţiuni de utilizare	81
RU	Инструкция по применению	85
SK	Návod na použitie	90
SL	Navodila za uporabo	94
SV	Bruksanvisning	97
TR	Kullanım kılavuzu	101
	Symbol Description	106

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EN Instructions for use

INTENDED USE

The PEEP-Valve is intended to generate a positive end-expiratory pressure.

Clinical benefit: Improvement of the oxygen supply

Patient target group: No restrictions

Place of use: Hospital and prehospital

INDICATIONS

- Can prevent a collapse of the alveoli and thus prevent atelectasis.

No further indications are known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Each time the product is used, visually inspect the product for damage (cracks, breakage, etc.) and carry out a functional check (see chapter "Checks"). A faulty product must be disposed of.
- Do not use oil, grease or any hydrocarbon-based substance on the product so as to avoid possible flammability during use.
- The reusable product must be processed completely prior to every use (see chapter "processing").
- Use original VBM parts only. Using parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.

- The product is MR unsafe.

- The product must not be used if the shelf life has elapsed.

USE

- ▶ To change the lower housing, disassemble the PEEP-Valve as shown in the diagram and then reassemble it with the other lower housing (see chapter "Assembly").
- ▶ Connect the PEEP-Valve (using the PEEP-Adapter if necessary) and set the desired PEEP-level by turning the adjustment knob on the PEEP-Valve.

For precise adjustment of the PEEP-Valve, a manometer can be connected to the ventilator for monitoring

PROCESSING (CLEANING / DISINFECTION, STERILISATION)

Notes

- The manufacturer is not responsible for damage resulting from improper processing.
- The user is obliged to validate their procedures or the devices and accessories and to comply with the validated parameters during every processing.
- The effectiveness has been verified by an independent test laboratory.
- To achieve effective reprocessing, heavy soiling must not be allowed to dry on the product and must be removed immediately after use.
- Silicone products must not come into contact with oil and grease.
- Prior to cleaning and sterilisation, set the PEEP-Valve to the lowest PEEP-level (turn the adjustment knob upwards).

AUTOMATED CLEANING / DISINFECTION

Notes

- Use de-ionised water.
- Use a washer disinfector which complies with ISO 15883.
- A Miele PG8581 with Cabinet Miele A 105 was used for validation.
- Do not use drying agents.

Procedure validated by the manufacturer

1. Disassemble products into three parts in accordance with the diagram (see section "Assembly").
2. Position products in suitable mesh trays in the injector carriage. Secure products with a cover net.

CAUTION

- Ensure effective rinsing of all hollow spaces by using a suitable injector jet.

3. Start programme with the following parameters:
 - a. Pre-rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 1 min.
 - b. Cleaning at 55 °C (131 °F) with an exposure time of 5 min using de-ionised water and a cleaning agent based on alkalinity donors, surfactants and enzymes.
The cleaning agent Neodisher® MediClean Forte (Dr. Weigert GmbH; dose: 0.5 % (5 ml/l)) was used for validation.
 - c. Rinse with de-ionised water (max. 100 CFU/ml) at 20 °C (68 °F) and an exposure time of 1 min.
 - d. Thermal disinfection at 90 °C (194 °F) - 95 °C (203 °F) for 5 min with de-ionised water.
 - e. Dry at 100 °C (212 °F) for 10 min.
4. Check for visible contamination. If necessary, repeat processing.
5. All parts must be dry and free of any damage (cracks, breakages, etc.) before reassembling. Reassemble products in accordance with the diagram (see section "Assembly").
6. Check products according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

Checks

- ▶ After disinfection, the product must be checked for damage (cracks, breakage, etc.) and correct function using a manometer. The following values must be checked:

88-13-006	88-13-007
5 cmH ₂ O	10 cmH ₂ O
10 cmH ₂ O	20 cmH ₂ O

The tolerance is ± 2 cmH₂O, and the values must be within this tolerance for 3 seconds. A faulty product must be disposed of (see chapter "Disposal").

PACKAGING

- Do not sterilise the assembled PEEP-Valve and PEEP-Adapter with other products. After the cleaning and disinfection procedure, the products must be packed into a package which is suitable for steam sterilisation.
- ▶ After the cleaning and disinfection procedure, the product must be packed into a sterile barrier system suitable for steam sterilisation. The sterile barrier systems must comply with ISO 11607-1. Striking See-Through Heat Sealable Rolls (K953776) were used for validation.

STERILISATION

Notes regarding sterilisation

- Product and sterilisation packaging must be protected from mechanical damage.
- Use a steriliser which complies with the requirements of DIN EN 285 or DIN EN 13060.

Sterilisation validated by the manufacturer

- Steam sterilise the product using a fractionated vacuum process.
- Exposure time is 5 min with a sterilisation temperature of 134 °C (273 °F).
- The drying time is 10 min.
- A Tuttnauer ELARA 11-D (K090783) autoclave was used for validation.

SHELF LIFE

WARNING

By using the products on patients with suspected prion disease there may be a high risk of transmission. In such case the physician decides if the product is either disposed of (see chapter "Disposal") or processed according to national directives.

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use. The shelf life of the product documented by the manufacturer is 5 years from the date of manufacture. The product can be processed up to 100 times within the shelf life. Every further use of the product is the responsibility of the user and the product must be inspected beforehand in accordance with chapter "Checks".

Use-by date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must have run through the complete processing process to exclude any risk for the manufacturer's staff and third parties. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

CAUTION

The product may be contaminated with potentially infectious substances of human origin.

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