


PEEP-Valve

single use

ASSEMBLY

PEEP-VALVE

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EN Instructions for use

INTENDED USE

The PEEP-Valve is intended to generate a positive end-expiratory pressure.

Clinical benefit: Improvement of the oxygen supply

Patient target group: No restrictions

Place of use: Hospital and prehospital

INDICATIONS

- Can prevent a collapse of the alveoli and thus prevent atelectasis.

No further indications are known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.
- This product must only be used by medically trained personnel.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before use, visually inspect the product for damage (cracks, breakage, etc.) and carry out a functional check (see chapter "Functional check"). A faulty product must be disposed of.
- Do not use oil, grease or any hydrocarbon-based substance on the product so as to avoid possible flammability during use.
- The product is intended for single use and must not be reused and/or processed. The function of the product is impaired by processing. Any reuse entails a potential infection hazard.
- Use original VBM parts only. Using parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.
- The product is MR unsafe.
- The product must not be used if the shelf life has elapsed.

VISUAL- AND FUNCTIONAL CHECK

- ▶ The product must be checked for damage and loose particles.
- ▶ If present, remove loose particles.
- ▶ Check correct function by using a manometer. The following values must be checked:

89-13-006	89-13-007
5 cmH ₂ O	10 cmH ₂ O
10 cmH ₂ O	20 cmH ₂ O

The tolerance is ± 2 cmH₂O, and the values must be within this tolerance for 3 seconds.

- ▶ A faulty product must be disposed of (see chapter "Disposal").

USE

- ▶ To change the lower housing, disassemble the PEEP-Valve as shown in the diagram and then reassemble it with the other lower housing (see chapter "Assembly").
- ▶ Connect the PEEP-Valve (using a suitable PEEP-Adapter if necessary) and set the desired PEEP-level by turning the adjustment knob on the PEEP-Valve.

For precise adjustment of the PEEP-Valve, a manometer can be connected to the ventilator for monitoring.

SHELF LIFE

The shelf life of the product is 5 years from the time of manufacture.

Use-by date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

CAUTION

- The product may be contaminated with potentially infectious substances of human origin.

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