
Tracheostomy Necktape

DE	TRACHEOSTOMIE-HALSBAND Gebrauchsanweisung 4 - 5	LT	TRACHEOSTOMINĖ APYKAKLĖ Naudojimo instrukcija 30 - 31
EN	TRACHEOSTOMY NECKTAPE Instructions for use 6 - 7	LV	TRAHEOSTOMIJAS KAKLA APSĒJS Lietošanas instrukcija 32 - 33
BG	ПОДДЪРЖАЩА ЛЕНТА ЗА ТРАХЕОСТОМИЯ Инструкции за употреба 8 - 9	NL	TRACHEOTOMIE-NEKBAND Gebruiksaanwijzing 34 - 35
CS	TRACHEOSTOMICKÁ KRČNÍ PÁSKA Návod k použití 10 - 11	NO	TRAKEOSTOMI-HALSBÅND Bruksanvisning 36 - 37
DA	TRACHEOSTOMIKRAVE Bruganvisning 12 - 13	PL	OPASKA MOCUJĄCA RURKĘ TRACHEOSTOMIJĄ Instrukcja użycia 38 - 39
EL	ΚΟΛΑΡΟ ΤΡΑΧΕΙΟΣΤΟΜΙΑΣ Οδηγίες χρήσης 14 - 15	PT	COLAR DE TRAQUEOSTOMIA Instruções de utilização 40 - 41
ES	COLLARÍN PARA TRAQUEOSTOMÍA Instrucciones de uso 16 - 17	RO	COLIER DE TRAHEOSTOMIE Instrucțiuni de utilizare 42 - 43
ET	TRAHHEOSTOOMIA KAELALINT Kasutamisujuhised 18 - 19	RU	ТРАХЕОСТОМИЧЕСКИЙ ФИКСАТОР Инструкция по применению 44 - 45
FI	TRAHEOSTOMIAKAULANAUHA Käyttöohje 20 - 21	SK	TRACHEOSTOMICKÁ FIXAČNÁ PÁSKA Návod na použitie 46 - 47
FR	COLLIER DE TRACHÉOTOMIE Manuel d'utilisation 22 - 23	SL	TRAHEOSTOMSKI OVRATNI TRAK Navodila za uporabo 48 - 49
HR	VRATNA TRAKA ZA TRAHEOTOMIJU Upute za uporabu 24 - 25	SV	TRAKEOSTOMIHALSBAND Bruksanvisning 50 - 51
HU	TRACHEOSTOMIÁS NYAKPÁNT Használati utasítás 26 - 27	TR	TRAKEOSTOMİ BOYUN BANDI Kullanım kılavuzu 52 - 53
IT	FASCIA DI FISSAGGIO PER CANNULA TRACHEOSTOMICA Istruzioni per l'uso 28 - 29		SYMBOL DESCRIPTION 54 - 57

ASSEMBLY

Deutsch

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REF	Patient target group
40-0005-011	Tracheostomy necktape, blue, infant / child
40-0005-033	Tracheostomy necktape, blue, child / adult
40-0005-044	Tracheostomy necktape, blue, adult, extra large

Clinical Benefit: The Tracheostomy Necktape prevents accidental extubation.

Place of use: Hospitals, nursing homes and home care

INDICATIONS

- Fixation of tracheostomy tubes for tracheostomised patients.

Further indications are not known.

CONTRAINDICATIONS

None known.

SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product and keep it for future reference.

- The product may only be used by medically trained personnel who have sufficient knowledge in handling the product.
- The user and / or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the respective country if an event occurs outside of the EU) in which the user and / or patient is located.
- Prior to use, a functional check must be performed (see section "Functional check"). Damaged products must not be used.

- No changes may be made to the product.
- The product is not sterile.
- The product is intended for single use and must not be reprocessed. Reprocessing influences the function of the product. Any reuse entails a potential risk of cross-contamination.
- The product must not be used if the packaging is damaged or the expiration date has expired.

FUNCTIONAL CHECK

- ▶ Check the product for damages and loose particles.
- ▶ Test the solidity of the product's two stitchings. Do not use excessive force when doing so.

Damaged products must not be used.

USE

Assembly

- ▶ Undo both fasteners from the Tracheostomy Necktape.
- ▶ Place the Tracheostomy Necktape centrally under the patient's neck.
- ▶ Guide both fasteners through the flange of the tracheostomy tube.

CAUTION

Do not place the Tracheostomy Necktape too tightly around the patient's neck because of the risk of decubitus formation or vasoconstriction.

- ▶ Fix both fasteners to the Tracheostomy Necktape.

Monitoring

Patient transfer, patient movement, moisture/humidity and prolonged use may cause the Tracheostomy Necktape to become loose.

- ▶ Regularly check that the Tracheostomy Necktape is in position and adjust if necessary.

If skin irritations occur in the area of the Tracheostomy Necktape, the Tracheostomy Necktape must be replaced with an alternative product.

Removal

The Tracheostomy Necktape must be removed / replaced after a maximum application period of 30 days or if damaged.

- ▶ Undo both fasteners from the Tracheostomy Necktape and remove them from the flange of the tracheostomy tube.

SHELF LIFE

Expiration date: See product label

STORAGE AND TRANSPORT CONDITIONS

	CAUTION
	<ul style="list-style-type: none">• Protect from heat and store in a dry place.• Keep away from sunlight and light sources.• Store and transport in the original packaging.
	

DISPOSAL

The used or damaged products must be disposed of in accordance with the applicable national and international statutory regulations.

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