
Tracheostomy Necktape

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Einsteinstrasse 1 | 72172 Sulz a.N. | Germany | Tel: +49 7454 9596-0

Fax: +49 7454 9596-33 | e-mail: info@vbm-medical.de | www.vbm-medical.de

DE	TRACHEOSTOMIE-HALSBAND Gebrauchsanweisung	5
EN	TRACHEOSTOMY NECKTAPE Instructions for use.....	9
BG	ПОДДЪРЖАЩА ЛЕНТА ЗА ТРАХЕОСТОМИЯ Инструкции за употреба	13
CS	TRACHEOSTOMICKÁ KRČNÍ PÁSKA Návod k použití.....	18
DA	TRACHEOSTOMIKRAVE Brugsanvisning.....	22
EL	ΚΟΛΑΡΟ ΤΡΑΧΕΙΟΣΤΟΜΙΑΣ Οδηγίες χρήσης.....	25
ES	COLLARÍN PARA TRAQUEOSTOMÍA Instrucciones de uso	30
ET	TRAHEOSTOOMIA KAELALINT Kasutamisejuhised	35
FI	TRAKEOSTOMIAKAULANAUHA Käyttöohje	39
FR	COLLIER DE TRACHÉOTOMIE Manuel d'utilisation	43
HR	VRATNA TRAKA ZA TRAHEOTOMIJU Upute za uporabu	47
HU	TRACHEOSTOMIÁS NYAKPÁNT Használati utasítás	51
IT	FASCIA DI FISSAGGIO PER CANNULA TRACHEOSTOMICA Istruzioni per l'uso.....	55

LT	TRACHEOSTOMINĖ APYKAKLĖ Naudojimo instrukcija.....	60
LV	TRAHEOSTOMIJAS KAKLA APSĒJS Lietošanas instrukcija	64
NL	TRACHEOTOMIE-NEKBAND Gebruiksaanwijzing	68
NO	TRAKEOSTOMI-HALSBÅND Bruksanvisning.....	72
PL	OPASKA MOCUJĄCA RURKĘ TRACHEOSTOMIJNĄ Instrukcja użycia.....	76
PT	COLAR DE TRAQUEOSTOMIA Instruções de utilização	80
RO	COLIER DE TRAHEOSTOMIE Instrucțiuni de utilizare	85
RU	ТРАХЕОСТОМИЧЕСКИЙ ФИКСАТОР Инструкция по применению.....	89
SK	TRACHEOSTOMICKÁ FIXAČNÁ PÁSKA Návod na použitie	93
SL	TRAHEOSTOMSKI OVRATNI TRAK Navodila za uporabo	98
SV	TRAKEOSTOMIHalsband Bruksanvisning.....	102
TR	TRAKEOSTOMI BOYUN BANDI Kullanma kılavuzu	106
	SYMBOL DESCRIPTION	110

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REF	Patient target group
40-0005-011	Tracheostomy Necktape, blue, infant / child
40-0005-033	Tracheostomy Necktape, blue, child / adult
40-0005-044	Tracheostomy Necktape, blue, adult, extra large

Clinical benefit: The Tracheostomy Necktape prevents accidental extubation.

Place of use: Hospitals, nursing homes and home care

INDICATIONS

- Fixation of tracheostomy tubes for tracheostomised patients.

Further indications are not known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep it for future reference.

- The product may only be used by medically trained personnel who have sufficient knowledge in handling the product.
- The user and / or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the respective country if an event occurs outside of the EU) in which the user and / or patient is located.
- Prior to use, a functional check must be performed (see section "Functional check"). Damaged products must not be used.
- No changes may be made to the product.
- The product is not sterile.
- The product is intended for single use and must not be reused and/or processed. The function of the product is impaired by processing. Any reuse entails a potential infection hazard.
- The product must not be used if the packaging is damaged or the expiration date has expired.

FUNCTIONAL CHECK

- ▶ Check the product for damages and loose particles.

- ▶ Test the solidity of the product's two stitchings. Do not use excessive force when doing so.

Damaged products must not be used.

USE

Assembly

- ▶ Undo both fasteners from the Tracheostomy Necktape.
- ▶ Place the Tracheostomy Necktape centrally under the patient's neck.
- ▶ Guide both fasteners through the flange of the tracheostomy tube.

CAUTION

Do not place the Tracheostomy Necktape too tightly around the patient's neck because of the risk of decubitus formation or vasoconstriction.

- ▶ Fix both fasteners to the Tracheostomy Necktape.

Monitoring

Patient transfer, patient movement, moisture/humidity and prolonged use may cause the Tracheostomy Necktape to become loose.

- ▶ Regularly check that the Tracheostomy Necktape is in position and adjust if necessary.

If skin irritations occur in the area of the Tracheostomy Necktape, the Tracheostomy Necktape must be replaced with an alternative product.

Removal

The Tracheostomy Necktape must be removed / replaced after a maximum application period of 30 days or if damaged.

- ▶ Undo both fasteners from the Tracheostomy Necktape and remove them from the flange of the tracheostomy tube.

SHELF LIFE

Expiration date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

The used or damaged products must be disposed of in accordance with the applicable national and international statutory regulations.

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