



Tourniquet Wipe Cuff

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Einsteinstrasse 1 | 72172 Sulz a.N. | Germany | Tel.: +49 7454 9596-0
Fax: +49 7454 9596-33 | e-mail: info@vbm-medical.de | www.vbm-medical.de

DE	TOURNIQUET WIPE CUFF Gebrauchsanweisung	5
EN	TOURNIQUET WIPE CUFF Instructions for use.....	8
BG	TOURNIQUET WIPE CUFF Инструкции за употреба	12
CS	TOURNIQUET WIPE CUFF Návod k použití.....	16
DA	TOURNIQUET WIPE CUFF Brugsanvisning.....	19
EL	TOURNIQUET WIPE CUFF Οδηγίες χρήσης.....	23
ES	TOURNIQUET WIPE CUFF Instrucciones de uso.....	27
ET	TOURNIQUET WIPE CUFF Kasutamisujuhised.....	30
FI	TOURNIQUET WIPE CUFF Käyttöohje	34
FR	TOURNIQUET WIPE CUFF Manuel d'utilisation	37
HR	TOURNIQUET WIPE CUFF Upute za uporabu	41
HU	TOURNIQUET WIPE CUFF Használati utasítás	45
IT	TOURNIQUET WIPE CUFF Istruzioni per l'uso.....	48

LT	TOURNIQUET WIPE CUFF Naudojimo instrukcija.....	52
LV	TOURNIQUET WIPE CUFF Lietošanas instrukcija	55
NL	TOURNIQUET WIPE CUFF Gebruiksaanwijzing.....	59
NO	TOURNIQUET WIPE CUFF Bruksanvisning.....	62
PL	TOURNIQUET WIPE CUFF Instrukcja użycia.....	66
PT	TOURNIQUET WIPE CUFF Instruções de utilização	69
RO	TOURNIQUET WIPE CUFF Instrucțiuni de utilizare	73
RU	TOURNIQUET WIPE CUFF Инструкция по применению	77
SK	TOURNIQUET WIPE CUFF Návod na použitie	81
SL	TOURNIQUET WIPE CUFF Navodila za uporabo	84
SV	TOURNIQUET WIPE CUFF Bruksanvisning.....	88
TR	TOURNIQUET WIPE CUFF Kullanma kılavuzu	91
	SYMBOL DESCRIPTION	95

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Instructions for use

INTENDED USE

The Tourniquet Wipe Cuff temporarily occludes the blood flow of a patient's upper or lower extremity in order to obtain a bloodless field.

Clinical benefit: The Tourniquet Wipe Cuff ensures a bloodless surgical field and reduces blood loss during extremity surgery.

Patient target group: Patients who require surgery at upper or lower extremities.

Place of use: Rooms which are suitable for surgical interventions.

INDICATIONS

The indication depends on the intervention and therefore on the selected Tourniquet Wipe Cuff.

Possible indications:

- Reduction of certain fractures
- Arthroscopy of knee, hand, finger or elbow
- Bone grafts
- Kirschner wire removal
- Traumatic or nontraumatic amputations
- Tumour and cyst excisions
- Subcutaneous fasciotomy
- Nerve injuries
- Tendon repair
- Replacement or revision of joints (knee, wrist or finger)
- Correction of a hammer toe
- Podiatry

No other indications are known.

CONTRAINDICATIONS

The contraindication depends on the intervention and therefore on the selected Tourniquet Wipe Cuff.

Possible contraindications:

- Open fractures of extremities
- Post-traumatic lengthy hand reconstruction
- Severe crushing
- Elbow surgery where there is excess swelling
- Severe hypertension
- Skin transplant
- Compromised circulation (e.g. peripheral artery disease)
- Diabetes mellitus

No other contraindications are known.

The physician must, in each individual case, assess the indications and contraindications in the light of his or her expert knowledge.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by a physician, or by medically trained personnel working under the instruction of a physician.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before every use, the product must be inspected visually for damages and a functional check must be performed (see "Visual inspection" and "Functional check"). A faulty product must not be used.
- It is not permitted to make any changes to the product.
- The Tourniquet Wipe Cuff is designed and tested for use with Tourniquets from VBM Medizintechnik GmbH. If the user uses a Tourniquet from another manufacturer, VBM Medizintechnik GmbH shall not accept any liability for the product.
- The Tourniquet Wipe Cuff is designed for a max. cuff pressure of 600 mmHg when in place on the patient. If exposed to a higher cuff pressure, the Tourniquet Wipe Cuff will be irreversibly damaged.
- Prior to every use, the product must be disinfected (see "Wipe disinfection").
- The product is not sterile.
- The product must not be used if the expiration date has expired.

PREPARATION PRIOR TO USE

- ▶ Choose the correct cuff size for the respective extremity (see "Product specifications").

VISUAL INSPECTION

- ▶ Check the product for damages (cracks, breaks, kinks) and loose particles.
- ▶ If present, remove loose particles.
- ▶ Check that the fastener is free of foreign bodies and holds sufficiently.

FUNCTIONAL CHECK

The functional check of the Tourniquet Wipe Cuff must be done with the Tourniquet.

- ▶ Follow the instructions for use of the Tourniquet.
- ▶ Wrap the Tourniquet Wipe Cuff tightly to create a counter-pressure when inflated.
- ▶ Connect the cuff tube to the Tourniquet.

The tube connection must snap firmly into place.

- ▶ Inflate the Tourniquet Wipe Cuff.

The entire system must be free of any air leakage.

A faulty product must be disposed of (see "Disposal").

USE

CAUTION

- The duration of the blood arrest is at the physician's discretion. A max. of 2 hours is normally recommended.
- In order to ensure reliable blood arrest while avoiding exposure of the patient to unnecessary stress, an appropriate cuff target pressure must be selected for the Tourniquet Wipe Cuff, depending on the cuff size, the extremity and the systolic blood pressure.
- The user must check that the selected cuff pressure is sufficient to ensure interruption of the blood flow during the entire procedure.
- Avoid kinking of cuff tube.

- ▶ Apply thin padding under the Tourniquet Wipe Cuff to avoid pressure sores and injury to the skin.
- In order to ensure that the cuff tube does not located in the operating field, it must be directed proximally (Image 1 + 2).

CAUTION

To ensure reliable blood arrest, the bladder of the Tourniquet Wipe Cuff must overlap on the extremity. However, too much overlap can result in cuff rolling.

- ▶ Apply the Tourniquet Wipe Cuff tightly on the extremity, but without using force. Two fingers should fit between the extremity and the Tourniquet Wipe Cuff.
- ▶ After applying the Tourniquet Wipe Cuff, additionally secure it with the colour-coded cuff ties (Image 1 + 2).

- ▶ Connect the cuff tube to the Tourniquet.
- The tube connection must snap firmly into place.
- ▶ Prior to cuff inflation, the extremity must be exsanguinated.
 - ▶ Prior to inflation, check for correct placement of the Tourniquet Wipe Cuff.
 - ▶ Inflate the cuff with the lowest necessary pressure required to occlude the blood flow.
- In order to avoid possible chemical burns, fluids must be prevented from getting under the Tourniquet Wipe Cuff.
- ▶ Disinfect the skin and drape the operating field.
 - ▶ Before commencing the surgical procedure, make sure by means of palpation or auscultation that the arterial blood flow has been occluded.

WARNING

IVRA (intravenous regional anaesthesia)

- Connect the blue cuff tube to the blue connection of the Tourniquet. Connect the red cuff tube to the red connection of the Tourniquet.
- After administration of the local anaesthetic, a minimum occlusion time of 20 minutes must be considered in order to prevent a toxic reaction.
- In the event of pressure loss in the inflated cuff chamber, the second cuff chamber must be inflated immediately.

CAUTION

- ▶ After use, completely deflate the Tourniquet Wipe Cuff and immediately remove cuff and padding from the extremity in order to prevent the risk of venous congestion.

- ▶ Disconnect the cuff tube from the Tourniquet.

WIPE DISINFECTION

CAUTION

- The product must not be mechanically or manually reprocessed or sterilised.
- In case of unsuccessful wipe disinfection or heavy contamination, the product must be disposed of as specified in "Disposal".

Wipe disinfection must be performed with commercially available surface disinfectants based on alcohol, QACs (quaternary ammonium compounds) or hydrogen peroxide. When choosing a product for disinfection, a disinfectant with appropriate ranges of action must be used: bactericidal, yeasticidal and virucidal. After wipe disinfection, inspect the product for visible contamination. If necessary, repeat wipe disinfection. After wipe disinfection, check the function of the product (see "Visual inspection" and "Functional check").

WARNING

By using the product on patients with suspected prion disease there may be a high risk of transmission. In such case the physician decides if the product is disposed of (see "Disposal").

SHELF LIFE

The shelf life of the product is 5 years and it can be reprocessed within the shelf life up to 100 times. The product must then be disposed of (see "Disposal"). Every further reuse is the responsibility of the user.

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must be cleaned and disinfected to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

DISPOSAL

The product must be disposed of in accordance with the applicable national and international legal regulations.

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