


LT[®]evo









CE 0123

Intended to be left blank.

Product description	4	IT Istruzioni per l'uso	30 - 31
Use	5	LT Naudojimo instrukcija	32 - 33
DE Gebrauchsanweisung	6 - 7	LV Lietošanas instrukcija	34 - 35
EN Instructions for use	8 - 9	NL Gebruiksaanwijzing	36 - 37
BG Инструкции за употреба	10 - 11	NO Bruksanvisning	38 - 39
CS Návod k použití	12 - 13	PL Instrukcja użycia	40 - 41
DA Brugsanvisning	14 - 15	PT Instruções de utilização	42 - 43
EL Οδηγίες χρήσης	16 - 17	RO Instrucțiuni de utilizare	44 - 45
ES Instrucciones de uso	18 - 19	RU Инструкция по применению	46 - 47
ET Kasutamisejuhise	20 - 21	SK Návod na použitie	48 - 49
FI Käyttöohje	22 - 23	SL Navodila za uporabo	50 - 51
FR Manuel d'utilisation	24 - 25	SV Bruksanvisning	52 - 53
HR Upute za uporabu	26 - 27	TR Kullanma kılavuzu	54 - 55
HU Használati utasítás	28 - 29	Symbol Description	56 - 59

PRODUCT DESCRIPTION

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INTENDED USE

The LT evi is a supraglottic airway to temporarily secure the airway.

Clinical benefit: ventilation of a patient.

Patient target group: child, adult

Place of use: hospital and prehospital including military use

INDICATIONS

- Difficult Airway Management and cardiopulmonary resuscitation.
 - Mechanical or spontaneous ventilation during routine anaesthesia procedures.
- No other indications are known.

CONTRAINDICATIONS

- Complete obstruction of the upper airway.
- In elective use: patients who have not fasted as well as patients with retained gastric content.

No other contraindications are known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel with adequate skills in the use of the product.

- Only use the product for the patient target group.

- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.

- Before use, the product must be inspected visually and a functional check must be performed (see chapters "visual inspection" and "functional check"). A faulty product must be disposed of.

- Application time: < 4 hours

- Do not damage the cuffs of the tube with sharp-edged objects.

- It is not permitted to make any changes to the product.

- The product is suitable for conducting oxygen, air, anaesthetic gases and their gas mixtures.

- Laser or electrocaustic equipment may cause the product to be flammable.

- Use original VBM parts only. Using parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.

- The product is MR conditional (see chapter "use of tube during an MRI examination").
- The product is intended for single use and must not be reused and/or processed. The function of the product is impaired by processing. Any reuse entails a potential infection hazard.

STERILE EO

- The tube is sterile (ethylene oxide).

- Visually inspect the sterile barrier system for breaches of integrity immediately prior to application.

- The product must not be used if the packaging is damaged or if the shelf life has elapsed.

- The syringe is not sterile.

PRODUCT DESCRIPTION

① 15 mm O.D. connector

② Inflation line for both cuffs

③ Pilot balloon with check valve and Luer connection

④ Valve opener

⑤ Ventilatory opening

⑥ Distal opening of drain tube

⑦ Distal cuff

⑧ Proximal cuff

⑨ Integrated biteblock

⑩ Teeth marks

⑪ Proximal opening of drain tube

ANATOMIC LANDMARKS

① Oesophagus

② Trachea

③ Vocal cords

④ Laryngeal inlet

⑤ Epiglottis

⑥ Tongue

⑦ Incisors

PREPARATION PRIOR TO USE

- ▶ Select appropriate tube size (see table "product specifications").

- ▶ Perform visual inspection, see chapter "visual inspection".

- ▶ Detach inflation line ② from connector ③ of tube and remove valve opener ④.

The valve opener ④ allows vented air to escape from the cuffs ⑨ / ⑩ during sterilization.

- ▶ Perform functional check, see chapter "functional check".

- ▶ Prepare syringe with the recommended volume (see table "product specifications").

- ▶ Lubricate tip and posterior part of the tube with a water-based lubricant.

- ▶ Have a spare tube ready.

If possible, pre-oxygenation is recommended.

VISUAL INSPECTION

- ▶ Check packaging for damage.

- ▶ Check product for damage and loose particles.

A faulty product must be disposed of (see chapter "disposal").

FUNCTIONAL CHECK

- ▶ Check for leaks:

- ▶ Deflate cuffs ⑨ / ⑩ completely (figure 1).

- ▶ Remove syringe and check that a vacuum is created and maintained.

A faulty product must be disposed of (see chapter "disposal").

USE**Tube insertion:****CAUTION**

Ensure that the level of anaesthesia or unconsciousness is sufficient before tube insertion.

- ▶ Position patient's head:

The ideal head position for tube insertion is an extended head and a flexed neck (figure 2).

- ▶ Hold tube in the area of the teeth marks ① with the ventilatory opening ⑤ facing towards the chin of the patient.

- ▶ Open patient's mouth.

The application of jaw thrust or chin lift manoeuvre (figure 3) helps to create a maximal pharyngeal space.

**CAUTION**

The cuffs ⑨ / ⑩ must be deflated completely to avoid damage during tube insertion.

- ▶ Introduce tube tip into the mouth in direction towards the hard palate. Do not use force to insert the tube to avoid injury to the patient.

- ▶ Advance tube into the hypopharynx (figure 4) until a definite resistance is felt. Maintain chin lift or jaw thrust during the whole insertion procedure.

The tip of the tube should be located in the oesophageal opening ①. The teeth marks ① are an indicator for the initial insertion depth. The incisors ⑦ should be positioned within this range.

Cuff inflation:

- ▶ Inflate cuffs ⑨ / ⑩ (figure 5).

**CAUTION**

- ▶ To ensure correct inflation, use the syringe for LTevo, see chapters "scope of delivery" and "accessories".
- ▶ Adjust cuff pressure soonest possible. Ensure that the cuff pressure is below 60 cmH₂O (figure 8).
- ▶ Monitor cuff pressure regularly during use.
The cuff pressure may change in the presence of oxygen or anaesthetic gases.
The cuff pressure may change during repositioning or transport of the patient.

Ventilation:

- ▶ Connect tube to a breathing system (figure 6).

**CAUTION**

- ▶ Use an appropriate procedure to monitor ventilation continuously.
- ▶ Fix tube soonest possible, but in any case before repositioning or transporting the patient (figure 8).
- ▶ Ensure suitable strain relief during ventilation.

Gastric tube insertion:

The drain tube (K) - (I) facilitates gastric pressure relief and the passage of a gastric tube to suction gastric fluids.

- ▶ To control correct positioning, it is recommended to insert a gastric tube via the drain tube (K) - (I) (figure 7) (see table "product specifications").
- ▶ Lubricate gastric tube with water-based lubricant.

If the gastric tube cannot be advanced this may indicate that the tip of the tube is positioned incorrectly.

- ▶ In this case remove gastric tube, deflate cuffs completely and reinsert tube.

Tips for optimization in case of inadequate ventilation:

- Tube inserted too shallow:
 - ▶ Remove gastric tube, deflate cuffs completely and reinsert tube. Ensure that jaw thrust or chin lift is maintained during the whole insertion procedure.
- Tube inserted too deep:
 - ▶ Withdraw tube with inflated cuffs between the teeth marks until ventilation is easy and free flowing (large tidal volume with minimal airway pressure).

Option: fiberoptic intubation

The ventilation channel (A) - (C) can be used as a conduit for fiberoptic intubation by introducing a tracheal tube under fiberoptic guidance.

- ▶ Before use, ensure that the tracheal tube is suitable for the tube size used (see table "product specifications") and is lubricated with a water-based lubricant.

Use of tube during an MRI examination:

- ▶ Fix pilot balloon (C) to the patient using adhesive tape.
- During MRI examination, the pilot balloon (C) can cause artefacts around the area where it is positioned.
- An MR field strength of up to 3 T is permitted.

Removal:

- ▶ Remove tube fixation.
- ▶ Remove gastric tube if used.
- ▶ Deflate cuffs completely.
- ▶ Remove tube where suction equipment and the ability for rapid intubation are present.

The tube is well tolerated until the return of protective reflexes.

SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture.

Use-by date: see product label

STORAGE AND TRANSPORT CONDITIONS**CAUTION**

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

PRODUCT SPECIFICATIONS

Tube size	#2	#3	#4	#5
Patient weight guidance	10 - 25 kg	25 - 50 kg	50 - 90 kg	> 90 kg
Minimum interden-tal gap	13 mm	17 mm	19 mm	19 mm
Recommended cuff inflation volume	30 ml	40 ml	45 ml	50 ml
Connector (A)	15 mm O.D.			
Gastric tube size (via drain tube (K) - (I))	≤ 12 Fr	≤ 16 Fr	≤ 18 Fr	≤ 18 Fr
Tracheal tube size (via ventilation channel (A) - (C))	O.D. ≤ 6.7 mm (I.D. 5 mm), uncuffed	O.D. ≤ 8.7 mm (I.D. 6.5 mm)	O.D. ≤ 10.7 mm (I.D. 8 mm)	O.D. ≤ 10.7 mm (I.D. 8 mm)

Conversion of pressure units: 1 hPa = 1.01973 cmH₂O = 0.75006 mmHg

Conversion of length units: 1 Fr = 0.33 mm

SCOPE OF DELIVERY

REF	Tube size	Syringe included	Packaging unit
32-10-002-1	#2	-	10
32-10-003-1	#3	-	10
32-10-004-1	#4	-	10
32-10-005-1	#5	-	10
32-10-102-1	#2	X	1
32-10-103-1	#3	X	1
32-10-104-1	#4	X	1
32-10-105-1	#5	X	1
32-10-209-1	#3 + #4 + #5	X	1 each

ACCESSORIES

REF	Description	Volume	Packaging unit
54-04-888	Syringe for LTevo	50 ml	10

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

**CAUTION**

The product may be contaminated with potentially infectious substances of human origin.

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