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VERSION A - RECTUS

VERSION B - NIST

## Tourniquet 2500

---

for single cuff operation

CE 0123

DE

EN

ES

FR

IT

G1001 - 26/2022-01

# Deutsch

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## 1 GENERAL INFORMATION

### Instructions for use

Prior to starting up the tourniquet device, carefully read and observe the instructions for use with the safety instructions.

This product must only be used by medically trained personnel.

### Intended use

The Tourniquet 2500 is a pneumatic tourniquet for maintaining a safe bloodless field during surgery on an extremity with a single cuff.

Place of use: Hospitals, operating theatres and medical practices.

### Indication / contraindication

The device only presents indications or contraindications in connection with other medical devices, e.g. tourniquet cuffs. These are listed with the corresponding medical device.

### Requirements for the medical device

The medical device meets the essential requirements of the Directive of the European Parliament and the Council about medical devices 93/42/EEC and 2007/47/EC.

## INFORMATION

- Disposal within the EU: Do not dispose of electrical or electronic equipment in household waste. Disposal within the EU must be carried out in accordance with Directive 2012/19/EU (WEEE Directive). The manufacturer will accept the electrical devices for disposal.
- Disposal outside of the EU: The devices and spare parts must be disposed of in accordance with the applicable local provisions.
- More extensive maintenance measures that are not described in these instructions may only be done by the manufacturer or the authorised service agent.

## 2 DEVICE DELIVERY

### Scope of delivery

#### Version A Pressure connection RECTUS:

- Hand inflator for manual use
- Cover
- Colour-coded coil tubing for tourniquet cuff
- 4 m (13 ft) pressure hose

### Scope of delivery

#### Version B Pressure connection NIST:

- Hand inflator for manual use
- Cover
- Colour-coded coil tubing for tourniquet cuff

### CAUTION

Does NOT contain a pressure hose.

### Mounting types:

#### REF xx-12-xxx

Table unit without cover

#### REF xx-13-xxx

Device with cover and universal clamp

#### REF xx-15-xxx

Device on mobile stand with basket and cover


#### REF xx-22-xxx

Table unit on mobile stand with basket, without cover

Tourniquet cuffs and other accessories are not included in the scope of delivery and must be ordered separately. Detailed product information can be obtained from the manufacturer.

### 3 ASSEMBLY AND PREPARATION FOR USE

#### Assembly of the mobile stand with basket REF xx-15-xxx

- Insert lockable castors (3) opposite each other into the base (1).
- Insert simple castors (2) into the base (1).
- Insert the stand tube (4) into the base (1).
- Insert the cylinder head screw M8x40 (7) with washer (6) and rotation lock (5) from below into the base (1) and tighten with Allen key (8).
- The bracket must be screwed to the basket with the two countersunk screws M5x12.


- The bracket with basket must be mounted to the upper part of the case with the two hexagon bolts M6x10 (insert nuts).

- Place the tourniquet device with the bracket and basket on the stand tube.
- Insert the enclosed cylinder head screw with internal hex M8x25 together with the lock washer through the hole in the bracket into the threaded hole of the stand tube. Tighten with the Allen key.

#### Assembly of the mobile stand with basket for table units REF xx-22-xxx

- Fasten the square rail to the stand tube with two Allen screws (a).
- Fasten the fixation plate to the stand tube with countersunk screw M8x20 (b).

- Place the device on the fixation plate. Fasten using cap nut M6 and spanner WS 10.


- Assemble the steel bracket with fixation plate (c) to the basket with 2 countersunk screws M5x12.

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### 3 ASSEMBLY AND PREPARATION FOR USE

#### Assembly of the universal clamp (REF xx-13-xxx)

Screw the universal clamp with the enclosed cylinder head screws with internal hex M6x25 either vertically or horizontally to the upper part of the case (insert nuts). The universal clamp thus enables the tourniquet device to be assembled on any tabletop, operating table, rail or stand tube.

#### Assembly and disassembly of the cover

The cover must be removed as follows:

- Tilt the cover slightly (approx. 30°).
- On one side of the cover recess, use your finger to press the cover from back to front until the pin comes out of the receptacle hole.

The cover must be assembled on the upper part of the case as follows:

- On one side, insert the guiding pin of the cover into the receptacle hole. Ensure that the blue coil tubing runs backwards through the cover recess.
- On the other side, position the pin in the guiding groove of the bearing eye and use your finger to press it backwards and lock it into place in the receptacle hole.

#### Compressed air connection

#### CAUTION

The tourniquet device may only be operated with compressed air.

#### Version A

##### Pressure connection RECTUS with pressure hose

Connect the pressure hose using the positive locking connector to the plug connection on the rear of the device.

Connect the other end of the pressure hose to a central compressed air supply using the corresponding plug connection.

When connected to the central compressed air supply (operating pressure of min. 2 to max. 10 bar), the pneumatic indicator eye on the front plate of the device changes to green. The tourniquet device is ready for operation.

#### Version B

##### Pressure connection NIST

Connect the pressure hose to the NIST connection on the rear of the device.

Connect the other end of the pressure hose to a central compressed air supply using the corresponding plug connection.

When connected to the central compressed air supply (operating pressure of min. 2 to max. 10 bar), the pneumatic indicator eye on the front plate of the device changes to green. The tourniquet device is ready for operation.

#### 4 TECHNICAL DATA

Weight (table unit without pressure hose)	2.26 kg	(80 oz)
Dimensions		
Height	140 mm	(6 inch)
Width	210 mm	(8 inch)
Depth	180 mm	(7 inch)
Operating pressure	2 - 10 bar	(29 - 145 psi)
Device's pressure connection		
Version A	RECTUS	
Version B	NIST AIR	
Adjustable range for tourniquet	0 - 600 mmHg ( $\pm 10\%$ )	
Adjustment accuracy	$\pm 3$ mmHg	
Set pressure limitation for tourniquet	600 mmHg ( $\pm 10\%$ )	
Display accuracy of the device manometer	$\pm 10$ mmHg	
Time alarm	Alarm sound for countdown to zero	
Connection of tourniquet system	Blue coil tubing with female positive locking connectors	

#### Environmental conditions

Transport	-10 ... +60 °C	(14 ... 140 °F)
	30 ... 95 % relative air humidity without condensation	
Storage / Operation	+10 ... +40 °C	(50 ... 104 °F)
	30 ... 95 % relative air humidity without condensation	

#### Shelf life

##### Tourniquet

The shelf life of the device is 8 years with proper use.

Date of manufacture: see rating plate.

##### Pressure hose

Expiry date: see print on the pressure hose.

## 5 DEVICE OVERVIEW

**a) Colour-coded coil tubing (blue)**

With positive locking connectors; fastened directly to the device, so no leaks

**b) Manometer**

Pressure range from 0 - 600 mmHg; always displays the actual cuff pressure

**c) Timer**

Countdown with alarm sound

**d) Indicator eye**

Compressed air / Manual

**e) Manual safety system with hand inflator**

In the event of compressed air failure, automatically switches to manual operation

**f) Deflate buttons (blue)**

Upon pressing, the cuff is immediately deflated, upon releasing it is inflated to the set value; for checking haemorrhaging after surgery; enables complete deflation in manual mode

**g) Precision pressure regulator knob**

Precise pressure regulation, automatic pressure compensation in the event of leaks, protective cap to stop accidental turning; pressure limitation to 600 mmHg ( $\pm 10\%$ )

## 6 OPERATING INSTRUCTIONS TOURNIQUET DEVICE

### CAUTION

In order to ensure a safe bloodless field while avoiding exposure of the patient to unnecessary stress, an appropriate cuff pressure must be selected, depending on the cuff size and the systolic blood pressure in the extremity. The duration of the bloodless field must not exceed two hours.

### CAUTION

A countdown of 2 hours is not possible! The timer can be preset to a max. of 99 min 59 sec. The countdown duration must be set again for a two-hour bloodless field (see "Instructions for use, 8 Operating Instructions Timer").

### After attaching and securing the single cuff, the tourniquet device is operated as follows:

1. Reset the precision pressure regulator knob to zero by turning counterclockwise as far as it will go.
2. Connect the connection hose of the single cuff and the coil tubing of the device with the positive locking connector.
3. After exsanguination, fill the single cuff to the desired cuff pressure or block by turning the knob. The current cuff pressure is continuously displayed on the right manometer and must be continuously monitored.
4. The cuff pressure can be regulated up (e.g. in the case of secondary haemorrhage) or down at any time with the precision pressure regulator knob.
5. After termination of surgery, the cuff must be deflated by slowly turning (counterclockwise to 0) the knob.

### Manual operation (safety system)

In the event of compressed air failure or disconnection from the central compressed air supply (e.g. for patient transport), the device automatically switches to manual operation mode (see indicator eye). The device continues to function after this. This means:

1. The currently set cuff pressure remains unchanged in the tourniquet system.
2. Inflation of the cuff now continues with the hand inflator.
3. The cuff chamber can be deflated by pressing the corresponding deflate button.

### CAUTION

#### Version A

#### Pressure connection RECTUS with pressure hose:

For devices that are operated with a wall gas probe with one-way valve (NOT recommended), it must be ensured that the pressure hose is disconnected from the device, as a sudden pressure drop in the cuff cannot be ruled out in the event of disconnection at the wall gas probe.

#### Version B

#### Pressure connection NIST:

For devices that are operated with a NIST connector and wall gas probe with one-way valve (NOT recommended), it must be ensured that the pressure hose is first disconnected from the wall gas probe and then from the device, as a sudden pressure drop in the cuff cannot be ruled out.

## 7 SAFETY INSTRUCTIONS

### Compressed air connection

The tourniquet device may only be operated with compressed air!

### Prior to use

Prior to each use of the tourniquet device, the functionality and tightness of the entire tourniquet system must be checked. To do so, the cuff must be attached tightly to a suitable object and filled to the maximum pressure of 500 mmHg. Then the pressure hose must be disconnected from the device. The indicator eye then changes from the green indicator "Compressed air" to the white indicator "Manual". Consistency of pressure must be checked on the analogue manometer. The pressure drop must not be greater than 30 mmHg in 5 minutes. If a greater pressure drop occurs, the device must be inspected (see "Service Manual, 1 Device Check").


### CAUTION

- Use the correct cuff size for the respective extremity. The manufacturer offers a complete selection of single cuffs.
- If cuffs from other manufacturers are used, compatibility must be ensured.
- Damaged cuffs and connections must not be used.
- Before connecting the cuff to the tourniquet device, firstly check that the knob has been reset to zero.
- The connecting tubes between the device and cuff must not be disconnected.
- While filling the cuff, carry out tactile inspection with your hand to check that the cuff chamber is being properly filled.
- During surgery, the tourniquet system must be continuously checked on the manometer. The manometer always displays the actual cuff pressure. Any pressure drop can be read on the manometer.

### Manual operation

- Manual operation and operation with the precision pressure regulator knob work independently from each other. Irrespective of the cuff pressure built up using the hand inflator, the cuff assumes the value set with the knob when connected to the central compressed air supply. This is why the knob's arrow position must be compared with the mmHg scale on the front of the device.
- When using manual operation, any leaks will not be offset automatically and can only be compensated for through manual inflation. Regular operation of the device should therefore occur with compressed air. Manual operation is only for short-time patient transport or as a backup in the event of compressed air failure.

## 8 OPERATING INSTRUCTIONS TIMER

### PRODUCT DESCRIPTION

- a Minute display
- b Seconds display
- c Start/pause button
- d Minute setting button (MIN)
- e Seconds setting button (SEC)
- f RESET button

Use the Timer to set the timing (count-up and count-down) as follows:

### USE

#### Setting the count-up

- ▶ Press the RESET button to set the display to 00:00.
- ▶ Press the Start/Pause button to start or temporarily stop the count-up.
- ▶ During count-up, ":" flashes in the display.

#### Setting count-down with alarm

- ▶ Press the RESET button to set the display to 00:00.
- ▶ Press the MIN or SEC button to set the time. The Timer can be preset to a maximum of 99 MIN 59 SEC.
- ▶ Press the Start/Pause button to start or temporarily stop the count-down.
- ▶ When the time has elapsed, an acoustic signal sounds.
- ▶ To switch off the alarm, press the Start/Pause button. At the same time, the preset time appears again in the display and can be started again by repeatedly pressing the Start/Pause button.

## 9 CLEANING INSTRUCTIONS

### Wipe disinfection

To disinfect the device, a wipe disinfection with non-flammable, commercially available surface disinfectants (e.g. Sani Cloth® Active from Ecolab Deutschland GmbH) must be performed. When choosing a product for disinfection, a disinfectant with appropriate ranges of action must be used: bactericidal, yeasticidal, tuberculocidal and virucidal. After the wipe disinfection, the products must be inspected for visible contamination. If necessary, repeat the wipe disinfection. After the wipe disinfection, the products must be checked in accordance with the section "Instructions for use, 1 Device Check".

The device cannot be mechanically reprocessed or sterilised.


# SERVICE MANUAL

## 1 DEVICE CHECK

### Functional and air-tight test

The following test must be carried out prior to each use and / or as needed (see "Instructions for use, 7 Safety Instructions, Prior to use", "Service Manual, 2 Troubleshooting"):

1. Plug the blue coil tubing with the plugs (see a).
2. Connect the device to compressed air supply. Indicator eye changes from white to green (see b).
3. Any pressure value must be adjustable with a precision of  $\pm 3$  mmHg with the pressure regulator (see c) and this must be able to be read on the manometer (see e).
4. Set the pressure to 500 mmHg.
5. Briefly press the blue deflate button (see d). Pressure drops to 0 mmHg; upon release, the preset value of 500 mmHg must be reached again (see e).
6. Set the pressure to 0 mmHg and disconnect the device from the compressed air supply.
7. Inflate to 500 mmHg using the hand inflator (see f) (can be read on the manometer, see e).
8. Set the timer (see g) to a 5 min countdown. After the 5 min have expired, the pressure drop on the manometer must not be greater than 30 mmHg. In the case of a greater pressure drop, see "Service Manual, 2 Troubleshooting, Leak in tourniquet system".

# 1 DEVICE CHECK

## Maintenance

The medical device may only be serviced by persons who have the technical knowledge, qualification and the required means to properly carry out the service.

After servicing the medical device, the structural and functional characteristics essential for safety and functionality must be checked, insofar as they can be influenced by the service procedure.

All guarantee or warranty claims are null and void if the user or an unauthorised body tries to remedy a defect themselves, which is not described in these instructions for use.

For a fast repair process, it is required that the medical device be returned with the most precise a fault description as possible.

To protect employees, every medical device must be thoroughly cleaned and disinfected before returning (cf. "Instructions for use, 9 Cleaning Instructions"). The manufacturer can refuse to repair soiled or contaminated medical devices for reasons of safety.

## Manufacturer inspection

It is recommended to return the device to the manufacturer for inspection after an operating period of 5 years.

## Metrological check

The accuracy  $\pm 10$  mmHg must be checked every 12 months.

If there is an indication that the tourniquet cannot keep the accuracy of measurement within the 12 months, or the metrological characteristics of the tourniquet have been affected, a metrological check has to be carried out immediately. This must only be carried out by a person who meets the requirements for metrological checks. The results of the check, as well as a minimum of the measured values and measurement procedures, have to be documented.

After successful check, the tourniquet has to be marked with a sign that indicates (clearly and traceable) the year of the next inspection and the person who carried out the metrological check.

In order to check the values, the female positive locking connector of the tourniquet has to be connected to a calibrated pressure gauge (measurement inaccuracy max.  $\pm 2$  mmHg).

Inflate the manometer with pressure range 600 mmHg to 100, 200, 300, 400, 500, 600 mmHg and note the values. The values must be within the tolerance. If this is not the case, the device must be returned to the manufacturer for repair. If the values are within the tolerance, the device can still be used.

If required measuring equipment is not available, the metrological check can be carried out by the manufacturer. In this case the tourniquet has to be sent to the manufacturer.

Note for conversion of mmHg to SI unit Pascal:


100 mmHg corresponds approximately to 13.33 kPa

## 2 TROUBLESHOOTING

Fault / failure	Cause / troubleshooting
<p><b>Indicator eye</b></p> <ul style="list-style-type: none"> <li>• Indicator eye does not change to green when the compressed air supply is connected.</li> <li>• Indicator eye stays green after disconnecting from the compressed air supply.</li> </ul>	<ul style="list-style-type: none"> <li>▶ Input pressure of the compressed air supply is missing or is too low. Provide input pressure of min. 2 bar to max. 10 bar.</li> <li>▶ Indicator eye is faulty. Repair by manufacturer.</li> <li>▶ Compressed air is only slowly released from the pressure hose. A wall gas probe with one-way valve should not be used.</li> <li>▶ Indicator eye is faulty. Repair by manufacturer.</li> </ul>
<p><b>Manometer</b></p> <ul style="list-style-type: none"> <li>• Pointer does not respond.</li> <li>• Pointer not exactly at zero (tolerance <math>\pm 10</math> mmHg).</li> <li>• Glass is damaged.</li> </ul>	<ul style="list-style-type: none"> <li>▶ Compressed air is missing (see indicator eye).</li> <li>▶ Pointer is stuck or loose. Replace manometer (see "Service Manual, 3 Replacement of Parts, Replacement of manometer").</li> <li>▶ Unit faulty. Completely replace manometer (see "Service Manual, 3 Replacement of Parts, Replacement of manometer").</li> <li>▶ Improper use. Replace glass (see "Service Manual, 3 Replacement of Parts, Replacement of manometer").</li> </ul>
<p><b>Pressure regulator</b></p> <ul style="list-style-type: none"> <li>• Pressure regulator cannot be turned anymore.</li> <li>• Pressure no longer reaches max. value (600 mmHg <math>\pm 10</math> %).</li> </ul>	<ul style="list-style-type: none"> <li>▶ Pressure regulator overwound. Loosen pressure regulator knob in accordance with instructions (see "Service Manual, 3 Replacement of Parts, Pressure regulator knob cannot be moved anymore").</li> <li>▶ Pressure regulator faulty. Repair by manufacturer.</li> </ul>
<p><b>Deflate button</b></p> <ul style="list-style-type: none"> <li>• Deflation with blue deflate button not possible.</li> </ul>	<ul style="list-style-type: none"> <li>▶ Deflate valve is jammed or faulty. Completely replace valve (see "Service Manual, 3 Replacement of Parts, Replacement of deflate valve").</li> </ul>
<p><b>Timer</b></p> <ul style="list-style-type: none"> <li>• Timer display is too weak or completely missing.</li> </ul>	<ul style="list-style-type: none"> <li>▶ Battery empty. Change battery (see "Service Manual, 3 Replacement of Parts, Replacement of battery").</li> </ul>
<p><b>Leak in tourniquet system</b></p>	<ul style="list-style-type: none"> <li>▶ Connections not properly connected. Reconnect connections.</li> <li>▶ Valves (deflate, one-way, 3/2-way valve) are not tight. Troubleshooting and repair by manufacturer.</li> </ul>

### 3 REPLACEMENT OF PARTS

#### Open the case


1. For the standing device REF xx-15-xxx, the device must firstly be removed from the stand tube. To do this, unscrew the cylinder head screw with internal hex M8x25 (see a).
2. Remove the basket by loosening the two countersunk screws with cross recess M5x12 (see b).

3. Tilt the cover slightly (approx. 30°): On one side of the cover recess, use your finger to press the cover from back to front until the pin comes out of the receptacle hole.


4. For table units with REF xx-22-xxx, the device must firstly be removed from the fixation plate. To do this, remove the two cap nuts.

5. Unscrew the four cylinder head screws with cross recess M4x10 on the right side (see a) and the four cylinder head screws on the left side.
6. Now pull the upper part of the case upwards and lower part downwards.


#### Replacement of timer

1. Open the case (see "Service Manual, 3 Replacement of Parts, Open the case").
2. Disconnect electrical connections (see a).
3. Remove the entire timer backwards by loosening the four crosshead screws (see b).


### 3 REPLACEMENT OF PARTS

#### Replacement of manometer


1. Open the case (see "Service Manual, 3 Replacement of Parts, Open the case").
2. Cut the hose on the manometer.
3. Loosen the nut (WS 11) on the bracket (see a). Remove the manometer forwards through the front plate.
4. To replace the glass, loosen the outer metal ring counterclockwise.

#### Replacement of battery

1. Open the battery compartment cover on the rear of the device using a screwdriver (see a).
2. Remove the battery using the screwdriver and replace with a new battery (AA 1.5 V) (check polarity).

#### Replacement of deflate valve


1. Open the case (see "Service Manual, 3 Replacement of Parts, Open the case").
2. Cut hose at deflate valve.

3. Using a screwdriver, press the flap on the valve body upwards (see a). Valve can be removed from the push-button.

#### Pressure regulator knob cannot be moved anymore

1. It must firstly be determined if the pressure regulator knob has been turned too far to the right or too far to the left. To do this, plug the coil tubing with plugs and connect the device to the compressed air supply.
2. If the manometer shows 0 mmHg pressure (see a), the pressure regulator knob must be turned clockwise (see b).

3. If the manometer shows 600 mmHg ( $\pm 10\%$ ) pressure (see c), the pressure regulator knob must be turned counterclockwise (see d).

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
# SPARE PARTS LIST

Version:  
REF 15-12-500  
REF 15-22-500


1. Case, upper part **REF 10-16-722**, (Box 1)
2. Case, lower part **REF 10-16-721**, (Box 1)
3. Mobile Stand, complete  
with fixation plate, square rail and basket  
**REF 10-22-555**, (Box 1)
4. Square Rail **REF 10-22-444**, (Box 1)
5. Basket with stainless steel hook  
**REF 51-20-007**, (Box 1)

Coil Extension Hose with spring,  
without connectors (Box 1 pc. in 3 m length)  
colour: blue  
**REF 20-10-444**


Version:  
REF 15-13-500  
REF 15-15-500


6. Cover transparent **REF 10-15-710**, (Box 1)
7. Case, lower part **REF 10-25-722**, (Box 1)
8. Case, upper part **REF 10-25-721**, (Box 1)
9. Mobile Stand, complete with bracket and basket  
**REF 10-15-555**, (Box 1)
10. Basket separate **REF 10-15-556**, (Box 1)
11. Universal Clamp **REF 51-53-555**, (Box 1)
12. Bracket with screws **REF 10-15-558**, (Box 1)

Coil Extension Hose with positive  
locking connectors (Box 1 pc. in 3 m length)  
colour: blue  
**REF 20-20-744**


# SPARE PARTS LIST

3/2-way-valve  
**REF 20-50-310**  
(Box 1)

Positive Locking Connector, female  
**REF SLZF40**  
(Box 5)


Deflate button, blue  
with valve  
**REF 20-50-550**  
(Box 1)

Plug for air-tight test  
of Tourniquets  
**REF 22-50-406**  
(Box 1)

Timer for external compartment  
complete with battery (Box 1 each)  
**REF 15-70-805**

Spare Battery **REF 15-50-999**

Hand Inflator  
with tube and spring  
**REF 20-10-333**  
(Box 1)




One-Way-Valve  
for hand inflator  
**REF 20-50-576-VK**  
(Box 1)

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