

Esmarch Bandage reusable

DE	ESMARCH BANDAGE WIEDERVERWENDBAR Gebrauchsanweisung 4 - 5	LT	ESMARCHO VYTURAS, DAUGKARTINIO NAUDOJIMO Naudojimo instrukcija..... 30 - 31
EN	ESMARCH BANDAGE REUSABLE Instructions for use..... 6 - 7	LV	ATKĀRTOTI IZMANTOJAMA ESMARHA BANDŽĀ Lietošanas instrukcija 32 - 33
BG	ЛЕНТА ЗА КРЪВОСПИРАНЕ (ПО ЕСМАРХ) ЗА МНОГОКРАТНА УПОТРЕБА Инструкции за употреба 8 - 9	NL	ESMARCH BANDAGE, HERBRUIKBAAR Gebruiksaanwijzing 34 - 35
CS	ESMARCHOVO ŠKRTIDLO, OPAKOVANĚ POUŽITELNÉ Návod k použití..... 10 - 11	NO	ESMARCH-BANDASJE TIL GJENBRUK Bruksanvisning..... 36 - 37
DA	GENANVENDELIG ESMARCH-BANDAGE Bruganvisning 12 - 13	PL	OPASKA ESMARCHA WIELOKROTNEGO UŻYTKU Instrukcja użycia..... 38 - 39
EL	ΕΠΙΔΕΣΜΟΣ ΕΣΜΑΡΧ ΕΠΑΝΑΧΡΗΣΙΜΟΠΟΙΟΥΜΕΝΟΣ Οδηγίες χρήσης..... 14 - 15	PT	LIGADURA ESMARCH REUTILIZÁVEL Instruções de utilização 40 - 41
ES	VENDA ESMARCH REUTILIZABLE Instrucciones de uso 16 - 17	RO	BANDAJ ESMARCH REUTILIZABIL Instrucțiuni de utilizare 42 - 43
ET	TAASKASUTATAV ESMARCHI BANDAAŽ Kasutamisujuhüis..... 18 - 19	RU	ЖГУТ ЭСМАРХА МНОГОРАЗОВЫЙ Инструкция по применению 44 - 45
FI	ESMARCH-SIDE, UUDELLEENKÄYTETTÄVÄ Käyttöohje 20 - 21	SK	OVÍNADLO ESMARCH OPAKOVANE POUŽITEĽNÉ Návod na použitie 46 - 47
FR	BANDE D'ESMARCH RÉUTILISABLE Manuel d'utilisation 22 - 23	SL	ESMARCHOVA PREVEZA, ZA PONOVRNO UPORABO Navodila za uporabo 48 - 49
HR	ESMARCHOV ZAVOJ ZA VIŠEKRAATNU UPORABU Upute za uporabu 24 - 25	SV	ESMARKBINDA, ÅTERANVÄNDBAR Bruksanvisning..... 50 - 51
HU	ESMARCH BANDAGE, ÚJRAFELHASZNÁLHATÓ Használati utasítás 26 - 27	TR	ESMARCH BANDAJ TEKRAR KULLANILABİLİR Kullanım kılavuzu 52 - 53
IT	BENDAGGIO ESMARCH RIUTILIZZABILE Istruzioni per l'uso..... 28 - 29		SYMBOL DESCRIPTION 54 - 56

Deutsch

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English

INTENDED USE

The Esmarch Bandage is used to exsanguinate a patient's upper or lower extremity before using a Tourniquet Cuff.

Clinical benefit: The Esmarch Bandage creates a bloodless field in combination with a Tourniquet Cuff.

Patient target group: Patients who require a surgery at upper or lower extremities using a Tourniquet Cuff.

Place of use: Rooms which are suitable for surgical interventions.

INDICATIONS

- Surgeries in bloodless field

No other indications are known.

CONTRAINDICATIONS

- Open fractures of the extremities
- Severe crushing injuries
- Severe hypertension
- Skin transplant - impaired circulation (e.g. peripheral arterial disease)

No other contraindications are known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by a physician, or by medically trained personnel working under the instruction of a physician.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU Member State (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Maximum application time: 15 min
- No changes may be made to the product.
- The product must not come into contact with sharp or pointed objects.
- The reusable product is supplied non-sterile and must be cleaned, disinfected and sterilised prior to every use.

PREPARATION PRIOR TO USE

- ▶ Select suitable Esmarch Bandage.

VISUAL INSPECTION

- ▶ Check the product for damage.

A faulty product must be disposed of (see "Disposal").

USE

Exsanguination is supported by elevation of the extremity.

- ▶ Wrap the Esmarch Bandage around the patient's extremity, from distal to proximal up to the Tourniquet Cuff (Fig. 1).
- ▶ Inflate the Tourniquet Cuff to occlude blood flow.
- ▶ Remove the Esmarch Bandage.

REPROCESSING

(CLEANING, DISINFECTION, STERILISATION)

GENERAL INFORMATION

- The manufacturer is not responsible for damages resulting from improper reprocessing.
- The user is obliged to validate their procedures or the devices and accessories and to comply with the validated parameters during every reprocessing.
- It is recommended to use an automated procedure due to a significantly higher efficiency of cleaning and disinfection.
- The effectiveness has been verified by an independent and accredited test laboratory.
- To achieve effective reprocessing, heavy soiling must not be allowed to dry on the product and must be removed immediately after use.

CLEANING / DISINFECTION

Automated cleaning / disinfection

Notes regarding the use of washer disinfectors

- ▶ Use de-ionised water.
- ▶ Use a washer disinfectant in compliance with ISO 15883. For the validation a Miele G7882 equipped with cabinet E435/3 was used.

Notes regarding the use of cleaning agents and disinfectants

- When alkaline cleaning agents are used, the agents must be subsequently neutralised.

CAUTION

Do not use any highly alkaline cleaning agents, as this may cause foaming.

- Do not use drying agents.

Procedure validated by the manufacturer:

1. Unroll and spread out the product.
2. Position the product in a suitable basket in the injector carriage. Secure product with a cover net.
3. Start program with the following parameters:
 - a. Pre-rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 1 min.
 - b. Cleaning at 55 °C (131 °F) with an exposure time of 5 min using de-ionised water and "Sekumatic® ProClean" cleaning agent (dose: 0.5% (5 ml/l)).
 - c. Neutralisation with "Sekumatic® FNZ" at 20 °C (68 °F) with an exposure time of 2 min (dose: 0.1% (1 ml/l)).
 - d. Rinse with de-ionised water (max. 100 CFU/ml) at 20 °C (68 °F) and an exposure time of 2 min.
 - e. Thermal disinfection at 93 °C (199 °F) for 5 min with de-ionised water.
 - f. Dry at 100 °C (212 °F).
4. Check for visible contamination. If necessary, repeat the reprocessing.
5. Check the product according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

Manual cleaning / disinfection

Notes regarding the use of cleaning agents and disinfectants

- Prepare the cleaning and disinfection solutions prior to each manual cycle.
- Use de-ionised water.

Procedure validated by the manufacturer:

1. Prepare a 2% (30 ml/l) cleaning and disinfection solution of "Sekusept® Aktiv" with de-ionised water at 20 °C (68 °F). Stir several times in the next 15 min. After 15 min the cleaning and disinfection solution is ready for use.
2. Unroll and spread out the product.
3. Clean the product with soft sponges in the cleaning and disinfection solution, the products must be completely immersed for at least 20 sec. Hard brushes and other materials which may damage the surface of the products must not be used.
4. Place the product into the cleaning and disinfection solution and leave for 15 min. Ensure all surfaces are covered with solution.
5. Remove cleaning and disinfection solution by rinsing thoroughly with de-ionised water (max. 100 CFU/ml). Remaining residues can shorten the product shelf life or can lead to material damage.
6. Dry the product. Avoid any accumulation of water.
7. Check for visible contamination. If necessary, repeat the reprocessing.
8. Check the product according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

CHECKS

After disinfection, the product must be inspected for damages (cracks, breakage, etc.).

A faulty product must be disposed of (see "Disposal").

PACKAGING

- ▶ Carefully roll up the product for sterilisation with a gauze bandage (Fig. 2).
- ▶ After the cleaning and disinfection procedure, the product must be packed into a sterile barrier system which is suitable for steam sterilisation. The sterile barrier system must be in compliance with ISO 11607-1. For the validation a single pouch sterile barrier system was used.

STERILISATION

Notes regarding sterilisation:

- Product or sterilisation packages must be protected from mechanical damage.
- Use a steriliser that meets the requirements of DIN EN 285 or DIN EN 13060.

Sterilisation validated by the manufacturer

- The product must be steam sterilised by a fractionated vacuum process.
- Exposure time is 5 min with a sterilisation temperature of 134 °C (273 °F).
- Drying time is 10 min.
- For the validation a Systec HX-320 autoclave was used.

REUSABILITY

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use. Considering the maximum shelf life of 5 years, reusable products can be reprocessed up to 100 times if appropriately used and under adherence of the cleaning, disinfection and sterilisation instructions.

Every further reuse is the responsibility of the user (see section "Checks").

WARNING

By using the products on patients with suspected prion disease there may be a high risk of transmission. In such case the physician decides if the product is either disposed of (see section "Disposal") or reprocessed according to national directives.

SHELF LIFE

Product shelf life is 5 years and can be reprocessed within the shelf life up to 100 times.

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must have run through the complete reprocessing process to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

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