**Cuff Cover, non-sterile**




**CE**

**Cuff Cover, sterile**

**CE** 0123

---

<b>DE</b>	<b>Schutzhülle</b> Gebrauchsanweisung .....	4 - 5	<b>NL</b>	<b>Beschermhoes</b> Gebruiksaanwijzing .....	12 - 13
<b>EN</b>	<b>Cuff Cover</b> Instructions for use .....	5 - 6	<b>PT</b>	<b>Manga de proteção</b> Instruções de utilização .....	13 - 14
<b>CS</b>	<b>Ochranný návlek</b> Návod k použití .....	6 - 7	<b>RU</b>	<b>Защитная оболочка</b> Инструкция по применению ..	14 - 15
<b>DA</b>	<b>Beskyttelseshylster</b> Bruganvisning .....	7 - 8	<b>SV</b>	<b>Manschettskydd</b> Bruksanvisning .....	15 - 16
<b>ES</b>	<b>Funda de protección</b> Instrucciones de uso .....	8 - 9	<b>TR</b>	<b>Kaf kilifi</b> Kullanım kılavuzu .....	16 - 17
<b>FR</b>	<b>Housse de brassard</b> Manuel d'utilisation .....	9 - 10		<b>Symbol Description</b> .....	18 - 20
<b>IT</b>	<b>Copertura bracciale tourniquet</b> Istruzioni per l'uso .....	11 - 12			

# Deutsch

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- LDPE (Low Density Polyethylen)

#### SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before use, the product must be inspected visually (see section "Visual inspection"). A faulty product must not be used.
- It is not permitted to make any changes to the product.

- The product is suitable for MRI.

- The product is intended for single use and must not be reused and/or reprocessed. The function of the product is impaired by reprocessing. Any reuse entails a potential infection hazard.

#### REF 20-11-111, 20-11-333, 20-11-222

- The product is not sterile.

- The product must not be used if the expiration date has expired.

#### REF 20-11-111-1, 20-11-333-1, 20-11-222-1

- The product is sterile (ethylene oxide).

- The product must not be used if the packaging is damaged or the expiration date has expired.

#### PREPARATION PRIOR TO USE

- ▶ Choose the correct size for the respective extremity (see section "Product specifications").

#### VISUAL INSPECTION

- ▶ Check the packaging for damage.
- ▶ Check the product for damage and loose particles.

A faulty product must be disposed of (see section "Disposal").

#### APPLICATION

- ▶ Slip the Cuff Cover over the extremity (Fig. 1).
- ▶ Place the Tourniquet Cuff (Fig. 2).
- ▶ Fold the distal part of the Cuff Cover over the Tourniquet Cuff. Seal the Cuff Cover fluid-tight at the distal edge to avoid skin damage (Fig. 3).
- ▶ After application, deflate the Tourniquet Cuff.
- ▶ Then immediately remove the Tourniquet Cuff with the Cuff Cover from the extremity to prevent the risk of venous congestion.

#### SHELF LIFE

The shelf life of the product is 5 years from the date of manu-

## English

#### INTENDED USE

The Cuff Cover protects Tourniquet Cuffs from soiling.

Clinical benefit: Hygienic protection of Tourniquet Cuffs.

Patient target group: Patients who require surgery at the upper or lower extremities in conjunction with a Tourniquet Cuff.

Place of use: Rooms which are suitable for surgical interventions.

#### INDICATIONS

- Surgeries in bloodless field

No other indications are known.

#### CONTRAINDICATIONS

- Open fractures of the extremities

No other contraindications are known.

facture.

Use-by date: see product label

## STORAGE AND TRANSPORT CONDITIONS

### CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

## DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

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