Manual Tourniquet

Intended to be left blank.

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English

INTENDED USE

The Manual Tourniquet is a manually operated tourniquet. It regulates the pressure of a Tourniquet Cuff which temporarily occludes the blood flow of a patient's upper or lower extremity in order to create a bloodless field. Clinical benefit: maintaining a bloodless field and reducing blood loss during procedures on extremities in conjunction with Tourniquet Cuffs.

Patient target group: patients who require surgery at the upper or lower extremities.

Place of use: rooms which are suitable for surgical interventions.

INDICATION / CONTRAINDICATION

Indications and contraindications depend on the selected Tourniquet Cuff.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- The product may only be used by a physician or by medically trained personnel working under the instruction of a physician.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before every use, the product must be inspected visually for damage (cracks, breakage, etc.) and subjected to a functional check (see "Visual and functional check"). A faulty product must not be used.
- It is not permitted to make any changes to the product.
- Any kinking / clamping of the connecting tubing between the Manual Tourniquet and the Tourniquet Cuff must be avoided.
- The cuff pressure must comply with the pressure specifications given in the Instructions for Use of the Tourniquet Cuffs.
- If tourniquet cuffs from other manufacturers are used, compatibility must be ensured.
- Prior to every use, the product must be disinfected (see "Wipe disinfection").
- The product is not sterile.
- The product is not suitable for MRI.

PRODUCT DESCRIPTION

- Manometer 0 - 600 mmHg
- Release screw for deflating the Tourniquet Cuff
- Blue / red coil tubing with female positive locking connectors
- Hand inflator for inflating the Tourniquet Cuff
- Hand pump for inflating the Tourniquet Cuff

VISUAL AND FUNCTIONAL CHECK

Prior to each use of the Manual Tourniquet, the functionality and tightness of the entire tourniquet system must be checked.

- ▶ Check the product for damage and loose particles.
- ▶ If present, remove loose particles.
- ▶ Roll-up the Tourniquet Cuff tightly to enable a counter-pressure to build up when the cuff is inflated.
- ▶ Close the release screw tightly.
- ▶ Inflate Tourniquet Cuff to a pressure of 600 mmHg.
- ▶ Reduce the pressure to 500 mmHg.
- ▶ Check the pressure on the analogue pressure gauge.

CAUTION

- The pressure drop must not be greater than 20 mmHg in 5 minutes.
- If the pressure drops, immediately label the Manual Tourniquet as defective and contact the manufacturer.

A faulty product must be disposed of (see "Disposal").

APPLYING THE TOURNIQUET CUFF

See corresponding Instructions for Use for Tourniquet Cuffs.

APPLICATION

WARNING

- While inflating the Tourniquet Cuff, perform tactile inspection with your hand to check that the cuff chamber is being properly inflated.
- During surgery, the cuff pressure of the tourniquet system must be continuously checked on the pressure gauge. The pressure gauge always displays the actual cuff pressure. Any pressure drop can be read on the pressure gauge.

- Connect the Tourniquet Cuff with the Manual Tourniquet using the positive locking connector.
- Then inflate the Tourniquet Cuff to the lowest required cuff pressure.
- To compensate for subsequent material expansion, first generate a pressure 50 mmHg higher and then correct it to the desired pressure using the release screw (e.g. for a desired pressure of 300 mmHg, first generate a pressure of 350 mmHg and then reduce it to 300 mmHg using the release screw).
- The actual cuff pressure is displayed continuously on the pressure gauge and must be monitored permanently.
- After completion of the surgery, deflate the Tourniquet Cuff using the release screw.

WIPE DISINFECTION

CAUTION

- The product must not be mechanically or manually reprocessed or sterilised.
- Do not immerse the product in liquids.
- Do not clean the pressure gauge with aggressive cleaning agents.
- In case of unsuccessful wipe disinfection or heavy contamination, the product must be disposed of as specified in "Disposal".

Wipe disinfection must be performed with commercially available surface disinfectants based on alcohol or QACs (quaternary ammonium compounds). When choosing a product for disinfection, a disinfectant with an appropriate spectrum of action must be used: bactericidal, yeasticidal and virucidal. After wipe disinfection, inspect the product for visible contamination. If necessary, repeat wipe disinfection. After wipe disinfection, check the function of the product (see "Visual and functional check").

METROLOGICAL CHECK

The metrological check may only be performed by qualified personnel in possession of suitable measuring equipment.

- The display accuracy ± 6 cmH₂O of the product must be verified every 24 months with a metrological check.
- Inflate the product to 100, 200, 300, 400, 500 and 600 mmHg with the hand inflator (d) or hand pump (e).
- The values must be within the indicated display accuracy.

If this is not the case, the manufacturer must be contacted.

If the operator does not have the required measuring equipment, the metrological check can be carried out by the manufacturer. In this case, send the product to the manufacturer.

If there are any signs before expiry of the 24 months indicating that the product fails to adhere to the display accuracy or that the metrological characteristics of the product have been affected, a metrological check must be performed immediately.

Conversion of pressure units: 1 hPa = 1.02 cmH₂O = 0.75 mmHg

SHELF LIFE

The shelf life of the product is 5 years.

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must be cleaned and disinfected to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

PRODUCT SPECIFICATIONS

REF	20-18-601	20-18-602	20-19-601	20-19-602
Display accuracy	± 6 mmHg			
Pressure range	0 - 600 mmHg			
Colour of coil connecting tubing	blue	red	blue	red
Weight (incl. packaging)	370 g		560 g	
Packaging unit	1	1	1	1

MATERIAL DATA

REF 20-18-601 20-18-602	Tubing: PUR (polyurethane) Inflation bulb: PVC (polyvinylchloride) Connector: POM (polyoxymethylene)
REF 20-19-601 20-19-602	Pump: aluminium Tubing: PUR (polyurethane) Connector: POM (polyoxymethylene) Hook and cylinder screw: stainless steel

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