Tourniquet Drape, non-sterile

CE

Tourniquet Drape, sterile

CE 0123

Intended to be left blank.

DE	Manschetten Drape Gebrauchsanweisung 4 - 5	FR	Drap pour brassard Manuel d'utilisation 7 - 8
EN	Tourniquet Drape Instructions for use 5 - 6	IT	Telino bracciale tourniquet Istruzioni per l'uso 9 - 10
ES	Cubierta para manguito Instrucciones de uso 6 - 7		Symbol Description 11 - 12

Deutsch

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English

INTENDED USE

The Tourniquet Drape creates a fluid-tight barrier between the surgical field and the Tourniquet Cuff.

Clinical benefit: Protection against accumulation of fluids underneath the Tourniquet Cuff.

Patient target group: Patients who require surgery at the upper or lower extremities in conjunction with a Tourniquet Cuff.

Place of use: Rooms which are suitable for surgical interventions.

INDICATIONS

- Surgeries in bloodless field
- No other indications are known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel.

- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.

- Before use, the product must be inspected visually (see section "Visual inspection"). A faulty product must not be used.

- It is not permitted to make any changes to the product.

- The product is intended for single use and must not be reused and/or reprocessed. The function of the product is impaired by reprocessing. Any reuse entails a potential infection hazard.

REF 24-00-100, 24-00-200, 24-00-300

- The product is not sterile.

- The product must not be used if the expiration date has expired.

REF 24-00-100-1, 24-00-200-1, 24-00-300-1

- The product is sterile (ethylene oxide).

- The product must not be used if the packaging is damaged or the expiration date has expired.

PREPARATION PRIOR TO USE

CAUTION

A too large Tourniquet Drape does not guarantee an adequate seal and bears the risk of fluid accumulation under the Tourniquet Cuff.

- Choose the correct size for the respective extremity (see section "Product specifications").

VISUAL INSPECTION

- Check the packaging for damage.
 - Check the product for damage and loose particles.
- A faulty product must be disposed of (see "Disposal").

APPLYING THE TOURNIQUET CUFF

See corresponding Instructions for Use for Tourniquet Cuffs.

APPLICATION

- Slip the Tourniquet Drape over the extremity to the distal edge of the Tourniquet Cuff (see cover image).
- Place the Tourniquet Drape proximally so that it does not protrude into the sterile field.

Ensure a good seal between the Tourniquet Cuff and the surgical field.

- ▶ After application, deflate the Tourniquet Cuff. Then immediately remove the Tourniquet Drape and the Tourniquet Cuff from the extremity to prevent the risk of venous congestion.

SHELF LIFE

Use-by date: see product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

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