Hand Inflator for Roll-On Cuff

Intended to be left blank.

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INTENDED USE

The Roll-On Cuff is inflated using a hand inflator.

Clinical benefit: enables the use of a roll-on cuff.

Patient target group: patients who require surgery at the upper or lower extremities.

Place of use: rooms which are suitable for surgical interventions.

INDICATIONS / CONTRAINDICATIONS

Indications and contraindications depend on the selected Roll-On Cuff.

SAFETY INSTRUCTIONS

• Read and follow the instructions for use carefully before using the product and keep them for future reference.

 • This product must only be used by medically trained personnel.

• The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.

• Before every use, the product must be inspected visually for damage (cracks, breakage, etc.) and subjected to a functional check (see "Visual and functional check"). A faulty product must not be used.

• It is not permitted to make any changes to the product.

• Any kinking / clamping of the connecting tubing between the Hand Inflator and the Roll-On Cuff must be avoided.

• The cuff pressure must comply with the pressure specifications given in the Instructions for Use of the Roll-On Cuff.

• Use original VBM individual parts only. Using individual parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.

• Prior to every use, the product must be disinfected (see "Wipe disinfection").

• The product is not sterile.

• The product is not suitable for MRI.

PRODUCT DESCRIPTION

a - Manometer 0 - 200 mmHg

b - Release screw for deflating the Roll-On Cuff

c - 1 m connecting tubing with quick-release coupling and hollow needle

d - Hand Inflator for inflating the Roll-On Cuff

VISUAL AND FUNCTIONAL CHECK

Prior to each use of the hand inflator, the functionality and tightness of the entire system must be checked.

▶ Check the product for damage and loose particles.

▶ If present, remove loose particles.

CAUTION

• Insert the needle vertically into the valve of the Roll-On Cuff. Insertion at an angle may damage the valve irreparably and thus cause it to leak.

▶ Close the release screw tightly.

▶ Inflate the Roll-On Cuff to a pressure of 150 mmHg.

▶ Reduce the pressure to 130 mmHg.

▶ Then check the pressure.

CAUTION

• The pressure drop must not be greater than 10 mmHg in 2 minutes.

• In case of a larger pressure drop, immediately label the device as defective and contact the manufacturer.

A faulty product must be disposed of (see "Disposal").

APPLICATION

See corresponding Instructions for Use for Roll-On Cuff.

WIPE DISINFECTION

CAUTION

• The product must not be mechanically or manually reprocessed or sterilised.

• Do not immerse the product in liquids.

• Do not clean the pressure gauge with aggressive cleaning agents.

• In case of unsuccessful wipe disinfection or heavy contamination, the product must be disposed of as specified in "Disposal".

Wipe disinfection must be performed with commercially available surface disinfectants based on alcohol or QACs (quaternary ammonium compounds). When choosing a product for disinfection, a disinfectant with appropriate ranges of action must be used: bactericidal, yeasticidal and virucidal. After wipe disinfection, inspect the product for visible contamination. If necessary, repeat wipe disinfection. After wipe disinfection, check the function of the product (see "Visual and functional check").

METROLOGICAL CHECK

The metrological check may only be performed by qualified personnel in possession of suitable measuring equipment.

- ▶ The display accuracy ± 6 mmHg of the product must be verified every 24 months with a metrological check.
- ▶ Inflate the product to 50, 100, 150 and 200 mmHg with the Hand Inflator (d).
- ▶ The values must be within the indicated display accuracy.

If this is not the case, the manufacturer must be contacted. If the operator does not have the required measuring equipment, the metrological check can be carried out by the manufacturer. In this case, send the product to the manufacturer.

If there are any signs before expiry of the 24 months indicating that the product fails to adhere to the display accuracy or that the metrological characteristics of the product have been affected, a metrological check must be performed immediately.

Conversion of pressure units:

$$1 \text{ hPa} = 1.02 \text{ cmH}_2\text{O} = 0.75 \text{ mmHg}$$

SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture.

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must be cleaned and disinfected to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

PRODUCT SPECIFICATIONS

REF	21-95-300
Display accuracy	± 6 mmHg
Pressure range	0 - 200 mmHg
Weight (incl. packaging)	380 g (0.84 lbs)
Packaging unit	1

MATERIAL DATA

REF 21-95-300	Inflation bulb: PVC (polyvinylchloride) Tubing from manometer to connector: PVC (polyvinyl chloride) Connector: POM (polyoxymethylene) O-ring: EPDM (ethylene propylene diene monomer) Tubing: silicone Hollow needle: stainless steel
------------------	---

ACCESSORIES

REF 21-95-302	Connecting tubing for Hand Inflator for Roll-On Cuff
------------------	--

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