



## Dispo Infusor

---

<b>DE</b>	<b>Dispo Infusor</b> Gebrauchsanweisung ..... 4 - 5	<b>LT</b>	<b>Dispo Infusor</b> Naudojimo instrukcija ..... 19 - 20
<b>EN</b>	<b>Dispo Infusor</b> Instructions for use ..... 5 - 6	<b>LV</b>	<b>Dispo Infusor</b> Lietošanas instrukcija ..... 21 - 22
<b>BG</b>	<b>Dispo Infusor</b> Инструкции за употреба ..... 6 - 7	<b>NL</b>	<b>Dispo Infusor</b> Gebruiksaanwijzing ..... 22 - 23
<b>CS</b>	<b>Dispo Infusor</b> Návod k použití ..... 7 - 8	<b>NO</b>	<b>Dispo Infusor</b> Bruksanvisning ..... 23 - 24
<b>DA</b>	<b>Dispo Infusor</b> Bruganvisning ..... 8 - 9	<b>PL</b>	<b>Dispo Infusor</b> Instrukcja użycia ..... 24 - 25
<b>EL</b>	<b>Dispo Infusor</b> Οδηγίες χρήσης ..... 10 - 11	<b>PT</b>	<b>Dispo Infusor</b> Instruções de utilização ..... 25 - 27
<b>ES</b>	<b>Dispo Infusor</b> Instrucciones de uso ..... 11 - 12	<b>RO</b>	<b>Dispo Infusor</b> Instrucțiuni de utilizare ..... 27 - 28
<b>ET</b>	<b>Dispo Infusor</b> Kasutamisujuhüis ..... 12 - 13	<b>RU</b>	<b>Dispo Infusor</b> Инструкция по использованию ..... 28 - 29
<b>FI</b>	<b>Dispo Infusor</b> Käyttöohje ..... 13 - 14	<b>SK</b>	<b>Dispo Infusor</b> Návod na použitie ..... 29 - 30
<b>FR</b>	<b>Dispo Infusor</b> Manuel d'utilisation ..... 14 - 16	<b>SL</b>	<b>Dispo Infusor</b> Navodila za uporabo ..... 30 - 31
<b>HR</b>	<b>Dispo Infusor</b> Upute za uporabu ..... 16 - 17	<b>SV</b>	<b>Dispo Infusor</b> Bruksanvisning ..... 32 - 33
<b>HU</b>	<b>Dispo Infusor</b> Használati utasítás ..... 17 - 18	<b>TR</b>	<b>Dispo Infusor</b> Kullanım kılavuzu ..... 33 - 34
<b>IT</b>	<b>Dispo Infusor</b> Istruzioni per l'uso ..... 18 - 19		<b>Symbol Description</b> ..... 35 - 36

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- Application time: maximum of 5 days, depending on the indication.
- It is not permitted to make any changes to the product.
- The product is not sterile.
- The product is not suitable for MRI.

- The product may be used repeatedly on a single patient. It must not be reprocessed or used on other patients. The function of the product is impaired by reprocessing. Any reuse entails a potential infection hazard.
- The product must not be used if the expiration date has expired.

## PRODUCT DESCRIPTION

- a - Hanger to hang up on the IV-stand
- b - Fleece disc
- c - Hook fastener for fixation of the fluid bag
- d - Netting
- e - Hook for fixation of the fluid bag
- f - Three-way stopcock
- g - Hand inflator

## VISUAL INSPECTION

- ▶ Check the product for damage.

A faulty product must be disposed of (see "Disposal").

## APPLICATION

### NOTE

- For optimal pressure transmission the fluid bag must be covered completely by the netting (d) of the Dispo Infusor.
- The material expands during inflation, which may initially cause a slight drop in pressure.
- To avoid damaging the Dispo Infusor, do not exceed the maximum pressure of 300 mmHg.
- The pressure drop due to emptying of the fluid bag must be monitored constantly and readjusted if necessary.

- ▶ Insert the fluid bag from below (Fig. 1).
- ▶ Only for 500 and 1000 ml: pass the hook fastener (c) through the eyelet of the fluid bag. Then attach the hook fastener to the fleece disc (b) according to the size of the fluid bag (Fig. 2).
- ▶ Only for 3000 ml: hang the fluid bag on the hook (e) (Fig. 3).
- ▶ Hang the Dispo Infusor on the IV-stand.
- ▶ Adjust the Dispo Infusor to the desired pressure using the hand inflator (g) and the three-way stopcock (f).

Connect the Dispo Infusor with Luer lock (REF 56-01-300) to the automatic pressure infusor for inflation.

### Function of the three-way stopcock (f) (Fig. 4):

- Inflation ①
- Deflation ②
- Airtight ③

## AFTER USE

- ▶ Deflate the Dispo Infusor via the three-way stopcock (f) or via the automatic pressure infusor (3000 ml only).
- ▶ Remove the fluid bag.

Before inserting the new fluid bag, the Dispo Infusor must be fully deflated.

## SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture.

Use-by date: See product label

## English

### INTENDED USE

To pressurise fluid bags by external pressure.

Clinical benefit:

- Rapid compensation of fluid/volume deficiency
- Irrigation of the operation field
- Support for continuous blood pressure monitoring

Patient target group: all patients with a medical indication.

Place of use: hospital and prehospital

### INDICATIONS

Pressure infusion, irrigation purposes, transfusion procedures and invasive pressure monitoring.

No other indications are known.

### CONTRAINDICATIONS

None known.

### SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before use, the product must be inspected visually (see Section "Visual inspection"). A faulty product must not be used.

## STORAGE AND TRANSPORT CONDITIONS

### CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

## DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

## PRODUCT SPECIFICATIONS

REF	Fluid bag	Model	Weight (incl. packaging)	Packaging unit
56-02-050	500 ml	with hand inflator, three-way stopcock and colour coded pressure indicator	840 g (30.0 oz)	5
56-02-100	1000 ml	with hand inflator, three-way stopcock and colour coded pressure indicator	870 g (30.7 oz)	5
56-01-300	3000 ml	with female Luer lock	780 g (27.5 oz)	5
56-02-300	3000 ml	with hand inflator, three-way stopcock and colour coded pressure indicator	1170 g (41.3 oz)	5

## MATERIAL DATA

REF	56-02-050 56-02-100	56-01-300	56-02-300
Fabric, piping	PES / cotton		
Connector	LDPE		
Loop	PET	-	-
Fleece disc, hook fastener	PA	-	-
Connector Luer lock, netting	PA		
Tubing	Silicone		
Hook	-	TPU	
Eyelet	-	Nickel-plated brass	
Inflation bulb	PVC	-	PVC
Connecting tubing	TPU	-	TPU
Three-way stopcock	PC / HDPE	-	PC / HDPE
Colour coded pressure indicator	ABS	-	ABS

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