





Clear Infusor

Intended to be left blank.

DE	Clear Infusor Gebrauchsanweisung	5 - 6	LT	Clear Infusor Naudojimo instrukcija	24 - 25
EN	Clear Infusor Instructions for use	6 - 7	LV	Clear Infusor Lietošanas instrukcija	26 - 27
BG	Clear Infusor Инструкции за употреба	8 - 9	NL	Clear Infusor Gebruiksaanwijzing	27 - 28
CS	Clear Infusor Návod k použití	9 - 10	NO	Clear Infusor Bruksanvisning	29 - 30
DA	Clear Infusor Brugsanvisning	11 - 12	PL	Clear Infusor Instrukcja użycia	30 - 31
EL	Clear Infusor Οδηγίες χρήσης	12 - 13	PT	Clear Infusor Instruções de utilização	32 - 33
ES	Clear Infusor Instrucciones de uso	14 - 15	RO	Clear Infusor Instrucțiuni de utilizare	33 - 34
ET	Clear Infusor Kasutamisinjuhised	15 - 16	RU	Clear Infusor Инструкция по использованию	35 - 36
FI	Clear Infusor Käyttöohje	17 - 18	SK	Clear Infusor Návod na použitie	36 - 37
FR	Clear Infusor Manuel d'utilisation	18 - 19	SL	Clear Infusor Navodila za uporabo	38 - 39
HR	Clear Infusor Upute za uporabu	20 - 21	SV	Clear Infusor Bruksanvisning	39 - 40
HU	Clear Infusor Használati utasítás	21 - 22	TR	Clear Infusor Kullanım kılavuzu	41 - 42
IT	Clear Infusor Istruzioni per l'uso	23 - 24		Symbol Description	43 - 44

Deutsch

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English

INTENDED USE

To pressurise fluid bags by external pressure.

Clinical benefit:

- Rapid compensation of fluid/volume deficiency
 - Irrigation of the operation field
 - Support for continuous blood pressure monitoring
- Patient target group: all patients with a medical indication.
Place of use: hospital and prehospital

INDICATIONS

Pressure infusion, irrigation purposes, transfusion procedures and invasive pressure monitoring.

No other indications are known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.
- This product must only be used by medically trained personnel.

- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before every use, the product must be inspected visually for damage (cracks, breakage, etc.) (see Section "Visual Inspection"). A faulty product must not be used.
- Application time: maximum of 5 days, depending on the indication.
- It is not permitted to make any changes to the product.
- Use original VBM individual parts only. Using individual parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.
- The product must be disinfected prior to use on any new patient (see Section "Wipe disinfection").
- The product is not sterile.
- The product is not suitable for MRI.

- The product must not be used if the expiration date has expired.

PRODUCT DESCRIPTION

- Hanger to hang up on the IV-stand
- Hook for fixation of the fluid bag
- Three-way stopcock
- Hand inflator

VISUAL INSPECTION

- ▶ Check the product for damage.

A faulty product must be disposed of (see "Disposal").

APPLICATION

NOTE

- For optimal pressure transmission the Clear Infuser must be wrapped completely around the fluid bag.
- The material expands during inflation, which may initially cause a slight drop in pressure.
- To avoid damaging the Clear Infuser, do not exceed the maximum pressure of 300 mmHg.
- The pressure drop due to emptying of the fluid bag must be monitored constantly and readjusted if necessary.

- ▶ Open the hook and loop fastener and unfold the Clear Infuser (Fig. 1).
- ▶ Hang the fluid bag on the hook (b) (Fig. 2).

- ▶ Wrap the Clear Infusor tightly around the fluid bag and close the hook and loop fastener (Fig. 3).
- ▶ Hang the Clear Infusor on the IV-stand.
- ▶ Adjust the Clear Infusor to the desired pressure using the hand inflator (d) and the three-way stopcock (c) (Fig. 4).

Connect the Clear Infusor with Luer lock (REF 57-01-300) to the automatic pressure infusor for inflation.

Function of the three-way stopcock (c) (Fig. 5):

- Inflation ①
- Deflation ②
- Airtight ③

AFTER USE

- ▶ Deflate the Clear Infusor via the three-way stopcock (c) or via the automatic pressure infusor (3000 ml only).
- ▶ Remove the fluid bag.

WIPE DISINFECTION

WARNING

The risk of contamination increases if the product is heavily soiled or if a bioburden is suspected. In this case the physician decides if the product is either disposed of (see Section "Disposal") or disinfected according to national directives.

CAUTION

- The product must not be mechanically or manually reprocessed or sterilised.
- Do not immerse the product/device in liquids.
- In case of unsuccessful wipe disinfection or heavy contamination, the product must be disposed of as specified in "Disposal".

Wipe disinfection must be performed with commercially available surface disinfectants based on alcohol, QACs (quaternary ammonium compounds) or hydrogen peroxide. When choosing a product for disinfection, a disinfectant with appropriate ranges of action must be used: bactericidal, yeasticidal and virucidal. After wipe disinfection, inspect the product for visible contamination. If necessary, repeat wipe disinfection. After the wipe disinfection, check the product (see "Visual inspection").

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must be cleaned and disinfected to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture and can be reprocessed 100 times within the period of the shelf life. Use-by date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

PRODUCT SPECIFICATIONS

REF	Fluid bag	Model	Weight (incl. packaging)	Packaging unit
57-02-050	500 ml	with hand inflator, three-way stopcock and colour coded pressure indicator	430 g (15.2 oz)	1
57-08-050	500 ml	with hand inflator, three-way stopcock and manometer	480 g (16.9 oz)	1

REF	Fluid bag	Model	Weight (incl. packaging)	Packaging unit
57-02-100	1000 ml	with hand inflator, three-way stopcock and colour coded pressure indicator	560 g (19.8 oz)	1
57-08-100	1000 ml	with hand inflator, three-way stopcock and manometer	610 g (21.5 oz)	1
57-02-300	3000 ml	with hand inflator, three-way stopcock and colour coded pressure indicator	820 g (28.9 oz)	1
57-08-300	3000 ml	with hand inflator, three-way stopcock and manometer	870 g (30.7 oz)	1
57-01-300	3000 ml	with female Luer lock	750 g (26.5 oz)	1

MATERIAL DATA

REF	57-02-050 57-02-100 57-02-300	57-08-050 57-08-100 57-08-300	57-01-300
Bladder, flap, hook, connector	TPU		
Tubing	Silicone		
Connector Luer lock	PA		
Hook and loop	PA / PU		
Inflation bulb	PVC	-	-
Connecting tubing	TPU	-	-
Three-way stopcock	PC / HDPE		
Colour coded pressure indicator	ABS	-	-
Manometer glass	-	PC	-
Manometer protective ring	-	EPDM	-

REPLACEMENT PARTS (FIG. 6)

REF	Description	Weight (incl. packaging)	Packaging unit
20-17-301	Manometer 300 mmHg for fixation to the pressure infusion cuff	75 g (2.7 oz)	1
51-02-666	Hand inflator with three-way stopcock and colour coded pressure indicator	82 g (2.9 oz)	1
51-03-666	Hand inflator with three-way stopcock	61 g (2.2 oz)	1

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