







Laryngeal Tube LTS II

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ISO connector
- colour-coded, for easy identification of different sizes

Valve opener

Inflation line with pilot balloon and Luer connector
- inflates and deflates both cuffs

Drain tube
- for gastric and suction catheter

Teeth marks
- for easy orientation and repositioning


Proximal cuff
- stabilises the tube and seals naso- and oropharynx

Ventilation

Ventilation orifices
- lie in front of the larynx

Pressure release

Distal cuff
- seals oesophagus
- reduces possibility of gastric insufflation

INTENDED USE

The LTS II is a supraglottic airway with drain tube.

INDICATIONS

- As alternative adjunct to temporarily secure the airway during difficult airway management and cardio-pulmonary resuscitation in clinical as well as preclinical environments (current guidelines of the European Resuscitation Council, American Heart Association and American Society of Anesthesiologists).
- For mechanical or spontaneous ventilation of patients during short anaesthesia procedures.

CONTRAINDICATIONS

- Complete obstruction of the upper airway.
 - In elective use: patients who have not fasted as well as patients with retained gastric contents.
- No other contraindications are known.

SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- The product may only be used by specialist medical personnel trained in the use of the product.
- During MRI treatment, the pilot balloon of the tube can cause artefacts around the area where it is positioned. A non-clinical test verified that the tube can be used in an MRI environment.
- The product is suitable for conducting oxygen, air and anaesthesia gases. The anaesthesia gases used must not interact with the product materials (see material data).
- Laser or electrocaustic equipment may cause the product to be flammable.
- The tube is supplied clean packed, ready for use and must be cleaned, disinfected and if necessary sterilised after first use and each subsequent use.
- The syringe is for single use and must not be reprocessed. Reprocessing leads to impairment of the function of the syringe. Any reuse entails a potential risk of cross-contamination.
- After disinfection and prior to each use, perform visual check for damages (cracks, ruptures etc.) and function check according to chapter "Control".

USE

- ▶ Choose the correct tube size (see table).
- Have a spare tube ready and prepared for immediate use.

- ▶ Deflate the cuffs completely (picture 1).
- ▶ Lubricate tip and posterior wall of the tube with water-soluble lubricant.
- ▶ Pre-oxygenation is advisable in all cases.
- ▶ Head position

The ideal head position for insertion is the "sniffing position". However, the tube's short length and its shape also allow it to be inserted with the head in a neutral position. For children, it is recommended (due to the anatomy) to have the head slightly extended with a chin lift manoeuvre.

Insertion:

CAUTION

Achieve an appropriate depth of anaesthesia before inserting the tube!

- ▶ Hold the tube like a pen in the area with the black lines (teeth marks) (picture 2).
- ▶ Open the mouth with the free hand (picture 3) and if possible, lift the chin. The tongue must not fold back during the downward movement of the tube.

CAUTION

Never use force to insert the tube! During insertion, the cuffs must not be damaged.

- ▶ Insert the tube either in the midline of the mouth or into the corner of the mouth while rotating the tube and advance to the upper teeth mark (picture 4).
- ▶ Use a syringe to inflate the cuffs with max. recommended volume (see table) (picture 5).

Thanks to the high-volume low-pressure cuffs, a seal can already be achieved with low cuff pressures (< 60 cmH₂O).

- ▶ Connect the tube to a ventilation system (picture 6).

CAUTION

- ▶ Ventilation must be checked each time the patient is repositioned.

Oxygen, air or anaesthesia gases can change the cuff volume or cuff pressure.

Excessive cuff pressure can lead to swelling of the tongue.

Reduce the cuff pressure to the Airway Leak Pressure and check regularly.

- ▶ Use an appropriate procedure to check ventilation and exclude gastric insufflation.
- ▶ If no sufficient ventilation can be verified, withdraw the tube with inflated cuffs between the teeth marks until ventilation is easy and free flowing (large tidal volume with minimal airway pressure).
- ▶ Check the cuff pressure again and adjust to the Airway Leak Pressure.
- ▶ Fix the tube.

From size 2, insertion with the optionally enclosed biteblock is recommended (picture 7).

The drain tube (see table) allows the insertion of a gastric tube or suction catheter.

Removal:

- ▶ Deflate the cuffs completely (picture 8).

Once the tube is in the correct position, it is well tolerated until the return of the protective reflexes.

- ▶ Remove the tube in an area where suction equipment and the ability for rapid intubation are present.

REPROCESSING (CLEANING, DISINFECTION, STERILISATION)

GENERAL INFORMATION

- Manufacturer is not responsible for damages resulting from improper reprocessing.
- It rests with the user to validate their procedures, including devices and accessories, and to comply with the validated parameters during every reprocessing.
- It is recommended to apply an automatic procedure due to a significantly higher efficiency of cleaning and disinfection.
- The effectiveness has been verified by an independent and accredited test laboratory.
- To achieve a more effective reprocessing rough contamination must not get dry on the product and must be removed immediately after use.
- Before cleaning and disinfection don't connect the red valve opener to the valve of the inflation line (only used for sterilization).
- Complete the Record Card each time the tube is reprocessed.
- Do not bring silicone products in contact with oil and grease.

Automatic Cleaning / Disinfection

Notes regarding the use of washer disinfectors

- Use de-ionised water.

Notes regarding the use of cleaning and disinfectants

- By the use of alkaline cleaning agents, neutralisation must be performed.
- Don't use drying agents.

Procedure validated by the manufacturer:

1. Connect ventilation lumen and drain tube to the appropriate adapters of the injector carriage. Place the biteblock in a suitable mesh tray into the injector carriage and secure with a cover net. **Ensure effective rinsing of all hollow spaces.**
2. Start program with following parameters:
 - a. Rinse with de-ionised water at 20 °C (68 °F) with an exposure time of 1 min.
 - b. Cleaning at 55 °C (131 °F) with an exposure time of 5 min. using de-ionised water and the cleaner „Sekumatic® ProClean“ (Dose: 0,5 % (5 ml/l))
 - c. Neutralisation with „Sekumatic® FNZ“ at 20 °C (68 °F), with an exposure time of 2 min. (Dose: 0,1 % (1 ml/l))
 - d. Rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 2 min.
 - e. Thermal disinfection at 93 °C (199 °F) for 5 min. with de-ionised water.
 - f. Drying at max. 100 °C (212 °F)
3. Examine the products for visible contamination. If necessary repeat the reprocessing.
4. Check the products acc. to the chapter "Control" and prepare them for sterilisation (see chapter "Packaging").

Manual Cleaning / Disinfection

Notes regarding the use of cleaning and disinfectants

- Prepare the cleaning- and disinfection solution before each manual cycle.
- Use de-ionised water.

Procedure validated by the manufacturer:

1. Prepare a 2 % (30 ml/l) cleaning- and disinfection solution with Sekusept® Aktiv and de-ionised water at 20 °C (68 °F). The cleaning- and disinfection solution is ready for use after 15 min.
2. Clean the products in the cleaning- and disinfection solution by using soft sponges. Areas difficult to access must be cleaned with soft brushes. Hard brushes and other materials which can damage the surface of the products must not be used.
3. Put the products into the cleaning- and disinfection solution with an exposure time of 15 min. Ensure complete wetting of all hollow spaces.
4. Remove the cleaning- and disinfecting solution removed by thoroughly rinsing with de-ionised water. Remaining residues can shorten the products shelf life or can lead to material damage.
5. Dry the products. Accumulation of water must be avoided.
6. Check for visible contamination. If necessary repeat the reprocessing.
7. Check the products acc. to the chapter "Control" and prepare them for sterilisation (see chapter "Packaging").

CONTROL

After disinfection and prior to each use the products must be checked acc. to the following visible and functional control:

- Check the products for damage (cracks, ruptures etc.) and loose particles.
- Check the inside of the tube for sufficient patency.
- Inflate the cuffs and check for leaks.

Damaged products must not be re-used.

PACKAGING

- Connect the red valve opener firmly to the inflation line which opens the valve of the pilot balloon for a free gas exchange. Avoid kinking of the inflation line.
- Alternatively the cuffs can be evacuated prior to packaging. Sterilisation must then follow immediately.

Note:

High temperatures and negative pressure in the autoclave lead to distention of the residual air in the cuff which can result in irreparable damage to the cuffs.

After the cleaning- and disinfection procedure the products must be packed into a package, which is suitable for steam sterilisation (acc. to the currently valid issue of DIN EN 868 and DIN EN ISO 11607).

STERILISATION

Notes regarding sterilisation:

- Protect the products and sterilisation pouch against mechanical damage.

Sterilisation validated by the manufacturer

- The products must be steam sterilised by a fractional vacuum process.
- Exposure time is 5 min. with a sterilisation temperature of 134 °C (273 °F).

REUSABILITY

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use. Considering the maximum shelf life of 5 years, reusable products can be reprocessed up to 50 times if appropriately used and under adherence of the cleaning-, disinfection and sterilisation instructions.

Every further reuse is in the responsibility of the user (see chapter "Control")

WARNING

By using the products on patients with suspected prion disease there might be a high risk of transmission. In such case the physician decides if the product is either disposed (see chapter "Disposal") or reprocessed acc. to the national directives.

PRODUCT SPECIFICATIONS

Tube size	#0	#1	#2	#2.5	#3	#4	#5
Weight/height	< 5 kg	5-12 kg	12-25 kg	125-150 cm	< 155 cm	155-180 cm	> 180 cm
Colour coding ISO connector 15 mm	trans-parent	white	green	orange	yellow	red	purple
Recommended cuff volume	10 ml	20 ml	35 ml	40-45 ml	50-60 ml	70-80 ml	80-90 ml
Maximum cuff pressure	60 cmH ₂ O						
Drain tube	8 Fr	10 Fr	14 Fr		16 Fr		
Bronchoscopy via ventilation lumen	< 3.0 mm		< 5.0 mm			< 6.0 mm	
External diameter of the tube	9 mm	11.5 mm	14 mm		17.5 mm		

SCOPE OF DELIVERY

Sets	Syringe	Biteblock
Single set, child #0, #1	20 ml	---
Single set, child #2, #2.5	60 ml	#5
Single set, adult #3	60 ml	#6
Single set, adult #4, #5	100 ml	#6

MATERIAL DATA

Tube Cuffs	Silicone
Tube connector	PSU (Polysulfone)
Luer connector	PP (Polypropylene), Silicone, Stainless Steel
Syringe	PP (Polypropylene), IR (Isoprene)
Biteblock	TPU (Thermoplastic Polyurethane)
Fixation Tape	PEF (Polyethylene foam), EL (Elasthane), PA (Polyamide)

SERVICE

Prior to return medical products for complaint/repair they must have run through the complete reprocessing process to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer preserves the right to refuse soiled or contaminated products.

STORAGE AND TRANSPORT CONDITIONS

The maximum shelf life of the product is 5 years.

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

Used or damaged products must be disposed in accordance with the applicable national and international legal regulations.

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**VBM Medizintechnik GmbH**