Cuff Controller

- DE** Gebrauchsanweisung und Serviceanleitung
- EN** Instructions for use and service manual
- FR** Manuel d'utilisation et Manuel de révisions

CE 0123

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Einsteinstrasse 1 | 72172 Sulz a.N. | Germany | Tel.: +49 7454 9596-0
Fax: +49 7454 9596-33 | e-mail: info@vbm-medical.de | www.vbm-medical.de

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TABLE OF CONTENTS

Instructions for use

1.	Intended use	26
2.	Indication	26
3.	Contraindication	26
4.	Requirements for the medical device	26
5.	Safety information	27
6.	Signal words and symbols	28
6.1	Signal words	28
6.2	Symbols on the device	28
6.3	Symbols on the labels	28
7.	Scope of delivery	29
8.	Product description	29
9.	Device overview	30
10.	Start-up	32
11.	Functional check - Use	32
12.	Operation	34
13.	Troubleshooting	35
14.	Cleaning	37
15.	Maintenance	37
16.	Technical data	37
17.	Environmental conditions	38
18.	Shelf life	38
19.	Consumables	38
20.	Disposal	38
21.	EMC table	39

Service manual

22.	Changing units for pressure display	41
23.	Servicing	42
23.1	Maintenance	42
23.1.1	Replacing the battery	42
23.2	Inspection	43
23.2.1	Functional check - Service	43
23.2.2	Leakage test	45
23.2.3	Calibration - Metrological check	46
23.3	Repairs	46

INSTRUCTIONS FOR USE

- ▶ Carefully read through and follow the instructions for use before starting up.

The instructions for use contain important information and directions that must be followed when using the device.

1. INTENDED USE

The Cuff Controller regulates and monitors the cuff pressure of tracheal tubes with high-volume low-pressure balloons. The Cuff Controller reduces the risk of pressure necroses, damage to the tracheal mucosa and aspiration. The Cuff Controller is intended for use in the operating theatre and in the intensive care unit.

2. INDICATION

- Monitoring the cuff pressure of tracheal tubes in intubated patients.
- Other indications are determined by the indication of the tracheal tube used.

3. CONTRAINDICATION

- There are no known contraindications for the Cuff Controller itself.
- Any contraindications result from the contraindications of the tracheal tube used.

4. REQUIREMENTS FOR THE MEDICAL DEVICE

The medical device meets the essential requirements of the Directive of the European Parliament and the Council about medical devices 93/42/EEC and 2007/47/EC.

5. SAFETY INFORMATION

- The product must only be used by trained medical personnel. The cuff pressure must be determined by a doctor. The state of the art in science and technology must be taken into account.
- The maximum length of use is determined by the tracheal tube.
- During defibrillation the connecting tube from the tracheal tube must be removed.
- The device contains ferromagnetic materials and must not be used during magnetic resonance imaging (MRI).
- A functional check must be done every time the device is started up.
- No changes may be made to the device.
- The device must be connected to the tracheal tube via the connecting tube.
- The device must be positioned in such a way that it can be quickly separated from the mains power supply.
- Check the device after a power failure and switch on again if necessary.
- Due to a risk of explosion, the device must not be used in the immediate vicinity (distance < 25 cm (< 9 ¾ inch) of flammable anaesthetic gas mixtures or in atmospheres with an oxygen concentration > 25 %.
- Protect the device including the power supply unit from splashes and moisture. The device must not be operated if liquids have penetrated the device.
- If the device will not be used for a longer period, the battery must be removed from the device (see service manual, section "Replacing the battery").
- In case of interference with other devices, proceed as follows:
 1. Increase the distance between the devices.
 2. Contact the manufacturer of the device.
- If any problems develop or damage occurs, immediately label the device as faulty and send back to the manufacturer for repairs.
- More extensive repairs that are not described in these instructions may only be done by the manufacturer or the authorised service agent of the manufacturer.

CAUTION

The connecting tube is intended for single use and must not be reprocessed. The function of the connection tube will be negatively affected by processing. Re-use involves the potential hazard of an infection.

6. SIGNAL WORDS AND SYMBOLS

6.1 SIGNAL WORDS

Symbol	Description
	CAUTION Indicates a hazardous situation with low risk which, if not avoided, could result in slight or moderate personal injury or property damage.
NOTE	NOTE helps to prevent damage to the device.
	Instruction: Prompt for the user to do something.

6.2 SYMBOLS ON THE DEVICE

Symbol	Description
	"On" / "Off" for a part of the device
	Temporarily interrupt the acoustic alarm
	Plus button: Increase value, confirm menu option in the service menu
	Minus button: Decrease value, select menu option in the service menu
	Charge display for the battery

6.3 SYMBOLS ON THE LABELS

Symbol	Description
	Manufacturer
	Date of Manufacture
	Expiry date
	Batch Code
	Article number
	Serial number
	Device type
	Non-sterile
	Do not use if package is damaged
	Keep away from sunlight and heat
	Keep dry

Symbol	Description
	Do not re-use
	Follow instructions for use
	Important
	Temperature limit
	Air humidity, limit
	Air pressure, limit
	Manufactured without the use of natural latex
	Caution: Federal law restricts this device to sale by or on the order of a physician. For USA and Canada only.
Clean packed	Non-sterile product, packaged in clean-room.
	CE marking
	CE marking with identification number of the notified body.
	Do not dispose of electrical and electronic devices in household waste
	Do not dispose of batteries in household waste
	Protection class II device
	Type BF applied part
	Positive polarity

7. SCOPE OF DELIVERY

- Cuff Controller incl. universal clamp
- External power supply unit Friwo FW8001M/12 A10 (12 V DC / max. 1.5 A)
- Connecting tube with female Luer lock connection and male Luer taper, 200 cm (78 inch) long

8. PRODUCT DESCRIPTION

The cuff pressure of tracheal tubes is monitored using the Cuff Controller and automatically kept constant. By selecting a suitable pressure, the Cuff Controller can prevent pressure necroses on the tracheal mucosa.


Current literature recommends a cuff pressure of 20 to 30 hPa (corresponds to 20 to 30 cmH₂O or 15 to 22 mmHg) for adults. The factory default setting for the Cuff Controller is therefore 25 hPa (corresponds to 25 cmH₂O or 18 mmHg). Depending on the patient, the ventilation method and the tracheal tube, a lower or higher cuff pressure may be necessary. The doctor must check the function of the tracheal tube and determine the ideal cuff pressure for the tracheal tube using conventional methods.

The Cuff Controller offers the following benefits and functions:


- precise pressure regulation, very consistent pressures and automatic pressure compensation of the cuff pressure using microprocessor-controlled pressure regulation
- enables blocking with the lowest pressures and reduces the risk of damage to the tracheal mucosa and pressure necroses
- compensates for slight leaks and reduces the risk of pneumonia as a result of aspiration
- set pressure is automatically regulated
- in case of brief separation from the mains power supply, the value set by the user is saved
- alarm system in case of larger leaks or disconnection
- auto-set function inflates the cuff of the tracheal tube to 25 hPa (corresponds to 25 cmH₂O or 18 mmHg) each time the device is switched on
- adjustable range 0 to 60 hPa (corresponds to 0 to 60 cmH₂O or 0 to 45 mmHg)
- change the pressure using buttons
- easily legible display of the current cuff pressure on the large LCD display
- choose between hPa, cmH₂O and mmHg as the units for the pressure display
- universal clamp to fix to a 25 mm (1 inch) standard rail

Conversion of pressure units: 1 hPa = 1.01973 cmH₂O = 0.75006 mmHg

9. DEVICE OVERVIEW

- 1 - LCD display
- 2 - +/- buttons - menu control and increase or reduce cuff pressure
- 3 - Battery charge display
- 4 - Alarm mute button - suppression of the acoustic alarm
- 5 - On/Off button
- 6 - Male Luer lock connector
- 7 - Connector for the power supply unit
- 8 - Universal clamp for 25 mm (1 inch) standard rail
- 9 - Connecting tube 200 cm (78 inch) for standard tracheal tubes: female Luer lock connector for the connection on the Cuff Controller, male Luer taper attachment for the connection on the tracheal tube
- 10 - External power supply unit with adapter


LCD display


- 1 - Battery and mains symbol
- 2 - Current cuff pressure
- 3 - Pressure units: hPa, cmH₂O or mmHg
- 4 - Target value
- 5 - Target value is currently being adjusted by the user
- 6 - Error message

Display for mains power supply

-  No connection to the mains power supply
-  Connection to the mains power supply


Display for battery status

-  Battery capacity > 75 %
-  Battery capacity 50 - 75 %
-  Battery capacity 25 - 50%
-  Battery capacity < 25 %
-  No battery or faulty battery
-  Battery is charging
-  Battery is not charging


10. START-UP

The device must be operated using the mains power supply.

- ▶ Insert the power supply pack into the plug (1).
- ▶ Insert the power cord into the strain relief (2).

- ▶ Attach the device with the universal clamp (4) to the standard rail (3).
- ▶ Check that the screw (5) is firm.
- ▶ Connect the power supply unit to the mains power supply.

CAUTION
Only use the power supply unit included in the delivery. The mains power supply must conform to the parameters of the power supply unit.

As soon as the device is connected to the mains power supply, the battery is charged. The yellow battery LED turns on.

- ▶ Carry out a functional check each time the device is used (see section "Functional control - Use").

11. FUNCTIONAL CHECK - USE

Each time the device is used, a functional check of the completely closed cuff system must be carried out.

- ▶ Connect the power supply unit to the device.
- ▶ Connect the power supply unit to the mains power supply.
- ▶ Connect the connecting tube and the tracheal tube to be used to the device. If the patient is already intubated with an tracheal tube, any tracheal tube can be used for the functional check.

NOTE

Do not compress, kink or disconnect the connecting tube.

Check automatic inflation

- ▶ Switch the device on using the On/Off button.

The auto-set function inflates the cuff of the tracheal tube to 25 hPa (corresponds to 25 cmH₂O or 18 mmHg) after the device is switched on. The display is illuminated.

NOTE

If the connecting tube or the cuff of the tracheal tube is faulty, the alarm "!!LO" or "!!HF" is triggered.

- ▶ Replace the connecting tube or the tracheal tube.

Check the cuff pressure

- ▶ Increase the cuff pressure with the plus button.

The device must adjust to the new target value.

- ▶ Reduce the cuff pressure with the minus button.

The device must adjust to the new target value.

Check the “!!LO” alarm

- ▶ Separate the connecting tube from the device.

After about 7 seconds the “!!LO” alarm is triggered. The display flashes and an acoustic alarm sounds.

- ▶ Connect the connecting tube to the device.

The alarm is automatically deactivated.

- ▶ Separate the connecting tube from the device.

After about 7 seconds the “!!LO” alarm is triggered. The display flashes and an acoustic alarm sounds.

- ▶ Press the alarm mute button.

The acoustic alarm is silenced for 30 seconds. The acoustic alarm then starts again. The visual alarm remains active. The error message and the background illumination flash on the display.

- ▶ Connect the connecting tube and the tracheal tube to the device.

The device adjusts to the target value and the alarm stops.

Check the connection to the mains power supply

- ▶ Separate the power supply unit from the mains power supply.

The plug symbol on the display disappears.

- ▶ Connect the power supply unit to the mains power supply.

The plug symbol on the display reappears.

NOTE

If the Cuff Controller does not pass the functional check, immediately label the device as faulty and send back to the manufacturer for repairs.

Check the battery capacity

Cause / fault resolution

Empty battery.

- ▶ Connect the device to the mains power supply.
- ▶ Charge the battery for at least 12 hours.

Battery polarity reversed.

- ▶ Insert battery correctly, see service manual section “Maintenance”.

No battery or faulty battery.

- ▶ Insert battery or replace faulty battery, see service manual section “Maintenance”.

- ▶ Switch device off with the On/Off button.

12. OPERATION

After switching on, the device automatically adjusts the cuff pressure to 25 hPa. The current cuff pressure is shown on the display.

If the pressure drops, the device immediately adjusts the pressure to the target value.

If the pressure increases, the device automatically adjusts to the target value with a 5 second delay.

NOTE

Frequent movements of the patient may cause fluctuations in the cuff pressure. The device can compensate these variations in pressure, but this can lead to a shortening of equipment life span.

The device triggers an acoustic and visual alarm if the target value is not reached if, for example, the tracheal tube is disconnected or there is a leak.

- ▶ Connect the power supply unit to the device.
- ▶ Attach the device with the universal clamp to the standard rail.
- ▶ Connect the power supply unit to the mains power supply.
- ▶ Connect the connecting tube to the device.
- ▶ Connect the intubated tracheal tube to the connecting tube.

CAUTION

To prevent loss of pressure, the connecting tube must not be separated from the device and / or the tracheal tube.

NOTE

Do not compress, kink or disconnect the connecting tube.

- ▶ Switch the device on using the On/Off button.

A short BEEP signal sounds.

CAUTION

The doctor must check the function of the tracheal tube and determine the ideal cuff pressure for the tracheal tube using conventional methods. When using higher ventilation frequencies and / or ventilation pressures, it may be that after setting the cuff pressure the seal with the trachea may not be tight. Check the ideal cuff pressure of the tracheal tube again after a few minutes.

The auto-set function inflates the cuff of the tracheal tube to 25 hPa (corresponds to 25 cmH₂O or 18 mmHg) after the device is switched on.

Current literature recommends the following cuff pressures:

- Adults: 20 to 30 hPa (corresponds to 20 to 30 cmH₂O or 15 to 22 mmHg)
- Children: 10 hPa (corresponds to 10 cmH₂O or 7 mmHg)
- ▶ Adjust cuff pressure using the plus or minus button.

Extubation

To prevent injuries to the tracheal mucosa, the following steps must be followed:

- ▶ Switch device off with the On/Off button.

A short BEEP signal sounds.

- ▶ Separate the connecting tube from the tracheal tube.

CAUTION

- ▶ Completely vent the tracheal tube with a syringe.

- ▶ Carry out extubation using conventional procedures.

Alarm system

If an alarm is triggered, a warning sounds and the display flashes. The error is shown on the display (see section "Troubleshooting").

- ▶ The acoustic alarm can be suppressed for 30 seconds using the alarm mute button.

If the fault is not resolved, the acoustic alarm sounds again.

13. TROUBLESHOOTING

Fault / failure	Cause / fault resolution
<p>No information on the display.</p> <p>Background illumination of the display does not light up.</p>	<p>Device is not connected to the mains power supply and / or the power supply unit.</p> <ul style="list-style-type: none"> ▶ Connect the device to the mains power supply and / or the power supply unit. ▶ Switch the device on using the On/Off button.
	<p>There is an interruption to the mains power supply:</p> <ul style="list-style-type: none"> ▶ Switch the device on using the On/Off button.
	<p>In battery mode the background illumination of the display is switched off.</p> <ul style="list-style-type: none"> ▶ Connect the device to the mains power supply.
	<p>Background illumination still does not light up despite connection to the mains power supply. The background illumination of the display is faulty.</p> <ul style="list-style-type: none"> ▶ Immediately label the device as faulty and send back to the manufacturer for repairs.
<p>!!!O</p>	<p>Leak alarm</p> <p>The device cannot establish pressure in the connecting tube and the tracheal tube. For 7 seconds the current cuff pressure is at least 2 hPa lower than the target value.</p> <ul style="list-style-type: none"> ▶ Check the pressure information on the display. ▶ Check the connections and connecting tube for damage and tight fit; reconnect connections if applicable. <div style="border: 1px solid black; padding: 5px; margin: 5px 0;"> <p>NOTE Do not compress, kink or disconnect the connecting tube.</p> </div> <p>Fault is still present. ▶ Replace the connecting tube.</p> <p>Fault is still present. ▶ Replace the tracheal tube.</p> <p>Fault is still present. ▶ Separate the connecting tube from the device. ▶ Manually adjust the cuff pressure of the tracheal tube. ▶ Immediately label the device as faulty and send back to the manufacturer for repairs.</p>

Fault / failure	Cause / fault resolution
<p>!!!HI</p>	<p>Pressure alarm</p> <p>Hose has become compressed during patient repositioning. For 7 seconds the current cuff pressure is at least 2 hPa larger than the target value.</p> <ul style="list-style-type: none"> ▶ If the device does not automatically adjust the pressure, check the connecting tube and tracheal tube for compression and kinks. <div style="border: 1px solid black; padding: 5px; margin: 5px 0;"> <p>NOTE Do not compress, kink or disconnect the connecting tube.</p> </div> <p>Fault is still present.</p> <ul style="list-style-type: none"> ▶ Separate the connecting tube from the device. ▶ Manually adjust the cuff pressure of the tracheal tube. ▶ Immediately label the device as faulty and send back to the manufacturer for repairs.
<p>!!!HF</p>	<p>Hardware error</p> <p>There is a hardware error.</p> <ul style="list-style-type: none"> ▶ A functional check must be done (see section "Functional check - Use"). <p>Fault is still present.</p> <ul style="list-style-type: none"> ▶ Immediately label the device as faulty and send back to the manufacturer for repairs.
<p>BAT</p> 	<p>Battery</p> <p>Empty battery.</p> <p>"BAT" is shown on the display.</p> <p>Device switches itself off after 60 seconds.</p> <ul style="list-style-type: none"> ▶ Connect the device to the mains power supply. ▶ Charge the battery for at least 12 hours.
	<p>Battery polarity reversed.</p> <ul style="list-style-type: none"> ▶ Insert battery correctly, see service manual section "Maintenance".
<p>No battery</p> 	<p>No battery or faulty battery.</p> <p>"No battery" appears on the display after switching on.</p> <p>A continuous tone sounds for 2 seconds and then the "No battery" message disappears.</p> <ul style="list-style-type: none"> ▶ Insert battery or replace faulty battery, see service manual section "Maintenance".

14. CLEANING

- ▶ Switch device off with the On/Off button.
- ▶ Unplug mains plug.
- ▶ Do not immerse the device in liquid or clean with a wet cloth.
- ▶ Clean device; see "Wipe disinfection".

Wipe disinfection

To disinfect the device, disinfect by wiping surfaces with conventional surface disinfectants based on alcohol (e. g., Incides® N or SaniCloth® Active from Ecolab Deutschland GmbH) is suitable. When choosing a product for disinfection, a disinfectant with appropriate ranges of action must be used: bactericidal, fungicidal, tuberculocidal and virucidal. After the wipe disinfection, the devices must be inspected for visible contamination. If necessary, repeat the wipe disinfection. After the wipe disinfection, the function of the device must be checked (see section "Functional check - Use").

CAUTION

The device must not be mechanically processed or sterilised.

15. MAINTENANCE

See service manual section "Maintenance".

16. TECHNICAL DATA

The product is tested for electrical safety in accordance with IEC 60601-1 (DIN EN 60601-1 3rd edition).

Protection class:	II, type BF
Power consumption:	10 VA
Operating pressure:	120 hPa
Adjustable pressure:	max. 60 hPa
Adjustable range:	0 to 60 hPa
Cuff pressure:	max. 60 hPa
Adjustment accuracy:	±1 hPa
Display accuracy:	±1 hPa
Alarm:	Acoustic and visual alarm signal (73 dB)
Connection to the device:	Male Luer lock connector
Connecting tube:	200 cm (78 inch) long with female Luer lock connection and male Luer taper connection
Weight incl. battery and universal clamp:	520 g (18 oz)
Battery type:	Varta 56722 NiMH battery 9 V 200 mAh type E-Block
Power supply unit Friwo FW8001M/12 A10:	Voltage input: 100 - 240 V AC / 50 - 60 Hz Secondary voltage: 12 V DC / max. 1.5 A
Dimensions of the device:	Height: 150 mm (6 inch) Width: 120 mm (4 ¾ inch) Depth: 105 mm (4 inch) Universal clamp: for 25 mm (1 inch) standard rail
Conversion of pressure units:	1 hPa = 1.01973 cmH ₂ O = 0.75006 mmHg

17. ENVIRONMENTAL CONDITIONS

Cuff Controller

Transport:	Temperature	-30 ... +60 °C (-22 ... +140 °F)
	Relative air humidity without condensation	5 ... 95 %
	Atmospheric pressure	70 - 106 kPa
Storage / operation:	Temperature	+10 ... +35 °C (+50 ... +95 °F)
	Relative air humidity without condensation	30 ... 95 %
	Atmospheric pressure	70 - 106 kPa

Connecting Tube

See label of the connecting tube.

18. SHELF LIFE

Cuff Controller


The shelf life of the device is seven years with proper use.

Date of manufacture: see rating plate.

Connecting Tube

Expiration date: see label of the connecting tube.

19. CONSUMABLES

REF	Length	Description
54-05-113	200 cm (78 inch)	Connecting tube with female Luer lock connection and male Luer taper.  male Luer female Luer lock
FTLLP	-	Luer lock plug, female
55-12-500-4	-	Varta 56722 NiMH battery 9 V 200 mAh type E-Block

20. DISPOSAL

The device and the battery must be disposed of separately.

► Remove the battery from the device (see service manual section "Maintenance").

Electrical and electronic devices

 Do not dispose of electrical and electronic devices in household waste. Disposal within the EU must be done in accordance with Directive 2012/19/EU (WEEE Directive). The manufacturer will accept the electrical devices for disposal.

Battery

The device contains a rechargeable battery that is needed for operation and for certain functions.

 Do not dispose of batteries in household waste. The battery must be disposed of in accordance with the relevant national and international legal regulations.

 **CAUTION** Do not open or short circuit the battery, do not immerse the battery in water or throw it in a fire.

Consumables

Used or damaged devices must be disposed of in accordance with the relevant national and international legal regulations.

21. EMC TABLE

Guidance and manufacturer's declaration – electromagnetic emissions		
The Cuff Controller is intended for operation in the electromagnetic environment indicated below. The customer or the user of the Cuff Controller should ensure that it is used in such an environment.		
Emission measurements	Compliance	Electromagnetic environment - Guidance
RF emissions as per CISPR 11	Group 1	The Cuff Controller uses RF energy only for its internal function. Its RF emissions are therefore very low and it is unlikely that this will cause interference to nearby electronic devices.
RF emissions as per CISPR 11	Class B	The Cuff Controller is intended for use in all facilities, including residential areas, and in any facilities connected directly to a public power supply also providing electricity to buildings that are used for residential purposes.
Harmonics as per IEC 61000-3-2	Class A	
Voltage fluctuations / flicker as per IEC 61000-3-3	Met	

Guidance and manufacturer's declaration – electromagnetic immunity			
The Cuff Controller is intended for operation in the electromagnetic environment indicated below. The customer or the user of the Cuff Controller should ensure that it is used in such an environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidance
Electrostatic discharge (ESD) as per IEC 61000-4-2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Floors should be wood, concrete or ceramic tile. If the floor is covered with synthetic materials, the relative humidity should be at least 30 %.
Electrical fast transient / bursts as per IEC 61000-4-4	±2 kV for mains power lines ±1 kV for input / output lines	±2 kV for mains power lines ±1 kV for input / output lines	The quality of the supply voltage should correspond to a typical business or hospital environment.
Surges as per IEC 61000-4-5	±1 kV voltage outer conductor - outer conductor ±2 kV voltage outer conductor - earth	±1 kV voltage outer conductor - outer conductor ±2 kV voltage outer conductor - earth	The quality of the supply voltage should correspond to a typical business or hospital environment.
Voltage dips, short interruptions and voltage variations in the supply voltage as per IEC 61000-4-11	< 5 % U_T for 0.5 cycles (> 95 % drop) 40 % U_T for 5 cycles (60 % drop) 70 % U_T for 25 cycles (30 % drop) < 5 % U_T for 5 sec (> 95 % drop)	< 5 % U_T for 0.5 cycles (>95 % drop) 40 % U_T for 5 cycles (60 % drop) 70 % U_T for 25 cycles (30 % drop) < 5 % U_T for 5 sec (>95 % drop)	The quality of the supply voltage should correspond to a typical business or hospital environment. If the user of the Cuff Controller requires continued function even if there are interruptions to the energy supply, it is recommended to supply the Cuff Controller from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field as per IEC 61000-4-8	3 A/m	3 A/m	If the images are disrupted, it may be necessary to move the Cuff Controller away from sources of mains frequency magnetic fields or to apply magnetic shielding: The mains frequency magnetic field should be measured at the planned location to ensure that it is sufficiently small.
Note: U_T is the mains alternating voltage before application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity

The Cuff Controller is intended for operation in the electromagnetic environment indicated below. The customer or the user of the Cuff Controller should ensure that it is used in such an environment.

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidance
Conducted RF disturbances as per IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	3 V _{eff} 150 kHz to 80 MHz	Portable and mobile radio equipment should not be used any closer to the Cuff Controller including cables than the recommended separation distance calculated from the equation relevant to the transmission frequency. Recommended separation distance: $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz with P as the maximum rated output power of the transmitter in Watt (W) according to the transmitter manufacturer and the recommended separation distance (d) in metres (m). The field strength of stationary radio transmitters, as determined by an electromagnetic site survey ^a , must be lower than the compliance level for all frequencies. ^b Interference is possible near devices that carry the sign "Non-ionising radiation".
Radiated RF disturbances as per IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

Note 1: At 80 MHz and 800 MHz the higher value applies.

Note 2: These guidelines may not apply to all situations. Transmitted electromagnetic waves may be absorbed and reflected by buildings, objects and humans.

^a The field strengths of stationary transmitters such as base stations of radio telephones and mobile terrestrial radio equipment, amateur radio stations, AM and FM radio, and TV broadcasters cannot be precisely predicted in theory. To assess the electromagnetic environment due to stationary RF transmitters, an electromagnetic site survey is recommended. If the measured field strength in the location in which the Cuff Controller is used exceeds the applicable compliance level above, the Cuff Controller must be observed to verify normal operation at each site of use. If abnormal performance is observed, additional measures, such as reorienting or relocating the Cuff Controller, may be necessary.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communication devices and the Cuff Controller

The Cuff Controller is intended for operation in an electromagnetic environment indicated below in which radiated RF disturbances are monitored. The customer or the user of the Cuff Controller can help to prevent electromagnetic interference by maintaining the minimum distances between portable and mobile RF communications equipment (transmitters) and the Cuff Controller as recommended below in accordance with the maximum output power of the communications equipment.

Rated maximum output power of the transmitter in W	Separation distance according to transmitter frequency in m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters with a rated maximum output power not stated in the above table, the recommended distance can be determined using the equation that applies to the relevant column where P is the rated maximum output power of the transmitter in Watt (W) according to the transmitter manufacturer's information.

Comment 1: To calculate the recommended separation distance from transmitters in the frequency range 80 MHz to 2.5 GHz, an additional factor of 10/3 has been used to reduce the probability that mobile or portable communications equipment unintentionally brought close to the patient causes a malfunction.

Note 2: These guidelines may not apply to all situations. Transmitted electromagnetic waves may be absorbed and reflected by buildings, objects and humans.

SERVICE MANUAL

Because the following sections of the service manual are identical to the instructions for use, these sections will not be repeated:

- SAFETY INFORMATION
- SIGNAL WORDS AND SYMBOLS
- DEVICE OVERVIEW
- TROUBLESHOOTING
- CLEANING
- DISPOSAL

22. CHANGING UNITS FOR PRESSURE DISPLAY

SERVICE MENU

```
leakage test
pressure check
X unit
```

UNIT SELECTION

```
unit: hPa
```

- ▶ Press the plus, minus and On/off buttons at the same time for 4 seconds.

“SERVICE MENU” is shown on the display.

- ▶ Select the menu option “Unit” using the minus button.
- ▶ Confirm the selection using the plus button.
- ▶ Select the preferred pressure units hPa, cmH₂O or mmHg using the minus button.

The units for the pressure display are saved when the device is switched off.

- ▶ Switch device off with the On/Off button.

23. SERVICING

The medical device may only be serviced by authorised persons or bodies.

After servicing the medical device, the structural and functional features essential for safety and functionality must be checked.

Only the tasks specified in these instructions for use and service manual may be carried out. If other tasks are carried out on the medical device, all claims under warranty or any other claims are rendered invalid.

23.1 MAINTENANCE

23.1.1 REPLACING THE BATTERY

Replace the battery at least once a year.

CAUTION

Comply with measures to prevent static electricity when changing the battery.

- ▶ Unplug mains plug before opening the housing.
- ▶ Unscrew both screws from the bottom of the device.
- ▶ Slide lower part of the housing from the upper part of the housing.

- ▶ Loosen hook-and-loop tape.
- ▶ Remove battery and replace with a new Varta 56722 NiMH battery 9 V 200 mAh type E-Block:
- ▶ When inserting the new battery, ensure the polarity is correct.
- ▶ Close hook-and-loop tape. The hook-and-loop tape must lie flush with the battery.
- ▶ Slide lower part of the housing onto the upper part of the housing. Do not kink or compress any of the hoses or cables in the interior.
- ▶ Screw both screws back into the bottom of the device. After each battery change, a complete functional check must be carried out (see section "Functional check - Service").
- ▶ Dispose of the old battery (see section "Disposal").
- ▶ Charge the battery for at least 12 hours.

23.2 INSPECTION

- The inspection of the device must be done annually.
- A device inspection involves carrying out the following steps / procedures.

23.2.1 FUNCTIONAL CHECK - SERVICE

- ▶ Connect the power supply unit to the device.
- ▶ Connect the power supply unit to the mains power supply.
- ▶ Connect the connecting tube and any preferred tracheal tube to the device.

NOTE

Do not compress, kink or disconnect the connecting tube.

Check automatic inflation

- ▶ Switch the device on using the On/Off button.

The auto-set function inflates the chosen tracheal tube automatically to 25 hPa after the device is switched on. The display is illuminated.

NOTE

If the connecting tube or the cuff of the tracheal tube is faulty, the alarm “!!LO” or “!!HF” is triggered.

- ▶ Replace the connecting tube or the tracheal tube.

Check the cuff pressure

- ▶ Increase the cuff pressure with the plus button.

The device must adjust to the new target value.

- ▶ Reduce the cuff pressure with the minus button.

The device must adjust to the new target value.

Check the “!!LO” alarm

- ▶ Separate the connecting tube from the device.

After about 7 seconds the “!!LO” alarm is triggered. The display flashes and an acoustic alarm sounds.

- ▶ Connect the connecting tube to the device.

The alarm is automatically deactivated.

- ▶ Separate the connecting tube from the device.

After about 7 seconds the “!!LO” alarm is triggered. The display flashes and an acoustic alarm sounds.

- ▶ Press the alarm mute button.

The acoustic alarm is silenced for 30 seconds. The acoustic alarm then starts again. The visual alarm remains active. The error message and the background illumination flash on the display.

- ▶ Connect the connecting tube and any preferred tracheal tube to the device.

The device adjusts to the target value and the alarm stops.

Check the “!!HI” alarm

- ▶ Gently compress the cuff of the tracheal tube by hand.

The cuff pressure increases and must be 2 hPa higher than the target value. After about 5 seconds the device compensates for the excessive pressure.

- ▶ Continuously increase the cuff pressure by hand.

The device readjusts and tries to compensate for the increasing cuff pressure. After about 2 seconds the “!!HI” alarm is triggered. The display flashes and an acoustic alarm sounds.

- ▶ Release the cuff of the tracheal tube.

The device adjusts to the target value and the alarm stops. The “!!LO” alarm may also be briefly triggered.

Check the connection to the mains power supply

- ▶ Separate the power supply unit from the mains power supply.

The plug symbol on the display disappears.

- ▶ Connect the power supply unit to the mains power supply.

The plug symbol on the display reappears.

NOTE

If the Cuff Controller does not pass the functional check, immediately label the device as faulty and send back to the manufacturer for repairs.

Check the battery capacity

Cause / fault resolution

Empty battery.

- ▶ Connect the device to the mains power supply.
- ▶ Charge the battery for at least 12 hours.

Battery polarity reversed.

- ▶ Insert battery correctly, see service manual section “Maintenance”.

No battery or faulty battery.

- ▶ Insert battery or replace faulty battery, see service manual section “Maintenance”.

- ▶ Switch device off with the On/Off button.

23.2.2 LEAKAGE TEST

The leakage test can be done with the Luer lock plug or any tracheal tube.

Leakage test with the Luer lock plug

- ▶ Connect the connecting tube to the device. Insert the Luer lock plug into the connecting tube.

Leakage test with any tracheal tube

- ▶ Connect the connecting tube and any preferred tracheal tube to the device.
- ▶ Switch the device on using the On/Off button.
- ▶ Increase the cuff pressure to 60 hPa with the plus button.

The cuff is inflated to 60 hPa.

- ▶ Wait about 1 minute until the cuff pressure in the tracheal tube of choice has stabilised.
- ▶ Switch device off with the On/Off button.

The following procedure is identical for the Luer lock plug and the tracheal tube of choice.

The leakage test must be carried out in the service menu.

SERVICE MENU

```
X leakage test
  pressure check
  unit
```

- ▶ Press the plus, minus and On/off buttons at the same time for 4 seconds.
- ▶ Select the menu option "Leakage test" using the minus button.
- ▶ Confirm the selection using the plus button.

LEAKAGE TEST

```
testpressure 60 hPa
diffpressure -1 hPa
.             .
```

The test starts automatically.

The test is complete after about 3 minutes. The display shows "Finished".

- ▶ Document the values displayed.
- ▶ Switch device off with the On/Off button.

The system has no leaks and is OK if the displayed pressure difference (see "diffpressure") is in the range -15 to +15 hPa. If the pressure difference is greater than 15 hPa, the connecting tube and the Luer lock plug (or the tracheal tube of choice) must be checked for damage. The test must be repeated with a new connecting tube and Luer lock plug (or tracheal tube of choice).

NOTE

If the Cuff Controller does not pass the leakage test, immediately label the device as faulty and send back to the manufacturer for repairs.

23.2.3 CALIBRATION - METROLOGICAL CHECK

The calibration is used to check if the measurement accuracy of the device is within the tolerances specified by the manufacturer.

NOTE

The device must not readjusted or changed.

- ▶ Connect the connecting tube to the device.
- ▶ Connect a calibrated reference measuring device (e. g., pressure calibrator Fluke 717 30G) to the connecting tube using the appropriate connectors.


The pressure check must be carried out in the service menu.

SERVICE MENU

```
leakage test
X pressure check
unit
```

- ▶ Press the plus, minus and On/off buttons at the same time for 4 seconds.
- ▶ Select the menu option "Pressure check" using the minus button.
- ▶ Confirm the selection using the plus button.

The current cuff pressure appears on the display, with a greater resolution, below the normal cuff pressure. This cuff pressure must be used when calibrating.

Several cuff pressures must be selected for the calibration. The complete pressure range (0 - 60 hPa) of the device must be covered.

- ▶ Adjust the selected cuff pressure with the plus / minus button.
- ▶ Read the cuff pressure with the higher resolution.
- ▶ Read the cuff pressure on the appropriate reference measuring device.

If the higher resolution pressure reading shown on the device deviates by more than ± 1 hPa, the device must be sent for readjustment to the manufacturer or an authorised body to adjust the display value.

- ▶ Repeat the process until all cuff pressures have been determined using the reference measuring device.
- ▶ Switch device off with the On/Off button.

If the adjustment behaviour of the device is unstable, increase the air volume to be regulated. The reason for the instability is that the adjustment algorithm is optimised for tracheal tubes and tends to become unstable with correspondingly smaller volumes.

23.3 REPAIRS

The medical device may only be repaired by the manufacturer or the authorised service agent of the manufacturer. Rapid turnaround of repairs presumes that the medical device is submitted with as precise a description of the fault as possible.

Prior to returning medical devices, the devices must have been thoroughly cleaned and disinfected (see section "Cleaning") to exclude any risk for the staff of the manufacturer. For safety reasons, the manufacturer preserves the right to refuse soiled or contaminated devices.

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