

POCKET

MONITOR









Cuff Pressure Gauge

CE

Intended to be left blank.

DE	POCKET & MONITOR CUFFDRUCKMESSGERÄT Gebrauchsanweisung 6 - 8	LT	POCKET & MONITOR MANŽETĖS SLĖGIO MATUOKLIS Naudojimo instrukcija..... 34 - 35
EN	POCKET & MONITOR CUFF PRESSURE GAUGE Instructions for use..... 8 - 10	LV	POCKET & MONITOR MANŠĒTES SPIE- DIENA MĒRĪERĪCE Lietošanas instrukcija 36 - 37
BG	POCKET & MONITOR МАНОМЕТЪР ЗА ИЗМЕРВАНЕ НА ВЪТРЕМАНШЕТНО НАЛЯГАНЕ Инструкции за употреба 10 - 12	NL	POCKET & MONITOR CUFFDRUKMETER Gebruiksaanwijzing..... 38 - 40
CS	POCKET & MONITOR MANŽETOVÝ MANOMETR Návod k použití..... 12 - 14	NO	POCKET & MONITOR MANSJETTRYKKMÅLER Bruksanvisning..... 40 - 42
DA	POCKET & MONITOR CUFFTRYKMÅLER Brugsanvisning..... 14 - 16	PL	POCKET & MONITOR MANOMETR DO MANKIETÓW Instrukcja użycia..... 42 - 44
EL	POCKET & MONITOR ΜΕΤΡΗΤΗΣ ΠΙΕΣΗΣ ΑΕΡΟΘΑΛΑΜΟΥ Οδηγίες χρήσης..... 16 - 19	PT	POCKET & MONITOR MANÓMETRO DE PRESSÃO DO BALÃO Instruções de utilização 44 - 46
ES	POCKET & MONITOR MANÓMETRO DE PRESIÓN Instrucciones de uso..... 19 - 21	RO	POCKET & MONITOR DISPOZITIV DE MĂSURARE A TENSIOMETRULUI Instrucţiuni de utilizare 46 - 48
ET	POCKET & MONITOR MANSETIRÕHU MÕÕTESEADE Kasutamisujuhised..... 21 - 23	RU	POCKET & MONITOR МАНОМЕТР ДЛЯ КОНТРОЛЯ ВНУТРИМАНЖЕТНОГО ДАВЛЕНИЯ Инструкция по применению 48 - 50
FI	POCKET & MONITOR KUFFIPAINEMITTARI Käyttöohje 23 - 25	SK	POCKET & MONITOR MANOMETER NA MERANIE TLAKU V MANŽETE CUFF Návod na použitie 51 - 53
FR	POCKET & MONITOR CONTRÔLEURS DE PRESSION DES BALLONNETS Manuel d'utilisation 25 - 27	SL	POCKET & MONITOR MANOMETER ZA MERJENJE TLAKA V TESNILKI Navodila za uporabo 53 - 55
HR	POCKET & MONITOR UREĐAJ ZA MJE- RENJE TLAKA U MANŽETI Upute za uporabu 27 - 29	SV	POCKET & MONITOR KUFFTRYCKSMÅTARE Bruksanvisning..... 55 - 57
HU	POCKET & MONITOR MANDZETTANYOMÁS-MÉRŐ Használati utasítás..... 29 - 31	TR	POCKET & MONITOR KAF BASINÇ ÖLÇER Kullanım kılavuzu 57 - 59
IT	POCKET & MONITOR MANOMETRO Istruzioni per l'uso..... 32 - 33		SYMBOL DESCRIPTION 60 - 63

Deutsch

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English

INTENDED USE

The Cuff Pressure Gauge is used to adjust and control the cuff pressure of high-volume low-pressure cuffs of tracheal tubes and cuffs of supraglottic airways.

Clinical benefit: Cuff Pressure Gauges are used for correct adjustment of the cuff pressure. This reduces the risk of tracheal injury or pharyngolaryngeal complications due to incorrect cuff pressure.

Patient target group: children and adults

Place of use: hospital, pre-hospital and nursing homes or any other places where medically trained personnel can use Cuff Pressure Gauges.

INDICATIONS

- Control and adjustment of the cuff pressure of tracheal tubes and supraglottic airways.
- No other indications are known.

CONTRAINDICATIONS

- None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.

- Before every use, the product must be inspected visually for damages (cracks, breakage, etc.) and a functional check must be performed (see "Functional check"). A faulty product must not be used.
- When using the connecting tubing, the system must always be disconnected at the control valve of the cuff, as otherwise the cuff pressure drops without control (picture 8).
- No changes may be made to the Cuff Pressure Gauge and/or connecting tubing.
- After every use, the Cuff Pressure Gauge must be disinfected (see "Wipe disinfection").
- The Cuff Pressure Gauge and connecting tubing are not sterile.
- The Cuff Pressure Gauge is not suitable for MRI.
- The connecting tubing is intended for single use and must not be reused and/or reprocessed. Reprocessing influences the function of the connecting tubing. Any reuse entails a potential infection hazard.
- The connecting tubing must not be used if the packaging is damaged or the expiration date has expired.

PRODUCT DESCRIPTION (PICTURE 1)

- a - Luer connection for inflation
- b - Hook for fixation
- c - Release valve for pressure adjustment
- d - Inflation bulb for pressure adjustment
- e - Vacuum valve for deflation

FUNCTIONAL CHECK

Leak test Cuff Pressure Gauge

- Seal the luer connection (a) with your finger (picture 2).
- Inflate to 40 cmH₂O using the inflation bulb (d).

The value must remain constant for 2 - 3 seconds.

If the pressure drops, immediately label the Cuff Pressure Gauge as defective and contact the manufacturer.

Leak test connecting tubing

- Connect the connecting tubing to the luer connection (a) of the already tested Cuff Pressure Gauge (picture 3).
- Seal the end of the connecting tubing with your finger (picture 4).
- Inflate to 40 cmH₂O using the inflation bulb (d).

The value must remain constant for 2 - 3 seconds.

If the pressure drops, the connecting tubing is leaking and must be replaced.

Vacuum valve (e)

- ▶ Inflate the Cuff to 40 cmH₂O using the Cuff Pressure Gauge (picture 5).
- ▶ Connect the luer connection of the vacuum valve (e) (picture 6) to the control valve of the cuff (f) (picture 7).
- ▶ Activate the inflation bulb (d) several times until it no longer reverts to its original shape and the cuff is completely empty (picture 7).

USE

- ▶ Connect the Cuff Pressure Gauge to the tracheal tube in place or the supraglottic airway.

For **tracheal tubes** a cuff pressure of 20 to 30 cmH₂O (in children max. 20 cmH₂O) is recommended. This reduces the risk of microaspiration and problems with ventilation caused by too low cuff pressure. Furthermore it reduces the risk of impairment of tracheal mucosal perfusion caused by too high cuff pressure.

For **supraglottic airways** a cuff pressure of 60 cmH₂O maximum is recommended.

- ▶ Check the cuff pressure regularly.
- ▶ In case of pressure drop, increase the cuff pressure using the inflation bulb (d). In case of pressure increase, reduce the cuff pressure using the release valve (c).

NOTE

The pressure drop caused by connection can be compensated using the inflation bulb (d) if necessary.

- ▶ After use, the cuff can be completely emptied using the vacuum valve (e).

WIPE DISINFECTION

CAUTION

- The product must not be mechanically or manually reprocessed or sterilised.
- Do not immerse the product in liquids.
- Do not clean the glass (plastic material) with aggressive cleaning agents.

Wipe disinfection must be performed with commercially available surface disinfectants based on alcohol or QACs (quaternary ammonium compounds). When choosing a product for disinfection, a disinfectant with appropriate ranges of action must be used: bactericidal, yeasticidal and virucidal. After wipe disinfection, inspect the product for visible contamination. If necessary, repeat wipe disinfection. After wipe disinfection, check the function of the product (see "Functional check").

METROLOGICAL CHECK

- ▶ The display accuracy ± 2 cmH₂O of the Cuff Pressure Gauge must be verified every 24 months with calibration.
- ▶ Inflate the Cuff Pressure Gauge to 30, 60 and 90 cmH₂O using the inflation bulb (d).
- ▶ The values must be within the indicated display accuracy. If this is not the case, the manufacturer must be contacted. If the operator does not have the required measuring equip-

ment, the metrological check can be carried out by the manufacturer. In this case, send the Cuff Pressure Gauge to the manufacturer.

If there are any signs before expiry of the 24 months indicating that the Cuff Pressure Gauge fails to adhere to the display accuracy or that the metrological characteristics of the Cuff Pressure Gauge are influenced, a metrological check must be carried out immediately.

Conversion of pressure units:

$$1 \text{ hPa} = 1.02 \text{ cmH}_2\text{O} = 0.75 \text{ mmHg}$$

SERVICE

Prior to the return of medical devices for complaints, the products must be cleaned and disinfected to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

SHELF LIFE

Cuff Pressure Gauge

The shelf life of the Cuff Pressure Gauge is 8 years.

Connecting Tubing

Use-by date: See connecting tubing label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

The Cuff Pressure Gauge and connecting tubing must be disposed of in accordance with the applicable national and international legal regulations.

PRODUCT SPECIFICATIONS

REF	54-04-000	54-05-000
Display accuracy	± 2 cmH ₂ O	
Pressure range	0 - 120 cmH ₂ O	
Weight (incl. packaging)	229 g (8.1 oz)	258 g (9.1 oz)

SCOPE OF DELIVERY

REF	Description	Packaging unit
54-04-000	POCKET Cuff Pressure Gauge with hook	1 pc
	Connecting Tubing for Cuff Pressure Gauge, Length 100 cm (40 inch)	1 pc
54-05-000	MONITOR Cuff Pressure Gauge with hook	1 pc
	Connecting Tubing for Cuff Pressure Gauge, Length 100 cm (40 inch)	1 pc

ACCESSORIES

REF	Description	Packaging unit
54-05-112	Connecting Tubing for Cuff Pressure Gauge, Length 100 cm (40 inch)	10 pc

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