



Tube Holder

Intended to be left blank.

	ASSEMBLY	4
DE	TUBUSHALTEBAND Gebrauchsanweisung	5 - 6
EN	TUBE HOLDER Instructions for use.....	7 - 8
CS	PÁSKA K FIXACI TUBUSU Návod k použití.....	9 - 10
ES	FIJADOR DE TUBO Instrucciones de uso.....	11 - 12
FR	LANIÈRE DE FIXATION Manuel d'utilisation	13 - 14
IT	FISSATUBO Istruzioni per l'uso.....	15 - 16
NL	TUBEVESTIGINGSBAND Gebruiksaanwijzing.....	17 - 18
	SYMBOL DESCRIPTION	19 - 20

ASSEMBLY

Deutsch

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English

INTENDED USE

The Tube Holder enables the tracheal tube to be securely fixed in place. It is suitable for both oral and nasal fixation.

REF	Patient target group
40-8010-033	Child / adult
40-8010-044	Adult

Clinical benefit: The Tube Holder prevents accidental extubation.

Place of use: clinic and preclinic

INDICATIONS

- Fixation of a tracheal tube.
- No other indications are known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU Member State (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Prior to use, a functional check must be performed (see section "Functional check"). A faulty product must not be used.

- No changes may be made to the product.
- The product is not sterile.
- The product is intended for single use and must not be reused and/or reprocessed. Reprocessing influences the function of the product. Any reuse entails a potential risk of cross-contamination.
- The product must not be used if the packaging is damaged or the expiration date has expired.

FUNCTIONAL CHECK

- ▶ Check the product for damages and loose particles.
- ▶ If present, remove loose particles.
- ▶ Test the solidity of the product's two stitchings. Do not use excessive force when doing so.

A faulty product must be disposed of (see "Disposal").

USE

Assembly

- ▶ Loosen the hook and loop fastener from the Tube Holder.
- ▶ Detach the self-adhesive hook tape (a) from the Tube Holder and remove the protective film (Fig. 1).
- ▶ Stick the self-adhesive hook tape around the tracheal tube (Fig. 2).

NOTE

The beige side of the Tube Holder must face outwards and the white side must be on the inside.

- ▶ Wrap the beige narrow side of the Tube Holder with its middle around the hook tape (Fig. 3).

- ▶ Place both ends of the Tube Holder around the patient's neck and fixate with the hook and loop fasteners (Fig. 4).

Monitoring

Patient transfer, patient movement, moisture/humidity and prolonged use may cause the Tube Holder to become loose.

- ▶ Regularly check the position of the Tube Holder and adjust if necessary.

If skin irritations occur in the area of the Tube Holder, the Tube Holder must be replaced with an alternative product.

Removal

The Tube Holder must be removed/replaced after a maximum application period of 30 days or if damaged.

- ▶ Undo both hook and loop fasteners from the Tube Holder. Then loosen the Tube Holder from the hook tape.

SHELF LIFE

Expiration date: see product label

STORAGE AND TRANSPORT CONDITIONS

 	<p>CAUTION</p> <ul style="list-style-type: none"> • Protect from heat and store in a dry place. • Keep away from sunlight and light sources. • Store and transport in the original packaging.
--	---

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

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