
## Tracheostomy Necktape

<b>ASSEMBLY</b>	3
<b>DE TRACHEOSTOMIE-HALS BAND</b>	
Gebrauchsanweisung .....	4 - 5
<b>EN TRACHEOSTOMY NECKTAPE</b>	
Instructions for use.....	6 - 7
<b>BG ПОДДЪРЖАЩА ЛЕНТА ЗА ТРАХЕОСТОМИЯ</b>	
Инструкции за употреба .....	8 - 9
<b>CS TRACHEOSTOMICKÁ KRČNÍ PÁSKA</b>	
Návod k použití.....	10 - 11
<b>DA TRACHEOSTOMIKRAVE</b>	
Bruganvisning .....	12 - 13
<b>EL ΚΟΛΑΡΟ ΤΡΑΧΕΙΟΣΤΟΜΙΑΣ</b>	
Οδηγίες χρήσης.....	14 - 15
<b>ES COLLARÍN PARA TRAQUEOSTOMÍA</b>	
Instrucciones de uso .....	16 - 17
<b>ET TRAHHEOSTOOMIA KAELALINT</b>	
Kasutamisujuhised.....	18 - 19
<b>FI TRAKEOSTOMIAKAULANAuha</b>	
Käyttöohje .....	20 - 21
<b>FR COLLIER DE TRACHÉOTOMIE</b>	
Manuel d'utilisation .....	22 - 23
<b>HR VRATNA TRAKA ZA TRAHEOTOMIJU</b>	
Upute za uporabu .....	24 - 25
<b>HU TRACHEOSTOMIÁS NYAKPÁNT</b>	
Használati utasítás .....	26 - 27
<b>IT FASCIA DI FISSAGGIO PER CANNULA TRACHEOSTOMICA</b>	
Istruzioni per l'uso.....	28 - 29
<b>LT TRACHEOSTOMINĖ APYKAKLĖ</b>	
Naudojimo instrukcija.....	30 - 31
<b>LV TRAHEOSTOMIJAS KAKLA APSĒJS</b>	
Lietošanas instrukcija .....	32 - 33
<b>NL TRACHEOTOMIE-NEKBAND</b>	
Gebruiksaanwijzing.....	34 - 35
<b>NO TRAKEOSTOMI-HALS BÅND</b>	
Bruksanvisning.....	36 - 37
<b>PL OPASKA MOCUJĄCA RURKĘ TRACHEOSTOMIJNĄ</b>	
Instrukcja użycia.....	38 - 39
<b>PT COLAR DE TRAQUEOSTOMIA</b>	
Instruções de utilização .....	40 - 41
<b>RO COLIER DE TRAHEOSTOMIE</b>	
Instrucţiuni de utilizare .....	42 - 43
<b>RU ТРАХЕОСТОМИЧЕСКИЙ ФИКСАТОР</b>	
Инструкция по применению...	44 - 45
<b>SK TRACHEOSTOMICKÁ FIXAČNÁ PÁSKA</b>	
Návod na použitie .....	46 - 47
<b>SL TRAHĚOSTOMSKI OVRATNI TRAK</b>	
Navodila za uporabo .....	48 - 49
<b>SV TRAKEOSTOMIHALS BAND</b>	
Bruksanvisning.....	50 - 51
<b>TR TRAKEOSTOMİ BOYUN BANDI</b>	
Kullanım kılavuzu .....	52 - 53
<b>SYMBOL DESCRIPTION</b>	
.....	54 - 57

## ASSEMBLY

# Deutsch

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REF	Patient target group
40-0008-011	Infants / children
40-0008-033	Children / adults
40-0008-044	Adult, extra long

Clinical benefit: The Tracheostomy Necktape prevents accidental extubation.

Place of use: Hospitals, nursing homes and home care

### INDICATIONS

- Fixation of tracheostomy tubes for tracheostomised patients.

No other indications are known.

### CONTRAINDICATIONS

None known.

### SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- The product may only be used by medically trained personnel who have sufficient knowledge in handling the product.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU Member State (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Prior to use, a functional check must be performed (see section "Functional check"). Damaged products must not be used.
- No changes may be made to the product.
- The product is not sterile.

- The product is intended for single use and must not be reused and/or reprocessed. Reprocessing influences the function of the product. Any reuse entails a potential risk of cross-contamination.

- The product must not be used if the packaging is damaged or the expiration date has expired.

### FUNCTIONAL CHECK

- ▶ Check the product for damages and loose particles.
- ▶ If present, remove loose particles.
- ▶ Test the solidity of the product's two stitchings. Do not use excessive force when doing so.

A faulty product must be disposed of (see "Disposal").

### USE

#### Assembly

- ▶ Undo both hook and loop fasteners from the Tracheostomy Necktape.
- ▶ Place the Tracheostomy Necktape centrally under the patient's neck.
- ▶ Guide both hook and loop fasteners through the slits in the flange of the tracheostomy tube.

#### CAUTION

Do not place the Tracheostomy Necktape too tightly around the patient's neck because of the risk of decubitus formation or vasoconstriction.

- ▶ Fix both hook and loop fasteners to the Tracheostomy Necktape.

### Monitoring

Patient transfer, patient movement, moisture/humidity and prolonged use may cause the Tracheostomy Necktape to become loose.

- ▶ Regularly check that the Tracheostomy Necktape is in position and adjust if necessary.

If skin irritations occur in the area of the Tracheostomy Necktape, the Tracheostomy Necktape must be replaced with an alternative product.

## Removal

The Tracheostomy Necktape must be removed/replaced after a maximum application period of 30 days or if damaged.

- ▶ Undo both hook and loop fasteners from the Tracheostomy Necktape and slide them out of the slits in the flange of the tracheostomy tube.

## SHELF LIFE

Expiration date: see product label

## STORAGE AND TRANSPORT CONDITIONS

 <b>CAUTION</b>	<ul style="list-style-type: none"><li>• Protect from heat and store in a dry place.</li><li>• Keep away from sunlight and light sources.</li><li>• Store and transport in the original packaging.</li></ul>
	

## DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

## PRODUCT SPECIFICATIONS

<b>REF</b>	<b>40-0008-011</b>	<b>40-0008-033</b>	<b>40-0008-044</b>
<b>Length</b>	340 mm (13 inch)	490 mm (20 inch)	590 mm (23 inch)
<b>Width of hook and loop fastener / side parts</b>	6 mm / 10 mm (¼ inch / ½ inch)		
<b>Weight</b> (incl. packaging)	100 g (3.5 oz)	140 g (5 oz)	150 g (5¼ oz)
<b>Packaging unit</b>	10	10	10

## MATERIAL DATA

<b>REF</b>	<b>40-0008-011</b>	<b>40-0008-033</b>	<b>40-0008-044</b>
<b>Foam</b>	Cotton, PU (polyurethane), PA (polyamide)		
<b>Hook and loop fastener</b>	PES (polyester), PP (polypropylene)		

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