



Endotape

ASSEMBLY	3
DE	ENDOTAPE Gebrauchsanweisung	4
EN	ENDOTAPE Instructions for use.....	5
CS	ENDOTAPE Návod k použití.....	6
DA	ENDOTAPE Bruganvising	7
ES	ENDOTAPE Instrucciones de uso	8
FR	ENDOTAPE Manuel d'utilisation	9
IT	ENDOTAPE Istruzioni per l'uso.....	10
RU	ENDOTAPE Инструкция по применению	11
	SYMBOL DESCRIPTION	12

ASSEMBLY

Deutsch

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Clinical benefit: The Endotape prevents accidental extubation.

Patient target group: children and adults

Place of use: clinic and preclinic

INDICATIONS

- Fixation of a tracheal tube, laryngeal tube or laryngeal mask.

No other indications are known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- The product may only be used by medically trained personnel who have sufficient knowledge in handling the product.

- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU Member State (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.

- Prior to use, a functional check must be performed (see section "Functional check"). Damaged products must not be used.

- No changes may be made to the product.

- The product is not sterile.

- The product is intended for single use and must not be reused and/or reprocessed. The function of the product is impaired by reprocessing. Any reuse entails a potential infection hazard.

- The product must not be used if the packaging is damaged or the expiration date has expired.

- The elastic tape contains natural rubber latex. Direct contact may cause allergic reactions.

FUNCTIONAL CHECK

- ▶ Check the product for damage and loose particles.
- ▶ If present, remove loose particles.
- ▶ Check the seams of the product for strength. Do not use excessive force when doing so.

A faulty product must be disposed of (see "Disposal").

USE

Assembly

- ▶ Loosen the wide hook and loop fastener from the Endotape.
- ▶ Apply Endotape with its narrow surface on the chin (oral) or the upper lip (nasal) (Fig. 1).
- ▶ Place the Endotape around the patient's neck and fixate the wide hook and loop fastener to the Endotape.
- ▶ Undo both hook and loop fasteners of the elastic tape.
- ▶ Place elastic tape around the ventilation tube (Figs. 2 and 3).
- ▶ Fix both hook and loop fasteners to the Endotape (Fig. 4).

Monitoring

Patient transfer, patient movement, moisture/humidity and prolonged use may cause the Endotape to become loose.

- ▶ Regularly check the position of the Endotape and adjust if necessary.

If skin irritations occur in the area of the endotape, the Endotape must be replaced with an alternative product.

Removal

The Endotape must be removed/replaced after a maximum application period of 30 days or if damaged.

- ▶ Loosen all hook and loop fasteners from the Endotape and remove the Endotape.

SHELF LIFE

Expiration date: see label on the product.

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international regulations.

PRODUCT SPECIFICATIONS

REF	40-8018-044
Length	540 mm (21 inch)
Width of hook and loop fastener / side parts	15 mm / 30 mm (½ inch / 1¼ inch)
Weight (incl. packaging)	190 g (6.7 oz)
Packaging unit	10

MATERIAL DATA

REF	40-8018-044
Foam	Cotton, PU (polyurethane), PA (polyamide)
Hook and loop fastener	PES (polyester), PP (polypropylene)
Elastic tape	PA (polyamide), rubber (natural latex)
Silicone tube	Silicone
Hook and loop fastener on the elastic tape	PP (polypropylene)

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