




## Water Trap

---

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Einsteinstrasse 1 | 72172 Sulz a. N. | Germany | Tel.: +49 7454 9596-0  
Fax: +49 7454 9596-33 | e-mail: [info@vbm-medical.de](mailto:info@vbm-medical.de) | [www.vbm-medical.de](http://www.vbm-medical.de)

<b>DE</b>	Gebrauchsanweisung .....	4
<b>EN</b>	Instructions for use.....	7
<b>BG</b>	Инструкции за употреба .....	10
<b>CS</b>	Návod k použití.....	14
<b>DA</b>	Brugsanvisning.....	17
<b>EL</b>	Οδηγίες χρήσης .....	21
<b>ES</b>	Instrucciones de uso .....	25
<b>ET</b>	Kasutusjuhend .....	28
<b>FI</b>	Käyttöohje .....	31
<b>FR</b>	Manuel d'utilisation .....	35
<b>HR</b>	Upute za uporabu .....	38
<b>HU</b>	Használati utasítás .....	42
<b>IT</b>	Istruzioni per l'uso.....	45
<b>LT</b>	Naudojimo instrukcija.....	48
<b>LV</b>	Lietošanas instrukcija .....	52
<b>NL</b>	Gebruiksaanwijzing.....	55
<b>NO</b>	Bruksanvisning.....	59
<b>PL</b>	Instrukcja użycia.....	62
<b>PT</b>	Instruções de utilização .....	66
<b>RO</b>	Instrucţiuni de utilizare .....	69
<b>RU</b>	Инструкция по применению.....	73
<b>SK</b>	Návod na použitie .....	77
<b>SL</b>	Navodila za uporabo .....	80
<b>SV</b>	Bruksanvisning.....	83
<b>TR</b>	Kullanım kılavuzu .....	87
	Symbol Description.....	91

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## EN Instructions for use

### INTENDED USE

To collect condensate during ventilation.

Clinical benefit: Supporting the ventilation of a patient.

Patient target group: Infant, child and adult

Place of use: Hospital

### INDICATIONS

- Ventilation of a patient.

No other indications are known.

### CONTRAINDICATIONS

None known.

### SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel with adequate skills in airway management.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.

- Before use, the product must be inspected visually (see section "visual inspection"). A faulty product must be disposed of.
- It is not permitted to make any changes to the product.
- Application time:  $\leq 30$  days
- The product is suitable for conducting oxygen, air, anaesthetic gases and their gas mixtures.
- The reusable product is supplied non-sterile and must be processed completely prior to every use (see chapter "processing").
- The product is MR unsafe.
- The product must not be used if the expiration date has expired.

## PREPARATION PRIOR TO USE

- ▶ Select a suitable water trap according to the procedure planned (see chapter "product specifications").

## VISUAL INSPECTION

- ▶ Check the packaging for damage.
- ▶ Check the product for damage, loose particles and patency.
- ▶ If present, remove loose particles.

A faulty product must be disposed of (see chapter "disposal").

## USE

### WARNING

Constantly check ventilation by standard techniques or observation of chest movements.

- ▶ Assemble the product as shown in the diagram:
  - ▶ Insert the sealing ring (d) into the corresponding groove of the cover (b).
  - ▶ Slide the tappet (c) into the opening of the cover (b) from below as shown in the diagram.
  - ▶ Attach the container (e) of the water trap to the cover (b) by rotation.
  - ▶ Attach the Y-connector (a) to the cover (b) as shown in the diagram.
- ▶ Connect the water trap to the breathing system via the Y-connector (a).  
Ensure correct and firm connection.

## Option

- ▶ Separate the container (e) of the water trap from the cover (b) by rotation, in order to empty collected condensate.
- ▶ Reattach the container (e) of the water trap to the cover (b) by rotation.  
Ensure that the sealing ring (d) is inserted correctly.

## AFTER USE

- ▶ Disconnect the water trap from the breathing system.

## PROCESSING (CLEANING, DISINFECTION, STERILISATION)

### Notes

- The manufacturer is not responsible for damage resulting from improper processing.
- The user is obliged to validate their procedures or the devices and accessories and to comply with the validated parameters during every processing.
- The effectiveness has been verified by an independent test laboratory.
- To achieve effective processing, heavy soiling must not be allowed to dry on the product and must be removed immediately after use.

## CLEANING / DISINFECTION

### Automated cleaning / disinfection

#### Notes

- ▶ Use de-ionised water.
- ▶ Use a washer disinfector which complies with ISO 15883.  
A Miele PG8581 with Cabinet Miele A 105 was used for validation.
- ▶ Do not use drying agents.

#### Procedure validated by the manufacturer

1. Disassemble the product as shown in the diagram.
2. Place the components in the injector carriage in a suitable mesh tray. Secure components with a cover net.

#### CAUTION

Ensure effective rinsing of all hollow spaces.

3. Start program with the following parameters:
  - a. Pre-rinse with water at 20 °C (68 °F) and an exposure time of 1 min.
  - b. Cleaning at 55 °C (131 °F) with an exposure time of 5 min using de-ionised water and a cleaning agent based on alkalinity donors, surfactants and enzymes.  
The cleaning agent Neodisher® MediClean Forte (Dr. Weigert GmbH; dose: 0.5 % (5 ml/l)) was used for validation.
  - c. Rinse with de-ionised water (max. 100 CFU/ml) at 20 °C (68 °F) and an exposure time of 1 min.
  - d. Thermal disinfection at 90 °C (194 °F) - 95 °C (203 °F) for 5 min with de-ionised water.
  - e. Dry at 100 °C (212 °F) for 10 min.
4. Check for visible contamination. If necessary, repeat processing.
5. Check the components according to chapter "checks".
6. Check the components for dryness.
7. Lubricate the sealing ring (d) using e. g. silicone oil.
8. Prepare the product for sterilisation (see chapter "packaging").

#### CHECKS

- ▶ After disinfection, the components must be inspected for damage (cracks, breakage, etc.). A faulty product must be disposed of (see chapter "Disposal").

#### PACKAGING

- ▶ After the cleaning and disinfection procedure, the components must be packed into a sterile barrier system suitable for steam sterilisation. The sterile barrier systems must comply with ISO 11607-1.  
Steriking See-Through Heat Sealable Rolls (K953776) were used for validation.

#### STERILISATION

##### Notes

- ▶ Product and sterilisation packaging must be protected from mechanical damage.
- ▶ Must not be connected to other connectors during sterilisation.
- ▶ Use a steriliser which complies with the requirements of DIN EN 285 or DIN EN 13060

#### Sterilisation validated by the manufacturer

1. Steam sterilise the components using a fractionated vacuum process.
2. Exposure time is 5 min with a sterilisation temperature of 134 °C (273 °F).
3. The drying time is 10 min.

A Tuttnauer ELARA 11-D (K090783) autoclave was used for validation.

## SHELF LIFE

 **WARNING**  
By using the products on patients with suspected prion disease there may be a high risk of transmission. In such case the physician decides if the product is either disposed of (see chapter "Disposal") or processed according to national directives.

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use. The shelf life of the product documented by the manufacturer is 5 years from the date of manufacture. The product can be processed up to 100 times within the shelf life.

Every further use of the product is the responsibility of the user and the product must be inspected beforehand in accordance with chapter "checks".

Use-by date: See product label

## STORAGE AND TRANSPORT CONDITIONS

 **CAUTION**

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.

## SERVICE

Prior to the return of medical devices for complaints/repairs, the products must have run through the complete processing process to exclude any risk for the manufacturer's staff and third parties. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

## DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

 **CAUTION**  
The product may be contaminated with potentially infectious substances of human origin.

## PRODUCT SPECIFICATIONS

<b>REF</b>	66-40-500	66-40-510
Connection to a breathing system via Y-connector (a)	15 mm I.D. / 22 mm O.D.	10 mm O.D.
Volume of container (e)	ca. 55 ml (1.86 fl oz)	ca. 55 ml (1.86 fl oz)
Packaging unit	6	6

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