Stylet

Description	REF	Size	Length	for Tracheal Tube	Packaging unit	Weight (incl. packaging)
Stylet with moveable silicone connector (15 mm I.D.)	90-10-001	small	34 cm	2.5 - 3.5 mm I.D.	10	280 g
	90-10-002	medium	38 cm	4.0 - 4.5 mm I.D.		330 g
	90-10-003	large	45 cm	5.0 - 6.5 mm I.D.		410 g
	90-10-004			7.0 - 11.0 mm I.D.		490 g
Stylet without silicone connector	90-10-111	small	34 cm	2.5 - 3.5 mm I.D.	25	400 g
	90-10-222	medium	38 cm	4.0 - 4.5 mm I.D.		530 g
	90-10-333	large	45 cm	5.0 - 6.5 mm I.D.		750 g
	90-10-444			7.0 - 11.0 mm I.D.		930 g

DE	EINFÜHRUNGSMANDRIN Gebrauchsanweisung3 - 5	IT	MANDRINO Istruzioni per l'uso.....21 - 23
EN	STYLET Instructions for use.....6 - 8	PT	ESTILETE Instruções de utilização24 - 26
CS	INTUBAČNÍ MANDRÉN Návod k použití.....9 - 11	RU	МАНДРЕН Инструкция по применению..27 - 30
DA	INDFØRINGSMANDRINEN Bruganvisning12 - 14	TR	TÜP STILELER Kullanım kılavuzu31 - 33
ES	ESTILETE Instrucciones de uso15 - 17		SYMBOL DESCRIPTION34 - 35
FR	MANDRIN Manuel d'utilisation18 - 20		

Deutsch

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English

INTENDED USE

The tube stylet is a malleable intubation stylet. The shape of the tracheal tube can be adjusted with the aid of the tube stylet.

Clinical benefit: To facilitate orotracheal intubation.

Patient target group: Neonates to adults

Place of use: Hospital and pre-hospital conditions.

INDICATION

The indications for use of the tube stylet arise from the indications that apply to the tracheal tube used.

CONTRAINDICATIONS

The contraindications for use of the tube stylet arise from the contraindications that apply to the tracheal tube used.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- The product must only be used by trained medical personnel with adequate skills in the use of the product.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- No changes may be made to the product.
- The product is not suitable for use when undergoing MRI.

USE

- ▶ Select appropriate size (see cover sheet).
- ▶ Lubricate the tube stylet with a water soluble lubricant.
- ▶ Insert the tube stylet into the tracheal tube and verify easy passage of the stylet in and out of the tracheal tube.

Optional: The version with silicone connector allows to connect the silicone connector with the 15 mm

connector of the tracheal tube to adjust the desired length by moving the tube stylet. Attach the silicone connector to the stylet beforehand (see cover image).

- ▶ Bend the tracheal tube to the desired shape and verify again the easy passage in and out.

CAUTION

The tip of the tube stylet must not protrude from the tracheal tube.

- ▶ While performing intubation, the tube stylet must be fixed in position by hand at the upper end of the tracheal tube.
- ▶ If during withdrawal of the tube stylet unusual resistance is felt, integrity of the tube stylet must be checked.
- ▶ If the tube stylet cannot be withdrawn from the tracheal tube, the tracheal tube with tube stylet as one must be removed.

REPROCESSING

(CLEANING, DISINFECTION, STERILISATION)

General information

- Products are supplied non-sterile and must be unpacked, cleaned, disinfected and sterilised prior to first use.
- Products must be reprocessed prior to each use.
- Manufacturer is not responsible for damages resulting from improper reprocessing.
- It rests with the user to validate their procedures, including devices and accessories, and to comply with the validated parameters during every reprocessing.
- It is recommended to apply an automatic procedure due to a significantly higher efficiency of cleaning and disinfection.
- The effectiveness has been verified by an independent and accredited test laboratory.
- To achieve a more effective reprocessing rough contamination must not get dry on the product and must be removed immediately after use.

Cleaning / Disinfection

Automatic Cleaning / Disinfection

Notes regarding the use of washer disinfectors

- Use de-ionised water.

Notes regarding the use of cleaning and disinfectants

- By the use of alkaline cleaning agents, neutralisation must be performed.
- Don't use drying agents.

Procedure validated by the manufacturer:

1. Disconnect the silicone connector from the tube stylet.
2. Place the products in an injector carriage that is equipped with suitable sieve trays. Secure the products with a cover net.
3. Start programme with the following parameters:
 - a. Pre-rinse with deionised water at 20 °C (68 °F) with an exposure time of 1 min.
 - b. Cleaning at 55 °C (131 °F) with an exposure time of 5 min using de-ionised water and the cleaner "Sekumatic® ProClean" (Dose: 0,5 % (5 ml/l)).
 - c. Neutralisation with "Sekumatic® FNZ" at 20 °C (68 °F) with an exposure time of 2 min (Dose: 0.1 % (1 ml/l)).
 - d. Rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 2 min.
 - e. Thermal disinfection at 93 °C (199 °F) for 5 min with de-ionised water.
 - f. Dry at 100 °C (212 °F) for 10 min.
4. Check for visible contamination. If necessary repeat the reprocessing.
5. Check the devices according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

Manual Cleaning / Disinfection

Notes regarding the use of cleaning and disinfectants

- Prepare the cleaning- and disinfection solution before each manual cycle.
- Use de-ionised water.

Procedure validated by the manufacturer:

1. Disconnect the silicone connector from the tube stylet.
2. Prepare a 2 % (30 ml/l) cleaning- and disinfection solution with "Sekusept® Aktiv" and de-ionised water at 20 °C (68 °F). After 15 min the cleaning and disinfection solution can be used.

3. Clean the products in the cleaning- and disinfection solution by using soft sponges. Areas difficult to access must be cleaned with soft brushes. Hard brushes and other materials which can damage the surface of the products must not be used.
4. Place the devices into the cleaning and disinfection solution and leave for 15 min. Ensure all internal cavities are covered with solution. Ensure effective rinsing of all hollow spaces.
5. Remove the cleaning- and disinfection solution by thoroughly rinsing with de-ionized water. Remaining residues can shorten the products shelf life or can lead to material damage.
6. Dry the products. Accumulation of water must be avoided.
7. Check for visible contamination. If necessary repeat the reprocessing.
8. Check the devices according to the section "Checks" and prepare them for sterilisation (see section "Packaging").

Control

After disinfection and prior to each use, perform the following inspection and function check on the products:

- Check the products for damages (cracks, breakage etc.).
- The metal core must be completely covered by the plastic coating.
- Bring the tube stylet to its original straight shape. To ensure movement inside the tracheal tube, the tube stylet must not be uneven.

Damaged products must not be re-used.

Packaging

- Disconnect the silicone connector from the stylet.

After the cleaning- and disinfection procedure the products must be packed together into a pouch which is suitable for steam sterilisation.

Sterilisation

Notes regarding sterilisation:

- Protect the products and sterilisation pouch against mechanical damage.

Sterilisation validated by the manufacturer

- The products must be steam sterilised by a fractional vacuum process.

- Exposure time is 5 min with a sterilisation temperature of 134 °C (273 °F).

Reusability

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use. Considering the maximum life span of 5 years, reusable devices can be reprocessed up to 10 times if used properly and if the cleaning, disinfection and sterilisation instructions are complied with.

Every further reuse is in the responsibility of the user (see chapter "Control")

WARNING

By using the products on patients with suspected prion disease there might be a high risk of transmission. In such case the physician decides if the product is either disposed (see chapter "Disposal") or reprocessed acc. to the national directives.

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Keep dry and cool.
- Keep away from sunlight and light sources.

SERVICE

Prior to return medical products for complaint/repair they must have run through the complete reprocessing process to exclude any risk for the manufacturer's staff. For safety reasons the manufacturer preserves the right to refuse soiled or contaminated devices.

DISPOSAL

Used or damaged products must be disposed acc. to the applicable national and international legal directives.

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