
Meconium Aspirator

DE EN ES FR IT SYM 8015

VBM Medizintechnik GmbH
Einsteinstrasse 1 | 72172 Sulz a. N. | Germany | Tel: +49 7454 9596-0
Fax: +49 7454 9596-33 | e-mail: info@vbm-medical.de | www.vbm-medical.de

DE MEKONIUM ABSAUGADAPTER

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Instructions for use

INTENDED USE

Adapter between suction line and tracheal tube/suction catheter for tracheal suctioning of aspirated meconium of the newborn.
Clinical benefit: tracheal suctioning of meconium.
Patient target group: newborn.
Place of use: hospital incl. out-of-hospital facilities.

INDICATION

Suspicion of prenatal or perinatal aspiration of meconium-containing amniotic fluid.
Further indications are not known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.
- The product may only be used by a physician or authorised and trained medical personnel who have sufficient knowledge in handling the product.
- Only use the product for the patient target group.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Prior to each use, the products must be inspected visually (see chapter "Visual check"). A faulty product must not be used.
- Application time: < 15 minutes
- The product is suitable for MRI.
- The product is intended for single use and must not be reused and/or processed. The function of the product is impaired by processing. Any reuse entails a potential infection hazard.
- The product is not sterile.
- The product must not be used if the shelf life has elapsed.

VISUAL CHECK

- ▶ Check the packaging for damage.
 - ▶ Check the product for damage (cracks, ruptures) and loose particles.
 - ▶ If present, remove loose particles.
- A faulty product must be disposed of (see "Disposal").

PREPARATION PRIOR TO USE

- ▶ Connect suitable suction line with the barbed end (a) of the Meconium Aspirator.

USE

- CAUTION**
- Execute suctioning with a maximum vacuum of -80 mmHg.
- Do not execute suctioning over a longer period of time.

Application with tracheal tube

- ▶ Connect the Meconium Aspirator to the connector of the already placed tracheal tube via the 15 mm I.D. connector (b).
- ▶ Regulate suction manually by blocking (suction) and unblocking (interruption) the fingertip port (d).

- CAUTION**
- The suctioning must be executed intermittently and with a maximum duration of 2 seconds.

Option: Application with suction catheter

- ▶ Connect the Meconium Aspirator to a suitable suction catheter via the funnel shaped connection (c). Insert suction catheter through the already placed tracheal tube or under visual control.
- ▶ Regulate suction manually by blocking (suction) and unblocking (interruption) the fingertip port (d).

- CAUTION**
- The suctioning must be executed intermittently and with a maximum duration of 2 seconds.

AFTER USE

- ▶ Disconnect Meconium Aspirator from tracheal tube/ suction catheter and proceed with original treatment.

SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture.

Use-by date: See product label

STORAGE AND TRANSPORT CONDITIONS

- CAUTION**
- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

- CAUTION**
- The product may be contaminated with potentially infectious substances of human origin.

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

PRODUCT SPECIFICATIONS

REF	49-30-000
Barbed end (a)	8.0 mm - 10.0 mm O.D. (0.31" - 0.39" O.D.)
Connector (b)	15 mm I.D. (0.59" I.D.)
Funnel shaped connection (c)	≥ 7.5 mm O.D. (≥ 0.3" O.D.)
Fingertip port (d)	8,0 mm I.D. (0.31" I.D.)
Length	57 mm (2.2")
Weight (incl. packaging)	150 g (5.3 oz)
Packaging unit	10

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