




Pocket Introducer

DE	POCKET INTRODUCER Gebrauchsanweisung	IT	POCKET INTRODUCER Istruzioni per l'uso
EN	POCKET INTRODUCER Instructions for use	PT	POCKET INTRODUCER Instruções de utilização
CS	POCKET INTRODUCER Návod k použití	RU	POCKET INTRODUCER Инструкция по применению
DA	POCKET INTRODUCER Brugsanvisning	SV	POCKET INTRODUCER Bruksanvisning
ES	INTRODUCTOR "POCKET" Instrucciones de uso	TR	POCKET INTRODUCER Kullanma kılavuzu
FR	POCKET INTRODUCER Mode d'emploi		SYMBOL DESCRIPTION

VBM Medizintechnik GmbH
Einsteinstrasse 1 | 72172 Sulz a.N. | Germany | Tel.: +49 7454 9596-0
Fax: +49 7454 9596-33 | e-mail: info@vbm-medical.de | www.vbm-medical.de

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MATERIALDATEN
• Polyetherblockamid

EN **POCKET INTRODUCER**
Instructions for Use

INTENDED USE
The Pocket Introducer facilitates orotracheal intubation during laryngoscopy, in particular in difficult airway situations, and is intended for single lumen tracheal tubes with an inner diameter $\geq 6,0$ mm. Clinical benefit: to facilitate orotracheal intubation.
Place of use: hospital and prehospital

INDICATIONS
• To facilitate orotracheal intubation, especially in patients where the visualisation of the glottis is inadequate.
No other indications are known.

CONTRAINDICATIONS
• Complete obstruction of the upper airways.
• Anatomical conditions that make tracheal access impossible (e.g. Cormack & Lehane IV).
No other contraindications are known.

SAFETY INSTRUCTIONS
• Read and follow the instructions for use carefully before using the product and keep them for future reference.
• The product may only be used by a physician or authorised and trained medical personnel who have sufficient knowledge in airway management.
• The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
• Before use, the product must be inspected visually (see "Visual check"). A faulty product must not be used.
• Maximum application time: < 60 minutes
• The product will be irreversibly damaged if the assembled hinges are disassembled. Therefore, disassembly and subsequent reassembly are prohibited.
• If resistance is encountered when advancing the product, the application must be discontinued. In this case, continue the application with a new product.
• The product does not contain any ferromagnetic materials. During MRI treatment, the product can cause artefacts around the area where it is positioned.
• It is not permitted to make any changes to the product.
• The product is intended for single use and must not be reused and/or reprocessed. The function of the product is impaired by reprocessing. Any reuse entails a potential infection hazard.
• The product is sterile (ethylene oxide).

• The product must not be used if the packaging is damaged or the expiration date has expired.

PREPARATION PRIOR TO USE
• Unfold the Pocket Introducer (Fig. 1), making sure not to damage the hinges "a". The plastic connections between the hinges must remain intact.
• Press the hinges together (Figure 2) until the locking mechanism of the hinges engage firmly "Click".
• If the hinges are not properly engaged and/or a hinge "b" is damaged and/or a plastic connection has become detached, the product must be disposed of (see "Disposal").

VISUAL INSPECTION
• Check the locking mechanism of the hinges.
• Check the product for damages (cracks, breaks) and loose particles.
A faulty product must be disposed of (see "Disposal").

APPLICATION
• Lubricate the Pocket Introducer with a water-soluble lubricant.

CAUTION
Do not advance the Pocket Introducer and tracheal tube beyond the Carina tracheae. The graduation along the length of the Pocket Introducer states the distance from the tip and enables the user to estimate the depth of insertion.

• Bend the tip (0 - 10 cm) of the Pocket Introducer to the desired shape if necessary (Fig. 3). Do not bend at the hinges to avoid damaging them.
CAUTION
If resistance is encountered when advancing the product, the application must be discontinued. In this case, continue the application with a new product.

• Using laryngoscopy, insert the Pocket Introducer with the curved tip through the glottis into the trachea (Fig. 4).
• Check the correct positioning of the Pocket Introducer prior to intubation of the tracheal tube.
• Hold the Pocket Introducer in position and advance the tracheal tube over the Pocket Introducer into the trachea (Fig. 5).
• Hold the tracheal tube in position and remove the Pocket Introducer.
• Check the correct positioning of the tracheal tube by common procedure and proceed according to the instructions of the tracheal tube manufacturer.

SHELF LIFE
Use-by date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION
• Protect from heat and store in a dry place.
• Keep away from sunlight and light sources.
• Store and transport in the original packaging.

DISPOSAL

CAUTION
The product may be contaminated when used with potentially infectious substances of human origin.

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

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