Jet-Ventilation Catheter

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English

INTENDED USE/INDICATION

The sterile Jet-Ventilation Catheter allows a safe and atraumatic transtracheal access for oxygenation/ventilation with a manual jet (Manujet III) or an automatic High Frequency Jet Ventilator (HFJV).

Further indications are not known.

Patient	Size
Infant	16 G
Child	14 G
Adult	13 G

CONTRAINDICATIONS

Not known.

SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- The product must only be used by a physician with adequate skills in the use of the product.
- Before use, the products must be inspected visually for damages (cracks, breakage, etc.), according to chapter "Visual check".
- The user must ensure that expiration takes place via the upper airways to reduce the risk of barotrauma.
- The Jet-Ventilation Catheter is made of teflon and not flammable by itself. If used with laser avoid direct contact of the laser beam with the catheter to prevent damage of the catheter.
- Laser impact of the catheter in the presence of elevated O₂ concentration may result in ignition of patient tissue. O₂ concentration in fresh gas should always be below 50 %.

- Ventilation values cannot be compared with the ventilation via a cuffed endotracheal tube due to the loss of air through the upper airway.
- The product is for single use and must not be reprocessed. Reprocessing influences the function of the product. Any reuse involves the potential risk for cross contamination.
- The product is sterile (ethylene oxide). The product must not be used if the packaging is damaged or the shelf life has elapsed.

VISUAL CHECK

- ▶ Check the product for damage and loose particles.

USE

- 1 ▶ Hyperextend the head of the patient. It might be helpful to place a pillow under the patient's neck.
 - ▶ Remove the Jet-Ventilation Catheter from the packaging.

For routine clinical use, we recommend to puncture and insert the catheter only under fiberoptic control.

- ▶ A syringe, attached to the catheter, serves as a handle and for identification for correct positioning of the catheter.
- ▶ Locate the cricothyroid membrane by palpation of the depression between the thyroid and cricoid cartilage. Stabilize this point with index finger and thumb for puncture.

WARNING

As soon as the catheter is located in the trachea, the metal needle must not be introduced further to avoid perforation of the rear tracheal wall.

▶ Puncture the cricothyroid membrane. Because of the sharp tip and conical shape of the needle, an incision is not necessary. After perforation of the cricothyroid membrane, a loss of resistance can be felt.

- 2** ▶ Aspirate air with the syringe to determine the position of the catheter.

If this is possible, the needle is in the trachea. The metal needle must not be introduced further.

- 3** ▶ Push the catheter forward along the needle until the flange rests on the neck.

The metal needle serves as a guide and prevents kinking of the catheter.

- 4** ▶ Now the metal needle can be removed. For safety reasons, check the correct positioning again by aspiration of air via a syringe. If applicable, fill the syringe with 2 - 3 ml sterile saline solution.

CAUTION

- ▶ Secure the catheter safely with the necktape to prevent kinking or dislocation.
- ▶ In case of kinking remove the catheter, reject and replace with a new one.

5 CAUTION

- ▶ Clinically confirm position on the catheter before jetting. If in doubt do not jet!

- ▶ Connect the patient via the Luer Lock connecting tube to the manual jet or the automatic jet ventilator.

The 15 mm connector of the catheter allows connection of a breathing circuit to oxygenate in emergencies (with adjusted flow close the valve according to the volume of the rebreathing bag of the ventilator).

SHELF LIFE

Expiry date: See product label.

STORAGE AND TRANSPORT CONDITIONS

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

Used or damaged products must be disposed in accordance with the applicable national and international regulations.

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