




Endoscopy Mask

CE 0123

Intended to be left blank.

DE	Endoskopiemaske Gebrauchsanweisung5 - 6	LT	Endoskopinė kaukė Naudojimo instrukcija24 - 25
EN	Endoscopy Mask Instructions for use6 - 7	LV	Endoskopijas maska Lietošanas instrukcija26 - 27
BG	Ендоскопска маска Инструкции за употреба8 - 9	NL	Endoscopiemasker Gebruiksaanwijzing27 - 28
CS	Endoskopická maska Návod k použití9 - 10	NO	Endoskopimaske Bruksanvisning29 - 30
DA	Endoskopimaske Bruganvisning11 - 12	PL	Maska endoskopowa Instrukcja użycia30 - 31
EL	Μάσκα ενδοσκόπησης Οδηγίες χρήσης12 - 13	PT	Máscara de endoscopia Instruções de utilização32 - 33
ES	Mascarilla de endoscopia Instrucciones de uso14 - 15	RO	Mască endoscopică Instrucțiuni de utilizare33 - 34
ET	Endoskoopimask Kasutamisujuhised15 - 16	RU	Эндоскопическая маска Инструкция по использованию35 - 36
FI	Endoskoopimaski Käyttöohje17 - 18	SK	Endoskopická maska Návod na použitie36 - 37
FR	Masque de fibroscopie Manuel d'utilisation18 - 19	SL	Maska za endoskopijo Navodila za uporabo38 - 39
HR	Endoskopska maska Upute za uporabu20 - 21	SV	Endoskopimask Bruksanvisning39 - 40
HU	Endoszkópos maszk Használati utasítás21 - 22	TR	Endoskopi maskesi Kullanım kılavuzu41 - 42
IT	Maschera per endoscopia Istruzioni per l'uso23 - 24		Symbol Description43 - 47

Deutsch

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English

INTENDED USE

Endoscopy Mask: To seal mouth and nose during ventilation/oxygenation of a patient and to allow simultaneous insertion of an endoscope.

Bronchoscope Airway: To guide and protect the endoscope/tracheal tube.

Fixation Tape: To hold a face mask in position.

Clinical Benefit: Prevention and control of hypoventilation and desaturation. Reduces interruptions during endoscopic procedures.

Patient target group: Children and adults

Place of use: Hospital

INDICATIONS

- Temporary mechanical ventilation/oxygenation during fibreoptic intubation, bronchoscopy and upper gastrointestinal endoscopies.

No other indications are known.

CONTRAINDICATIONS

- In non-fasted patients with increased risk of aspiration, mask ventilation must not be employed.

No other contraindications are known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- This product must only be used by medically trained personnel.

- Only use the product for the patient target group.

- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.

- Before use, the product must be inspected visually (see section "Visual inspection"). A faulty product must not be used.

- Application time: ≤ 2 hours

- It is not permitted to make any changes to the product.

- The product is suitable for conducting oxygen, air, anaesthetic gases and their gas mixtures.

- The Endoscopy Mask and the Fixation Tape are not sterile.

- The product is intended for single use and must not be reused and/or reprocessed. The function of the product is impaired by processing. Any reuse entails a potential infection hazard.

- The Endoscopy Mask is not suitable for MRI.

- The Bronchoscope Airway and the Fixation Tape are suitable for MRI.

- The Bronchoscope Airway is sterile (ethylene oxide).

- The product must not be used if the packaging is damaged or the expiration date has expired.

PRODUCT DESCRIPTION (FIG. 1)

- a - Membrane
- b - Plug
- c - One-way valve with female Luer connector
- d - Mask cushion
- e - Extension tubing
- f - Bronchoscope Airway (accessory)
- g - Fixation Tape (accessory)

PREPARATION PRIOR TO USE

- ▶ Select the suitable mask size.

VISUAL INSPECTION

- ▶ Check the packaging for damage.
- ▶ Check the product for damage and loose particles.
- ▶ If present, remove loose particles.

A faulty product must be disposed of (see "Disposal").

- ▶ Adjust the volume of the mask cushion (d) with a common syringe via the one-way valve (c).
- ▶ Position the opening of the membrane (a) depending on the oral or nasal endoscopic procedure.
- ▶ Optional: Attach the extension tubing (e) to the Endoscopy Mask.

Optional: For oral use, it is recommended to use a bite protection for the endoscope e. g. VBM Bronchoscope Airway (f) (Fig. 2).

Use of the Bronchoscope Airway is the responsibility of the user.

The Bronchoscope Airway prevents the tongue from blocking the airways in unconscious patients.

- ▶ Select the suitable size.
- ▶ Perform visual inspection (see section "Visual inspection").

Do not insert or remove the Bronchoscope Airway with excessive force.

- ▶ Insert the Bronchoscope Airway.

APPLICATION

CAUTION

The user must ensure efficient ventilation/oxygenation during the entire procedure.

- ▶ Connect the Endoscopy Mask to the breathing system.
- ▶ Apply the Endoscopy Mask firmly over the patient's mouth and nose so that it forms an airtight seal of the upper airways.
- ▶ Hold the Endoscopy Mask in place.
- ▶ Optional: Fix the Endoscopy Mask adequately with a Fixation Tape (g) (Fig. 3).
- ▶ Do not remove the membrane during the procedure.
- ▶ Remove the plug (b) of the membrane to obtain access to the patient's mouth or nose (Fig. 4).
- ▶ Insert the lubricated endoscope through the opening of the membrane (Fig. 5).
- ▶ After the procedure remove the endoscope and close the membrane with the plug.

AFTER USE

- ▶ Remove and disconnect the Endoscopy Mask and the accessories used.

SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture.

Use-by date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

CAUTION

The product may be contaminated with potentially infectious substances of human origin.

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

PRODUCT SPECIFICATIONS

REF	30-40-111	30-40-333	30-40-555	30-40-777
Size of Endoscopy Mask	small	medium	large	large
Connector / Extension tubing	15 mm O.D.	15 mm O.D.	15 mm O.D.	15 mm O.D.
Length of extension tubing	155 mm (6.1 inch)			
Opening of the membrane	Ø 2.0 mm (0.08 inch)	Ø 3.0 mm (0.12 inch)	Ø 5.0 mm (0.20 inch)	Ø 10.0 mm (0.39 inch)
Weight (incl. packaging)	620 g (21.9 oz)	840 g (29.6 oz)	960 g (33.9 oz)	960 g (33.9 oz)
Packaging unit	10	10	10	10

ACCESSORIES

REF	30-40-420-1	30-40-440-1	35-30-255
Description	Bronchoscope Airway		Fixation Tape
#	2	4	-
Length	65 mm (2 1/2 inch)	85 mm (3 3/8 inch)	-
CE Marking			
Weight (incl. packaging)	172 g (6.1 oz)	250 g (8.8 oz)	480 g (16.9 oz)
Packaging unit	10	10	10

MATERIAL DATA

REF	Description	Material
30-40-111, 30-40-333, 30-40-555, 30-40-777	Endoscopy Mask	PVC
	Membrane	Silicone
	Extension tubing	PP, EVA (tubing), PC (connector)
30-40-420-1, 30-40-440-1	Bronchoscope Airway	EVA
35-30-255	Fixation Tape	EPDM

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