








Quicktrach II

CE 0123

DE	Gebrauchsanweisung 4	LT	Naudojimo instrukcija..... 18 - 19
EN	Instructions for use 6	LV	Lietošanas instrukcija 19 - 20
BG	Инструкции за употреба 7 - 8	NL	Gebruiksaanwijzing 20 - 21
CS	Návod k použití 8 - 9	NO	Bruksanvisning..... 21 - 22
DA	Brugsanvisning 9 - 10	PL	Instrukcja użycia..... 22 - 23
EL	Οδηγίες χρήσης..... 10 - 11	PT	Instruções de utilização 23 - 24
ES	Instrucciones de uso 11 - 12	RO	Instrucțiuni de utilizare 25
ET	Kasutamisejuhised 12 - 13	RU	Инструкция по применению. 26 - 27
FI	Käyttöohje 13 - 14	SK	Návod na použitie 27 - 28
FR	Manuel d'utilisation 14 - 15	SL	Navodila za uporabo 28 - 29
HR	Upute za uporabu 15 - 16	SV	Bruksanvisning..... 29 - 30
HU	Használati utasítás 16 - 17	TR	Kullanım kılavuzu 30 - 31
IT	Istruzioni per l'uso 17 - 18		Symbol Description..... 32 - 35

Deutsch

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English

INTENDED USE

Percutaneous cricthyrotomy set for securing the airway in case of upper airway obstruction.

Clinical benefit: ventilation of a patient

Patient target group: adults

Place of use: hospital and prehospital including military use

INDICATION

Upper airway obstruction or ultima ratio if all other attempts to ventilate the patient have failed.

Further indications are not known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them with the product.
- The product must only be used by a physician or by authorized and trained medical personnel with adequate skills in the use of the product.
- Only use the product for the patient target group.
- The needle tip is sharp.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Prior to use, the product must be inspected visually, and a functional check must be performed (see chapter "Visual and functional check"). A faulty product must be disposed of.
- Maximum application time: 4 hours
- Avoid damage to the cuff by sharp-edged objects.
- Constantly check ventilation by standard techniques or chest movements. In case of expiration problems reduce ventilation frequency to avoid the risk of barotrauma.
- It is not permitted to make any changes to the product.
- The product is MR unsafe.

STERILE

- The product is intended for single use and must not be reused and/or processed. The function of the product is impaired by processing. Any reuse entails a potential infection hazard.
- The product is sterile (ethylene oxide).
- Visually inspect sterile barrier system for breaches of integrity immediately prior to application.
- The product must not be used if the packaging is damaged or the expiration date has expired.

VISUAL AND FUNCTIONAL CHECK

- Check the product for damage and loose particles.
 - Remove valve opener.
 - Inflate the cuff with the syringe and check for leaks.
- A faulty product must be disposed of (see "Disposal").
- Evacuate the cuff with the syringe.

USE

- Hyperextend the neck of the patient.
- Evacuate the cuff with the inserted syringe completely. Remove and fill syringe for inflating the cuff with 10 ml air.
- Locate the cricothyroid membrane by palpation of the depression between the thyroid and cricoid cartilage. Stabilise this point with index finger and thumb for puncture.
- Puncture the cricothyroid membrane in a 90° angle.

Due to the sharp and conical needle tip, a prior incision is not necessary. The opening of the trachea is obtained by dilation which also reduces the risk of bleeding.

- Insert the Quicktrach II caudally in a 45° angle further towards the trachea up to the stopper. The stopper helps to prevent the needle from being inserted too deep, helping to avoid a possible perforation of the rear tracheal wall.

- Aspirate air with the syringe to determine the position of the cannula.

If this is possible, the needle tip is in the trachea.

Option: Should aspiration of air be impossible because of an obese neck, then remove the stopper and carefully advance the cannula with the metal needle until it is possible to aspirate air.

- Remove the stopper from the cannula.
- Push the plastic cannula forward with the thumb until the safety clip audibly clicks into position.

This indicates that the tip of the metal needle is covered by the plastic cannula to prevent trauma.

- Further insert the Quicktrach II until the flange rests on the neck.
- Remove metal needle from the cannula.

- Inflate the cuff with the prepared syringe (10 ml air). Check for proper inflation by palpation of the pilot balloon.

- Secure the plastic cannula with the foam necktape.
- Ventilate the patient directly via the 15 mm connector or via the included extension tubing for enhanced flexibility.

EXTUBATION

- Prior to extubation deflate the cuff completely with a syringe.

SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture.

Use-by date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

CAUTION

- The product may be contaminated with potentially infectious substances of human origin.
- The needle tip is sharp.

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

PRODUCT SPECIFICATIONS

REF	30-10-004-1	30-10-904-1
Patient	Adult	
Cannula	4.0 mm I.D.	
Cannula connection	15 mm O.D.	
Connection extension tubing to breathing system	22 mm I.D.	

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