# Pelvic Sling

---

Intended to be left blank.

<b>DE</b>	<b>BECKENSCHLINGE ZUM EINMALGEBRAUCH</b> Gebrauchsanweisung .....4 - 5	<b>LT</b>	<b>VIENREIZLIETOJAMA IEGURŅA JOSTA</b> Lietošanas instrukcija .....23 - 24
<b>EN</b>	<b>PELVIC SLING FOR SINGLE USE</b> Instructions for use .....5 - 6	<b>LV</b>	<b>VIENKARTINIO NAUDOJIMO DUBENS ĮTVARAS</b> Naudojimo instrukcija .....25 - 26
<b>BG</b>	<b>ТАЗОВ ИМОБИЛИЗАТОР ЗА ЕДНОКРАТНА УПОТРЕБА</b> Инструкции за употреба .....7 - 8	<b>NL</b>	<b>BEKKENBAND VOOR EENMALIG GEBRUIK</b> Gebruiksaanwijzing .....26 - 27
<b>CS</b>	<b>PÁNEVNÍ FIXÁTOR NA JEDNO POUŽITÍ</b> Návod k použití .....8 - 9	<b>NO</b>	<b>BEKKENSTABILISATOR TIL ENGANGSBRUK</b> Bruksanvisning .....28 - 29
<b>DA</b>	<b>BÆKKENSLYNGE TIL ENGANGSBRUG</b> Brugsanvisning .....10 - 11	<b>PL</b>	<b>PAS STABILIZUJĄCY MIEDNICĘ DO JEDNORAZOWEGO UŻYCIA</b> Instrukcja użycia .....29 - 30
<b>EL</b>	<b>ΝΑΡΘΗΚΑΣ ΑΚΙΝΗΤΟΠΟΙΗΣΗΣ ΠΥΕΛΟΥ ΜΙΑΣ ΧΡΗΣΗΣ</b> Οδηγίες χρήσης .....11 - 12	<b>PT</b>	<b>CINTA PARA A BACIA PARA UTILIZAÇÃO ÚNICA</b> Instruções de utilização .....31 - 32
<b>ES</b>	<b>CINTURÓN PÉLVICO DE UN SOLO USO</b> Instrucciones de uso .....13 - 14	<b>RO</b>	<b>CENTURĂ PELVIANĂ DE UNICĂ FOLOSINȚĂ</b> Instrucțiuni de utilizare .....32 - 33
<b>ET</b>	<b>VAAGNARIHM ÜHEKORDSEKS KASUTAMISEKS</b> Kasutamisinjuhhis .....14 - 15	<b>SK</b>	<b>PANVOVÁ SLUČKA NA JEDNORAZOVÉ POUŽITIE</b> Návod na použitie .....34 - 35
<b>FI</b>	<b>LANTIOTUKIVYÖ KERTAKÄYTTÖÖN</b> Käyttöohje .....16 - 17	<b>SL</b>	<b>MEDENIČNI PAS ZA ENKRATNO UPORABO</b> Navodila za uporabo .....35 - 36
<b>FR</b>	<b>SANGLE PELVIENNE À USAGE UNIQUE</b> Mode d'emploi .....17 - 18	<b>SV</b>	<b>BÄCKENGÖRDEL FÖR ENGÅNGSBRUK</b> Bruksanvisning .....37 - 38
<b>HR</b>	<b>STABILIZATOR ZDJELICE ZA JEDNOKRATNU UPORABU</b> Upute za uporabu .....19 - 20	<b>TR</b>	<b>TEK KULLANIMLIK KALÇA VE BEL SABİTLEYİCİ</b> Kullanım kılavuzu .....38 - 39
<b>HU</b>	<b>EGYSZER HASZNÁLATOS MEDENCERÖGZÍTŐ ÖV</b> Használati utasítás .....20 - 21		<b>SYMBOL DESCRIPTION</b> .....40 - 41
<b>IT</b>	<b>CINTURA PELVICA MONOUSO</b> Istruzioni per l'uso .....22 - 23		<b>PICTURE</b> .....42

# Deutsch

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## English

### INTENDED USE

The Pelvic Sling is a non-invasive, external pelvis stabiliser for compression and stabilisation of unstable pelvic fractures.

Clinical benefit:

- Reduction of hemorrhage by reducing the intra-pelvic volume.
- External stabilisation of the pelvic ring to restore the anatomic conditions.

Patient target group: patients with a hip circumference between 70 and 140 cm (between 28 and 53 inch)

Place of use: clinic and preclinic

### INDICATION

- Unstable pelvic fracture

No other indications are known.

### CONTRAINDICATION

- Pregnancy
- Open pelvic fractures

No other contraindications are known.

### SAFETY INSTRUCTIONS


- Read and follow the instructions for use carefully before using the product and keep them for future reference.

- The product may only be used by medically trained personnel who have sufficient knowledge in handling the product.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU Member State (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- If a lateral compression fracture is suspected, the Pelvic Sling must be used with caution, as there is a risk that external pressure may cause malformation of the fracture ends.
- Before every use, the products must be inspected visually for damages (cracks, breakage, etc.) (see "Visual inspection").
- The product is intended for single use and must not be reused and/or reprocessed. Reprocessing influences the function of the product. Any reuse entails a potential infection hazard.

### PRODUCT DESCRIPTION

- a - Slide splint
- b - Cuffs
- c - Manometer
- d - Stopcock
- e - Hand inflator

### VISUAL INSPECTION

- ▶ Check for damages (cracks, breakage, etc.).

The faulty product must be disposed of (see "Disposal").

## USE

### Placement of the Pelvic Sling

- ▶ Remove all objects from the patient's pockets in the pelvic region.
- ▶ Choose the correct size (see "Product specifications").
- ▶ Prior to placement of the Pelvic Sling, perform an axial traction and internal rotation on the legs (exception: leg fracture).
- ▶ Use the slide splint (a) to slide the Pelvic Sling underneath the patient (Fig. 1). An assistant can help to lift the pelvis slightly. The printing "Inside" must face towards the patient.
- ▶ Then remove the slide splint (a).

#### WARNING

During placement consider correct positioning of the lateral cuffs (b) as a malposition of the cuffs (b) could lead to a hypoperfusion of the lower extremities.

- ▶ Place the Pelvic Sling so that the lateral cuffs (b) are positioned centrally on the Trochanter Major.
- ▶ The cuff markings (b) should be equally visible on the left and right side (Fig. 2).

The Pelvic Sling must be fitted tightly to the patients' body but without changing the position of the pelvis.

Close the Pelvic Sling similar to a belt.

- ▶ Slide the fixation tape through the buckle and pull it slightly (Fig. 3).
- ▶ Connect the hand inflator (e) to the Pelvic Sling (Fig. 4).
- ▶ Inflate the Pelvic Sling with the hand inflator (e) to the beginning the green zone on the manometer (c) (Fig. 4).

A pressure increase can lead to further significant reduction of the intrapelvic volume. The required pressure increase is at the discretion of the user, but should not exceed the green zone on the manometer (c).

#### WARNING

Ensure that the cuffs (b) are inflated during transport (check the pressure on the manometer (c)).

1. After inflation, turn the stopcock (d) to "ZU - CLOSED" (Fig. 5).
2. Remove the hand inflator (e) (Fig. 5).

- ▶ Fix the knees using the slide splint (a) in order to maintain internal rotation of the lower extremities (Fig. 5).

#### WARNING

The hand inflator (e) must be removed from the Pelvic Sling for X-ray, CT and MRI, as it may affect the image quality in the imaged area.

### Removing the Pelvic Sling

#### WARNING

Before removing the Pelvic Sling, surgical stabilisation of the pelvic ring must be ensured.

## PRODUCT SPECIFICATIONS

REF	34-10-006	34-10-004	34-10-005
Size	Junior S	Standard M	Large L
Hip circumference	70 - 90 cm (28 - 35 inch)	90 - 110 cm (35 - 43 inch)	110 - 140 cm (43 - 53 inch)
Weight (incl. packaging)	530 g (18.7 oz)	640 g (22.6 oz)	680 g (24.0 oz)
Packaging unit	1	1	1

## SHELF LIFE

Expiration date: See product label

## STORAGE AND TRANSPORT CONDITIONS

#### CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

## DISPOSAL

The used or defective products must be disposed of in accordance with the applicable national and international statutory regulations.

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