
SYMBOL DESCRIPTION

|    |                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                              |
|----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| MD | DE - Medizinprodukt<br>EN - Medical Device<br>BG - Медицинско изделие<br>CS - Zdravotnický prostředek<br>DA - Medicinsk udstyr<br>EL - Ιατροτεχνολογικό προϊόν<br>ES - Producto sanitario<br>ET - Meditsiineseade<br>FI - Lääkinnällinen laite                                                                                         | FR - Dispositif médical<br>HR - Medicinski proizvod<br>HU - Orvostechnikai eszköz<br>IT - Dispositivo medico<br>LT - Medicinos priemonė<br>LV - Mediciniska ierice<br>NL - Medisch hulpmiddel<br>NO - Medisinsk utstyr<br>PL - Wyrób medyczny                                                                                                          | PT - Dispositivo médico<br>RO - Dispozitiv medical<br>SK - Zdravotnícka pomôcka<br>SL - Medicinski pripomoček<br>SV - Medicinteknisk produkt<br>TR - Tibbi cihaz                                                                             |
|    | DE - Hersteller<br>EN - Manufacturer<br>BG - Производител<br>CS - Výrobce<br>DA - Producent<br>EL - Κατασκευαστής<br>ES - Fabricante<br>ET - Tootja<br>FI - Valmistaja                                                                                                                                                                 | FR - Fabricant<br>HR - Proizvođač<br>HU - Gyártó<br>IT - Fabbricante<br>LT - Gamintojas<br>LV - Ražotājs<br>NL - Fabrikant<br>NO - Produsent<br>PL - Producent                                                                                                                                                                                         | PT - Fabricante<br>RO - Producător<br>SK - Výrobca<br>SL - Proizvajalec<br>SV - Tillverkare<br>TR - Üretici                                                                                                                                  |
|    | DE - Verwendbar bis<br>EN - Use-by date<br>BG - Да се използва преди<br>CS - Datum expirace<br>DA - Kan anvendes indtil<br>EL - Ημερομηνία λήξης<br>ES - Fecha de caducidad<br>ET - Kõlblik kuni<br>FI - Viimeinen käyttöpäivämäärä                                                                                                    | FR - À utiliser jusqu'au<br>HR - Uporabiti do<br>HU - Lejárati dátum<br>IT - Data di scadenza<br>LT - Naudoti iki<br>LV - Izlietot līdz<br>NL - Te gebruiken tot<br>NO - Utløpsdato<br>PL - Data ważności                                                                                                                                              | PT - Válido até<br>RO - A se utiliza până la<br>SK - Použitelné do<br>SL - Rok uporabnosti<br>SV - Utgångsdatum<br>TR - Son kullanma tarihi                                                                                                  |
|    | DE - Artikelnummer<br>EN - Catalogue number<br>BG - Каталоген номер<br>CS - Číslo výrobku<br>DA - Artikelnummer<br>EL - Αριθμός προϊόντος<br>ES - Número de artículo<br>ET - Artikli number<br>FI - Tuotenumero                                                                                                                        | FR - Numéro d'article<br>HR - Broj artikla<br>HU - Cikkszám<br>IT - Numero articolo<br>LT - Gaminio numeris<br>LV - Artikula numurs<br>NL - Artikelnummer<br>NO - Artikkelnummer<br>PL - Numer artykułu                                                                                                                                                | PT - N.º do artigo<br>RO - Număr articol<br>SK - Číslo výrobku<br>SL - Številka izdelka<br>SV - Artikelnummer<br>TR - Ürün numarası                                                                                                          |
|    | DE - Charge<br>EN - Batch code<br>BG - Партиден код<br>CS - Číslo šarže<br>DA - Batchkode<br>EL - Κωδικός παρτίδας<br>ES - Código de lote<br>ET - Partii kood<br>FI - Erätunnus                                                                                                                                                        | FR - Numéro de lot<br>HR - Kod serije<br>HU - Tételszám<br>IT - Numero di lotto<br>LT - Partijos numeris<br>LV - Partijas kods<br>NL - Batchcode<br>NO - Batchnummer<br>PL - Kod partii                                                                                                                                                                | PT - Código do lote<br>RO - Cod lot<br>SK - Kód šarže<br>SL - Koda serije<br>SV - Batchkod<br>TR - Parti kodu                                                                                                                                |
|    | DE - Gebrauchsanweisung beachten<br>EN - Consult instructions for use<br>BG - Направете справка с инструкциите за употреба<br>CS - Řiďte se návodem k použití<br>DA - Overhold brugsanvisningen<br>EL - Τηρείτε τις οδηγίες χρήσης<br>ES - Véanse las instrucciones de uso<br>ET - Järgige kasutamishuht<br>FI - Noudata käyttöohjeita | FR - Respecter le manuel d'utilisation<br>HR - Slijediti upute za uporabu<br>HU - Kövesse a használati utasítást<br>IT - Rispettare le istruzioni per l'uso<br>LT - Laikyti šis naudojimo instrukcijos<br>LV - Ievērot lietošanas instrukciju<br>NL - Gebruiksaanwijzing opvolgen<br>NO - Følg bruksanvisningen<br>PL - Przestrzegać instrukcji użycia | PT - Cumpra as instruções de utilização<br>RO - Respectați instrucțiunile de utilizare<br>SK - Prečítajte si návod na použitie<br>SL - Upoštevajte navodila za uporabo<br>SV - Läs bruksanvisningen<br>TR - Kullanım kılavuzunu dikkate alın |
| !  | DE - Achtung<br>EN - Caution<br>BG - Внимание<br>CS - Pozor<br>DA - OBS<br>EL - Προσοχή<br>ES - Atención<br>ET - Tähelepanu<br>FI - Huomautus                                                                                                                                                                                          | FR - Attention<br>HR - Oprez<br>HU - Figyelem<br>IT - Attenzione<br>LT - Dėmesio<br>LV - Uzmanību<br>NL - Let op<br>NO - OBS<br>PL - Uwaga                                                                                                                                                                                                             | PT - Atenção<br>RO - Atenție<br>SK - Upozornenie<br>SL - Pozor<br>SV - Observera<br>TR - Dikkat                                                                                                                                              |
|    | DE - Nicht wiederverwenden<br>EN - Do not re-use<br>BG - Да не се използва повторно<br>CS - Nepoužívejte opakovaně<br>DA - Må ikke genbruges<br>EL - Μην επαναχρησιμοποιείτε<br>ES - No reutilizar<br>ET - Ärge taaskasutage<br>FI - Ei saa käyttää uudelleen                                                                          | FR - Ne pas réutiliser<br>HR - Nije za višekratnu uporabu<br>HU - Tilos újra felhasználni<br>IT - Non riutilizzare<br>LT - Nenaudokite pakartotinai<br>LV - Vienreizējai lietošanai<br>NL - Niet opnieuw gebruiken<br>NO - Må ikke gjenbrukes<br>PL - Nie używać ponownie                                                                              | PT - Não reutilizar<br>RO - A nu se reutiliza<br>SK - Nepoužívať opakovane<br>SL - Ni za ponovno uporabo<br>SV - Får ej återanvändas<br>TR - Tekrar kullanmayın                                                                              |

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English

INTENDED USE

The Pelvic Sling is a non-invasive, external pelvic stabiliser for compression and stabilisation of unstable pelvic fractures.  
Clinical benefit:

- Reduction of hemorrhage by reducing the intra-pelvic volume. Minimises internal bleeding by reducing intra-pelvic volume.
  - External stabilisation of the pelvic ring to restore the anatomic conditions.
- Patient target group: patients with a hip circumference between 70 and 140 cm (between 28 and 53 inch)  
Place of use: clinic and preclinic

INDICATION

- Unstable pelvic fracture
- No other indications are known.


CONTRAINDICATION

- Pregnancy
  - Open pelvic fractures
- No other contraindications are known.

SAFETY INSTRUCTIONS

- 
  - Read and follow the instructions for use carefully before using the product and keep them for future reference.
- 
  - The product may only be used by medically trained personnel who have sufficient knowledge in handling the product.
  - The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU Member State (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
  - If a lateral compression fracture is suspected, the Pelvic Sling must be used with caution, as there is a risk that external pressure may cause malformation of the fracture ends.
  - Before every use, the products must be inspected visually for damages (cracks, breakage, etc.) (see "Visual inspection").
  - The product is intended for single use and must not be reused and/or reprocessed. Reprocessing influences the function of the product. Any reuse entails a potential infection hazard.
- 

PRODUCT DESCRIPTION

- a - Slide splint
- b - Cuffs
- c - Manometer
- d - Stopcock
- e - Hand inflator

VISUAL INSPECTION

- ▶ Check for damages (cracks, breakage, etc.).  
The faulty product must be disposed of (see "Disposal").

USE

Placement of the Pelvic Sling

- ▶ Remove all objects from the patient's pockets in the pelvic region.
- ▶ Choose the correct size (see "Product specifications").
- ▶ Prior to placement of the Pelvic Sling, perform an axial traction and internal rotation on the legs (exception: leg fracture).
- ▶ Use the slide splint (a) to slide the Pelvic Sling underneath the patient (Fig. 1). An assistant can help to lift the pelvis slightly. The printing "Inside" must face towards the patient.
- ▶ Then remove the slide splint (a).

**WARNING**

During placement consider correct positioning of the lateral cuffs (b) as a malposition of the cuffs (b) could lead to a hypoperfusion of the lower extremities.

- ▶ Place the Pelvic Sling so that the lateral cuffs (b) are positioned centrally on the Trochanter Major.
- ▶ The cuff markings (b) should be equally visible on the left and right side (Fig. 2).

The Pelvic Sling must be fitted tightly to the patients' body but without changing the position of the pelvis.  
Close the Pelvic Sling similar to a belt.

- ▶ Slide the fixation tape through the buckle and pull it slightly (Fig. 3).
- ▶ Connect the hand inflator (e) to the Pelvic Sling (Fig. 4).
- ▶ Inflate the Pelvic Sling with the hand inflator (e) to the beginning the green zone on the manometer (c) (Fig. 4).

A pressure increase can lead to further significant reduction of the intrapelvic volume. The required pressure increase is at the discretion of the user, but should not exceed the green zone on the manometer (c).

**WARNING**

Ensure that the cuffs (b) are inflated during transport (check the pressure on the manometer (c)).

1. After inflation, turn the stopcock (d) to "ZU - CLOSED" (Fig. 5).
2. Remove the hand inflator (e) (Fig. 5).

- ▶ Fix the knees using the slide splint (a) in order to maintain internal rotation of the lower extremities (Fig. 5).

**WARNING**

The hand inflator (e) must be removed from the Pelvic Sling for X-ray, CT and MRI, as it may affect the image quality in the imaged area.

Removing the Pelvic Sling

**WARNING**

Before removing the Pelvic Sling, surgical stabilisation of the pelvic ring must be ensured.

SHELF LIFE

Expiration date: See product label

STORAGE AND TRANSPORT CONDITIONS

**CAUTION**

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

The used or defective products must be disposed of in accordance with the applicable national and international statutory regulations.

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