
PVC Resuscitator Set

for single use

Intended to be left blank.

DE	Gebrauchsanweisung	4 - 5
EN	Instructions for use	6 - 7
BG	Инструкции за употреба	8 - 9
CS	Návod k použití	10 - 11
DA	Bruganvisning	12 - 13
EL	Οδηγίες χρήσης	14 - 15
ES	Instrucciones de uso	16 - 17
ET	Kasutamisejuhise	18 - 19
FI	Käyttöohje	20 - 21
FR	Manuel d'utilisation	22 - 23
HR	Upute za uporabu	24 - 25
HU	Használati utasítás	26 - 27
IT	Istruzioni per l'uso	28 - 29
LT	Naudojimo instrukcija	30 - 31
LV	Lietošanas instrukcija	32 - 33
NL	Gebruiksaanwijzing	34 - 35
NO	Bruksanvisning	36 - 37
PL	Instrukcja użycia	38 - 39
PT	Instruções de utilização	40 - 41
RO	Instrucțiuni de utilizare	42 - 43
RU	Инструкция по использованию	44 - 45
SK	Návod na použitie	46 - 47
SL	Navodila za uporabo	48 - 49
SV	Bruksanvisning	50 - 51
TR	Kullanım kılavuzu	52 - 53
	Symbol Description	54 - 56

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INTENDED USE

For manual ventilation of patients.
Clinical benefit: Ventilation and oxygen supply of patients.
Patient target group:

Patient	 Body weight
Infant	< 5 kg
Child	5 - 20 kg
Adult	> 20 kg

Place of use: Hospital and prehospital

INDICATIONS

Any emergency with limited respiratory function or signs of hypoxia.
No other indications are known.

CONTRAINDICATIONS

None known.

SAFETY INSTRUCTIONS

• Read and follow the instructions for use carefully before using and keep them for future reference.

- This product must only be used by medically trained personnel.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before use, visually inspect the product for damage (cracks, breakage, etc.) and carry out a functional check (see chapter "Functional check"). A faulty product must be disposed of.
- Do not disassemble the pressure relief valve.
- Use original VBM individual parts only. Using parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this..

• The product is intended for single use and must not be reused and/or processed. The function of the product is impaired by processing. Any reuse entails a potential infection hazard.

• The product must not be used if the shelf life has elapsed.

• The product is not sterile.

• The product is MR unsafe.

FUNCTIONAL CHECK

The valve functions must be checked as follows to ensure proper functionality of the resuscitator:
A reservoir bag is needed to complete the test procedures described below:

1.0 Intake valve

a) ▶ Firstly compress the resuscitator with one hand and then close the Non-Rebreathing Valve of the resuscitator with the other hand. Release your grip on the resuscitator.
Rapid re-expansion of the resuscitator confirms efficient air intake via the intake valve.

b) ▶ Close the Non-Rebreathing Valve of the resuscitator and try to compress the resuscitator.
If the resuscitator cannot be compressed with reasonable force, or if resuscitator compression forces the air out between your hand and the neck of the resuscitator, the intake valve efficiently prevents backward leakage of air.

2.1 Non-Rebreathing Valve (pressure relief valve in closed position)
▶ Hold a reservoir bag over the Non-Rebreathing Valve and use your thumb to press on the reservoir bag connector.
▶ Ensure a tight seal between the Non-Rebreathing Valve and reservoir bag.
▶ Compress the resuscitator with your other hand several times. Ensure that the lip valve opens during compression.

Filling of the reservoir bag confirms that the Non-Rebreathing Valve efficiently directs air to the patient.

2.2 Non-Rebreathing Valve (pressure relief valve in opened position)
▶ Close the Non-Rebreathing Valve with your thumb and compress the resuscitator several times. The resuscitator can be compressed and the visible and audible opening of the pressure relief valve confirms its proper operation.

3.0 Reservoir valve (in the intake valve)
a) ▶ Proceed as described above in 2.1 in order to fill the reservoir bag with ambient air.
▶ Attach the filled reservoir bag to the intake valve and press on the reservoir bag.

The visible rise of the diaphragm when compressing the reservoir bag confirms that the reservoir bag efficiently vents excessive gas to the atmosphere.

b) ▶ Proceed as described above in 2.1 in order to fill a reservoir bag with ambient air.
▶ Attach the filled reservoir bag to the intake valve.
▶ Perform several compression-release cycles on the resuscitator until the reservoir bag is flat and empty. Rapid re-expansion of the resuscitator confirms efficient air intake.

USE

- ▶ Optional: Connect reservoir bag and oxygen tubing to the resuscitator.
- ▶ Connect the mask to the Non-Rebreathing Valve.
- ▶ Connect the oxygen tubing to a regulated oxygen source.
- ▶ Adjust the oxygen flow rate. The reservoir bag expands completely during inflow and nearly collapses as the compressed resuscitator refills during exhalation.
- ▶ Prior to connecting to a patient, check the function of the resuscitator and make sure that all connections are in proper position. Observe the intake valve, reservoir bag and Non-Rebreathing Valve in all phases of ventilation; leakages must not occur.

- Apply the Mask firmly over the patient's mouth and nose so that it forms an airtight seal of the upper airway. If the mask does not fit properly, the ventilation may be affected.
- Hold the Mask in position during ventilation.
- Compress the resuscitator bag to deliver a breath. Observe the patient's chest wall rise to confirm inspiration.
- Release the resuscitator bag to allow patient exhalation. Observe the chest wall fall to confirm exhalation.

SHELF LIFE

The shelf life of the product is 5 years from the time of manufacture.

Use-by date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

CAUTION

- The product may be contaminated with potentially infectious substances of human origin.

PRODUCT SPECIFICATIONS

REF	Patient	Resuscitator bag	Mask	Reservoir bag	Oxygen tubing	Packaging unit
84-10-195	Adult	1500 ml	#5	2500 ml	200 cm	6
84-10-295	Child	450 ml	#3	2500 ml	200 cm	6
84-10-395	Infant	280 ml	#1	600 ml	200 cm	6

TECHNICAL DATA

Patient	Infant	Child	Adult
Inspiratory resistance	≤ 5 cmH ₂ O @ 60 l/min		
Expiratory resistance	≤ 5 cmH ₂ O @ 60 l/min		
Pressure relief valve	40 cmH ₂ O		
Dead space	≤ 5 ml + 10 % of the supplied volume	≤ 5 ml + 10 % of the supplied volume	≤ 5 ml + 10 % of the supplied volume
Delivered volume	≥ 20 ml	≥ 150 ml	≥ 600 ml

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