
Intended to be left blank.

Intended to be left blank.

Surgicric II

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English

INTENDED USE

Surgical cricothyrotomy set for securing the airway in case of upper airway obstruction or as ultima ratio if all other attempts to ventilate the patient have failed. Surgicric II can be used with resuscitation bags or other standard ventilation systems.

Clinical benefit: Ventilation of a patient

Patient target group: Adults

Location of use: Hospital and prehospital including military use

INDICATION

- Emergency airway access via the cricothyroid membrane. Life-threatening dyspnoea that cannot be controlled in any other way

No other indications are known.

CONTRAINDICATION

- Not known.

SAFETY INFORMATION

 Read and follow the instructions for use carefully before using the product and keep them with the product.

 The product must only be used by a physician or by authorized and trained medical personnel with adequate skills in the use of the product.

- Only use the product for the patient target group.
- Only use the products for patients for whom a 6 mm tube is suitable for cricothyrotomy.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before use, the product must be inspected acc. to a functional check (see "Functional check"). A faulty product must not be used.
- The product serves to temporarily secure the airway. If possible it must be replaced by a permanent airway. Maximum recommended period of application: 4 hours
- Tube and extension tubing are suitable for conducting oxygen and air.
- The enclosed extension tubing is not suitable for use with mechanical ventilation systems.
- Avoid damage to the cuff of the tube by sharp-edged objects.
- Do not kink the inflation line of the tube.
- Use the scalpel, tracheal hook, speculum and tube with appropriate force
- Lock the safety mechanism of the scalpel for disposal only after use
- It is not permitted to make any changes to the product.
- The product is not suitable for MRI.

STERILE

- The product is intended for single use and must not be reused and/or processed. The function of the product is impaired by processing. Any reuse entails a potential infection hazard.
- The product is sterile (ethylene oxide).
- Visually inspect sterile barrier system for breaches of integrity immediately prior to application.
- The product must not be used if the packaging is damaged or the expiration date has expired.

FUNCTIONAL CHECK

- ▶ Check the products for damage and loose particles.
- ▶ Inflate the cuff with the syringe and check for leakage.

A faulty product must be disposed of (see "Disposal").

- ▶ Remove all air from the cuff with the syringe.

USE

- ▶ Hyperextend the patient's head.
- ▶ Locate the cricothyroid membrane by the palpable depression between the thyroid cartilage and the cricoid cartilage and stabilise this site with index finger and thumb.

NOTE

- For clear identification of the cricothyroid membrane, an additional incision through skin and soft tissue may be necessary in obese patients or patients with neck trauma.

- ▶ Make a midline vertical skin incision using the scalpel.
- ▶ If the cricothyroid membrane is not identified after skin incision, scissors can be used to bluntly dissect the soft tissue down to the cricothyroid membrane.
 - ▶ To keep the incision open, the tracheal hook can be used as an additional aid if necessary.
- ▶ Make a horizontal incision through the cricothyroid membrane.
 - ▶ Remove the scalpel.
- ▶ Insert the closed speculum horizontally into the incision and advance it into the trachea.
 - ▶ Rotate the speculum 90° in the caudal direction.
 - ▶ Lift the closed speculum slightly so that the tip of the speculum points caudally.
- ▶ Open the speculum.
 - ▶ Insert the tube into the trachea via the speculum until the flange abuts the speculum.
 - ▶ Remove the speculum.
 - ▶ Insert the tube deeper. The flange must be flush with the patient's neck.
 - ▶ For a possible position check, connect the syringe to the dilator and aspirate air.
 - ▶ If necessary, the tube can be repositioned after opening the flange. The centimetre markings on the tube indicate the insertion depth.
- ▶ Then tighten the flange to prevent axial dislocation of the tube.
- ▶ Inflate the cuff to the minimum pressure required to seal the trachea using the syringe and check the cuff at regular intervals.
- ▶ Unlock and remove the dilator. Do not change the position of the tube when doing so.

8

CAUTION

Refrain from excessive rotating movements and traction. Avoid accidental disconnection or kinking of the inserted tube.

- ▶ Connect the tube to a ventilation system.
 - If necessary, use the enclosed extension tubing when using manual ventilation systems.
 - During mechanical ventilation ensure for proper strain-relief.
- ▶ Ventilate the patient. Verify correct placement and successful ventilation at regular intervals by a suitable method.
- ▶ Place the necktape around the patient's neck and connect it firmly to the flange.
- ▶ Check the tube for secure hold.
- ▶ Tracheal suctioning is possible through the connector of the extension tubing with a suction catheter \leq 16 Fr.
- ▶ Prior to extubation deflate the cuff completely with a syringe. The tube is radiopaque.

SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture.

Expiration date: See product label

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

CAUTION

- The product may be contaminated with potentially infectious substances of human origin.
- Injury hazard, the scalpel is sharp.

The used or faulty product must be disposed of in accordance with the applicable national and international regulations.

WEIGHT

REF	30-08-117-1
Weight (incl. packaging)	380 g (13.4 oz)
Packaging unit	1

MATERIAL DATA

Scalpel	Blade: Carbon steel; handle: PA, SAN	
Tracheal Hook, Scissors, Speculum	Stainless steel	
Tube	Tube: PVC; connector: PSU; flange: PP	
Syringe	PP, IR	
Extension tubing	PP	
Swivel connector	Connector: PC; suction cap: SJ; o-ring: EPDM	
Necktape	Foam	Inside: PU; outside: PA
	Hook and loop fastener	PES, PP

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