Surgicric III

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English

INTENDED USE

Percutaneous cricothyrotomy set for securing the airway in case of upper airway obstruction or as last resort if all other attempts to ventilate the patient have failed. Surgicric III can be used with resuscitation bags or other standard systems.

INDICATION

- Emergency airway access via the cricothyroid membrane. Life-threatening dyspnoea that cannot be controlled in any other way.

CONTRAINDICATION

- Patients who are not suitable for use of a 6 mm tracheal tube during cricothyrotomy.
- Further contraindications are not known.

SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- The product must only be used by a physician with adequate skills in the use of the product.
- Before use, the products must be inspected for damages (cracks, breakage, etc.), according to chapter "Functional check".
- The product serves to temporarily secure the airway. If possible it must be replaced by a permanent airway. Maximum recommended period of application: 4 hours.
- Avoid damage to the cuff of the tracheal tube by sharp-edged objects.
- Introducer needle 13 G, tracheal tube and extension tubing are suitable for conducting oxygen and air.
- Only introducer needle 13 G, tracheal tube, syringe, extension tubing and necktape do not contain any ferromagnetic materials.
- The product is for single use and must not be reprocessed. Reprocessing influences the function of the product. Any reuse involves the potential risk for cross contamination.
- The product is sterile (ethylene oxide). The product must not be used if the packaging is damaged or the shelf life has elapsed.

FUNCTIONAL CHECK

- ▶ Check the products for damage and loose particles.
 - ▶ Inflate the cuff with the syringe and check for leakage.
- Damaged or leaking products must not be used.
- ▶ Remove all air from the cuff with the syringe.

USE

- ▶ Hyperextend the patient's head.

NOTICE

For obese patients or patients with neck trauma, an incision through the skin and soft tissue may be necessary in order to clearly identify the cricothyroid membrane.

- ▶ A syringe, attached to the catheter, serves as a handle and for identification for correct positioning of the catheter.
- 1 ▶ Locate the cricothyroid membrane by palpation of the depression between the thyroid and cricoid cartilage. Stabilize this point with index finger and thumb for puncture.

2 WARNING

As soon as the catheter is located in the trachea, the metal needle must not be introduced further to avoid perforation of the rear tracheal wall.

- ▶ Puncture the cricothyroid membrane.
- 3 ▶ Check for correct position of the introducer needle by aspirating air with the syringe. If this is possible, the needle tip is in the trachea.
 - ▶ Do not further advance the metal needle into the trachea.
 - 4 ▶ Insert the catheter along the metal needle in caudal direction until the flange rests on the neck. The metal needle serves as a guide and prevents kinking of the catheter.
 - ▶ Remove the metal needle. It is possible to check for correct position by aspirating air with the syringe through the catheter.

- ▶ For temporary oxygenation of the patient the catheter can be connected to a breathing system. Connector of the catheter: 15 mm O.D. The catheter must first be fixed with the necktape.
- 5** ▶ Attach the guidewire to the catheter. For depth orientation there are graduation marks every 10 cm on the guidewire.
- 6** ▶ Advance the guidewire to the 20 cm graduation mark.
- 7** ▶ Hold the guidewire in place and remove the sleeve of the guidewire.
 - ▶ Open the necktape at one side of the catheter.
- 8** ▶ Hold the guidewire in place and remove the catheter with the necktape.
- 9** ▶ Use the safety scalpel and make a 5 mm horizontal incision on each side of the guidewire.
- 10** ▶ Railroad the pre-assembled dilator with tracheal tube over the guidewire into the trachea.
- 11** ▶ The pre-adjusted flange must be flush with the patient's neck. During this procedure make sure to hold the guidewire at the proximal end.
 - ▶ If required by the patient's anatomy the tracheal tube must be inserted deeper into the trachea after loosening the flange. Then tighten the flange to prevent axial dislocation of the tracheal tube.
 - ▶ Remove the guidewire from the tracheal tube.
 - ▶ Inflate the cuff to the minimum pressure required to seal the trachea using the syringe and check the cuff at regular intervals.
 - ▶ It is possible to check for correct position by aspirating air with the syringe through the dilator.
- 12** ▶ Unlock and remove the dilator. Make sure that the position of the tracheal tube is not changed.

CAUTION

Refrain from excessive rotating movements and traction. Avoid accidental disconnection or kinking of the inserted tracheal tube.

- ▶ Attach the extension tubing to the tracheal tube and connect to a manual breathing system. During mechanical ventilation ensure for proper strain-relief.

- ▶ Ventilate the patient. Verify correct placement and successful ventilation at regular intervals by a suitable method.
- ▶ Place the necktape around the patient's neck and connect it firmly to the flange.
- ▶ Check the tracheal tube for secure hold.
- ▶ Tracheal suctioning is possible through the connector of the extension tubing with a suction catheter ≤ 16 Fr.
- ▶ Prior to extubation deflate the cuff completely with a syringe.
- ▶ The tracheal tube is radiopaque.

SHELF LIFE

Expiry date: see product label.

STORAGE AND TRANSPORT CONDITIONS

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

Used or damaged products must be disposed in accordance with the applicable national and international regulations.

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