
Intended to be left blank.

Intended to be left blank.

## Surgicric III

---

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# English

## INTENDED USE

Percutaneous cricothyrotomy set for securing the airway in case of upper airway obstruction or as ultima ratio if all other attempts to ventilate the patient have failed. Surgicirc III can be used with resuscitation bags or other standard ventilation systems.

Clinical benefit: Ventilation of a patient.

Patient target group: Adult

Location of use: Hospital and prehospital including military use.

## INDICATIONS

- Emergency airway access via the cricothyroid membrane. Life-threatening dyspnoea that cannot be controlled in any other way.
- No other indications are known.

## CONTRAINDICATIONS

None known.

## SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them with the product.
- The product must only be used by a physician or by authorized and trained medical personnel with adequate skills in the use of the product.
- Only use the product for the patient target group.
- Only use the product on patients for whom a 6 mm I.D. tube is suitable for cricothyrotomy.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before use, the product must be inspected acc. to a functional check (see chapter "functional check"). A faulty product must not be used.
- The product serves to temporarily secure the airway. As soon as possible it must be replaced by a permanent airway. Maximum application time: 4 hours.
- Avoid damage to the cuff of the tube by sharp objects.
- Constantly check ventilation by standard techniques or chest movements.
- It is not permitted to make any changes to the product.
- The tube and the extension tubing are suitable for conducting oxygen and air.
- The extension tubing enclosed is not suitable for use with mechanical ventilation systems.
- Do not kink the inflation line of the tube.
- Use the scalpel, introducer catheter and tube with appropriate force.
- Lock the safety mechanism of the scalpel for disposal only after use.
- The product is not suitable for MRI.

STERILE/EO

- The product is intended for single use and must not be reused and/or processed. The function of the product is impaired by processing. Any reuse entails a potential infection hazard.
- The product is sterile (ethylene oxide).
- Visually inspect sterile barrier system for breaches of integrity immediately prior to application.
- The product must not be used if the packaging is damaged or the expiration date has expired.

## FUNCTIONAL CHECK

- ▶ Check the product for damage and loose particles.
  - ▶ Inflate the cuff with the syringe (b) and check for leakage.
- A faulty product must be disposed of (see chapter "disposal").
- ▶ Deflate the cuff completely using the syringe.

## USE

- ▶ Hyperextend the patient's head.
  - ▶ Attach the syringe to the introducer catheter (f), which serves as a handle and for position control.
- ▶ Locate the cricothyroid membrane by palpation of the depression between the thyroid and cricoid cartilage. Stabilise this point with index finger and thumb for puncture.

### NOTICE

For obese patients or patients with neck trauma, an incision through the skin and soft tissue may be necessary in order to clearly identify the cricothyroid membrane.

2

### WARNING

As soon as the introducer catheter is located in the trachea, the needle must not be inserted further to avoid perforation of the rear tracheal wall.

3

- ▶ Puncture the cricothyroid membrane with the introducer catheter.
- ▶ Check for correct position of the introducer catheter by aspirating air with the syringe. If this is possible, the needle tip is in the trachea.

4

- ▶ Do not further advance the needle into the trachea.
- ▶ Insert the cannula along the needle in caudal direction until the flange rests on the neck. The needle serves as a guide and prevents kinking of the catheter.

5

- ▶ Remove the needle. To check the position again if necessary, attach the syringe to the cannula and aspirate air.

6

- ▶ Attach the guidewire (g) to the cannula. For depth orientation there are graduation marks every 10 cm on the guidewire.

7

- ▶ Advance the guidewire to the 20 cm graduation mark.

8

- ▶ Remove the sleeve of the guidewire by holding the guidewire without removing it from the cannula.

9

- ▶ Hold the guidewire in place and remove the cannula.

10

- ▶ Use the scalpel (a) to make a horizontal incision of approx. 5 mm on each side of the guidewire.

11

- ▶ Railroad the tube (c) with the pre-assembled dilator over the guidewire into the trachea.

12

- ▶ The flange must be flush with the patient's neck. During this procedure make sure to hold the guidewire at the proximal end.

- ▶ If required by the patient's anatomy the tube may need to be inserted deeper into the trachea after loosening the flange. Then tighten the flange to prevent axial dislocation of the tube.

- ▶ Remove the guidewire from the tube.
- ▶ Inflate the cuff to the minimum pressure required to seal the trachea using the syringe and check the cuff at regular intervals.

- ▶ It is possible to check for correct position by aspirating air with the syringe through the dilator.

- ▶ Unlock and remove the dilator. Make sure that the position of the tube is not changed.

### CAUTION

Refrain from excessive rotating movements and traction. Avoid accidental disconnection or kinking of the inserted tube.

- ▶ Connect the tube to a ventilation system.
- ▶ If necessary, use the enclosed extension tubing (e) when using manual ventilation systems.
- ▶ During mechanical ventilation ensure for proper strain-relief.
- ▶ Ventilate the patient. Verify correct placement and successful ventilation at regular intervals by a suitable method.
- ▶ Place the necktape (d) around the patient's neck and connect it firmly to the flange.

- ▶ Check the tube for secure hold.
- Tracheal suctioning is possible through the connector of the extension tubing with a suction catheter  $\leq$  16 Fr.
- ▶ Prior to extubation evacuate the cuff completely with a syringe.
- The tube is radiopaque.

#### SHELF LIFE

The shelf life of the product is 5 years from the date of manufacture.

Expiration date: See product label

#### STORAGE AND TRANSPORT CONDITIONS

##### CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

#### DISPOSAL

##### CAUTION

- The product may be contaminated with potentially infectious substances of human origin.
- Risk of injury, needle tip and scalpel are sharp.

The used or faulty product must be disposed of in accordance with the applicable national and international regulations.

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